



**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION  
DECISION SUMMARY  
ASSAY AND INSTRUMENT**

**I Background Information:**

**A 510(k) Number**

K203035

**B Applicant**

PerkinElmer Inc

**C Proprietary and Established Names**

Eonis SCID-SMA kit

**D Regulatory Information**

Product Code(s)	Classification	Regulation Section	Panel
PJI	Class II	21 CFR 866.5980 - Newborn Screening Test for Severe Combined Immunodeficiency Disorder (SCID)	IM - Immunology

**II Submission/Device Overview:**

**A Purpose for Submission:**

New Device

**B Measurand:**

T-cell receptor excision circles (TRECs) and Kappa-deleting element recombination circles (KRECs)

**C Type of Test:**

Multiplex real-time polymerase chain reaction-based assay using target sequence specific primers and probes

### **III Intended Use/Indications for Use:**

#### **A Intended Use(s):**

See Indications for Use below.

#### **B Indication(s) for Use:**

The Eonis™ SCID-SMA kit is intended for the semi-quantitative determination of TREC (T-cell receptor excision circle) as an aid in screening newborns for Severe Combined Immunodeficiency (SCID) and for the semi-quantitative determination of KREC (Kappa-deleting recombination excision circle) as an aid in screening newborns for X-linked agammaglobulinemia (XLA). The test is intended for DNA from blood specimens dried on a filter paper and for use on the QuantStudio™ Dx Real-Time PCR instrument.

This test is not intended for screening of SCID-like Syndromes, such as DiGeorge Syndrome, or Omenn Syndrome. It is also not intended to screen for less acute SCID syndromes such as leaky-SCID or variant SCID. The test is not indicated for screening B-cell deficiency disorders other than XLA, such as atypical XLA, or for screening of XLA carriers.

This test is not intended for use as a diagnostic test and a positive screening result should be followed by confirmatory testing.

#### **C Special Conditions for Use Statement(s):**

Rx - For Prescription Use Only

In vitro diagnostic use

#### **D Special Instrument Requirements:**

QuantStudio™ Dx Real-Time PCR Instrument (ThermoFisher, k123955)

Liquid Handling Instruments: JANUS G3 Automated Workstation Extraction instrument (labeled for IVD use) or JANUS Extraction Mini instrument (labeled for IVD use) and the JANUS PCR Mastermix instrument (labeled for IVD use).

### **IV Device/System Characteristics:**

#### **A Device Description:**

The Eonis SCID-SMA kit contains reagents to detect three biomarkers: TREC, KREC and exon 7 in the SMN1 gene. Detection of SMN1 gene was submitted and reviewed in DEN200044. This Decision Summary only discusses the validation data in support of TREC and KREC.

The newborn screening workflow for the Eonis SCID-SMA kit includes:

- Two liquid handling platforms (one for DNA extraction and one for PCR master mix setup)
- QuantStudio™ Dx Real-Time PCR instrument (k123955)
- Eonis Analysis Software

Each Eonis SCID-SMA kit contains reagents for up to 384 reactions or 1152 reactions including kit controls. The kit contents are listed in Table 1. Materials required but not provided include the Eonis DNA Extraction Kit, Eonis Analysis Software, and consumables.

**Table 1. Eonis SCID-SMA Kit content**

Component	Quantity
SCID-SMA Kit Controls	2 filter paper cassettes containing 4 sets of dried blood spots for 384 reaction kit
	4 filter paper cassettes containing 8 sets of dried blood spots for 1152 reaction kit
C1 Analyte-negative (TREC/KREC/SMN1) control	
C2 Low TREC/KREC and normal SMN1 control	
C3 High TREC/KREC and normal SMN1 control	
PCR Reagent 1	1 vial, 2.7 mL for 384 reaction kit 3 vials, 2.7 mL for 1152 reaction kit
PCR Reagent 2	1 vial, 2.7 mL for 384 reaction kit 3 vials, 2.7 mL for 1152 reaction kit
Lot-specific quality control certificate	1 pc

## B Principle of Operation:

The Eonis SCID-SMA kit is a multiplex real-time PCR based assay that uses target sequence-specific primers and TaqMan™ probes to amplify and detect TREC, KREC, and RPP30 in DNA extracted from newborn dried blood spot (DBS) using the Eonis DNA Extraction kit in a single PCR reaction. Each TaqMan probe has a unique dye linked to their terminal 5' end, allowing for the simultaneous detection of the four targets, if present. The amount of each target present in the DNA is determined by the intensity of fluorescence emitted by each dye released from the degraded probe during amplification and detected by the QuantStudio Dx Real-Time PCR instrument. The instrument measures fluorescence signals and converts them into comparative quantitative readouts which are expressed as a function of cycle threshold (Ct) values. RPP30 is used as an internal amplification control for DNA quality and to relatively quantitate TREC and KREC (copies/10<sup>5</sup> cells) using TREC ΔCt and KREC ΔCt. TREC and KREC concentrations are calculated using the following formulas to generate results expressed as copies/10<sup>5</sup> cells:

- TREC:  $2 \times 2^{-(\Delta C_t \text{ TREC})} \times 117000$
- KREC:  $2 \times 2^{-(\Delta C_t \text{ TREC})} \times 254000$

The interpretation of results based on copy number is shown in Tables 2 and 3.

**Table 2. TREC result interpretation matrix**

RPP30 Ct values		TREC Ct values		TREC conc. (copies/10 <sup>5</sup> cells)	Result interpretation for TREC
Ct < 15.0				Values not considered when RPP30 Ct < 15.0	Invalid*
15.0 ≤ Ct ≤ 32.00	AND	Ct < 18.0		Values not considered when TREC Ct < 18.0	Invalid*
	AND	Ct ≥ 18.0	AND	TREC conc. ≥ CutoffTREC	Presumptive normal
			AND	TREC conc. < CutoffTREC	Presumptive positive*
	AND	No Ct	N/A		Presumptive positive*
Ct > 32.00 or no Ct				Values not considered when RPP30 Ct > 32.00 or no Ct	Invalid*

\*Samples with invalid results should be repeated in singlicate, and samples with presumptive positive results should be repeat tested in duplicate. If one of the duplicate samples is positive, the result is presumptive positive.

**Table 3. KREC result interpretation matrix**

RPP30 Ct values		KREC Ct values		KREC conc. (copies/10 <sup>5</sup> cells)	Result interpretation for KREC
Ct < 15.0				Values not considered when RPP30 Ct < 15.0	Invalid*
15.0 ≤ Ct ≤ 32.00	AND	Ct < 18.0		Values not considered when TREC Ct < 18.0	Invalid*
	AND	Ct ≥ 18.0	AND	KREC conc. ≥ CutoffKREC	Presumptive normal
			AND	KREC conc. < CutoffKREC	Presumptive positive*
	AND	No Ct	N/A		Presumptive positive*
Ct > 32.00 or no Ct				Values not considered when RPP30 Ct > 32.00 or no Ct	Invalid*

\*Samples with invalid results should be repeated in singlicate, and samples with presumptive positive results should be repeat tested in duplicate. If one of the duplicate samples is positive, the result is presumptive positive.

### *Result Reporting*

Samples that result in values below 242 copies/10<sup>5</sup> cells are recommended to be reported as "<242 copies/10<sup>5</sup> cells". Samples that result in values above 4316 copies/10<sup>5</sup> cells are recommended to be reported as ">4320 copies/10<sup>5</sup> cells" and are considered to be TREC screen negative.

Samples that result in values below 459 copies/10<sup>5</sup> cells are recommended to be reported as "<459 copies/10<sup>5</sup> cells". Samples that result in values above 24343 copies/10<sup>5</sup> cells are recommended to be reported as ">24300 copies/10<sup>5</sup> cells" and are considered to be KREC screen negative.

## **C Instrument Description Information:**

### 1. Instrument Name:

QuantStudio Dx Real-Time PCR instrument (The information below is specific to the TREC and KREC analytes; refer to the decision summary for additional information regarding the features in this section K123955)

JANUS G3 Automated Workstation (for IVD use) consisting of the JANUS Extraction Instrument and JANUS PCR Mastermix Instrument

### 2. Specimen Identification:

Specimens are barcoded and entered into the testing queue consistent with the laboratory procedures.

### 3. Specimen Sampling and Handling:

The required sample is a 3.2 mm specimen punch from a dried blood spot (DBS) into a 96-well plate. DNA is extracted and purified using the Eonis DNA Extraction kit. The extraction can be completed manually or using an automated liquid handling platform, the JANUS G3 Workstation (IVD), for processing up to four 96-well plates simultaneously.

### 4. Calibration:

The device (Eonis SCID-SMA kit) calibrates the TREC results utilizing the internal control RPP30 gene (assumes two copies) and a mathematical algorithm and reports out the results in the unit of copies/10<sup>5</sup> cells.

### 5. Quality Control:

Quality control is addressed for each separately cleared specific assay to be run on the instrument. For the Eonis SCID-SMA assay, quality is maintained by the process control, used during sample processing and amplification in the assay. Additionally, there are external positive and negative controls available that are used in accordance with the assay instructions for use.

Controls are used throughout the entire procedure from sample preparation to PCR as a quality check for performance of the Eonis test. The Eonis SCID-SMA kit contains three multi-level dried blood spot controls: C1, C2 and C3. The DBS controls are filter paper cassettes with each control level spotted in duplicate.

- C1 is an analyte-negative control and does not contain TREC, KREC analytes.
- C2 and C3 are analyte-positive controls that have different levels of TREC and KREC analytes (lower level in C2 and a higher level in C3).
- All three kit controls have the same level of RPP30 analyte.
- Blank reaction wells without sample disks (NTC [no template control]) are used to monitor the absence of contamination. All controls are prepared from human whole blood with a hematocrit value of 48-55% and contain purified salmon-sperm DNA and synthetic TREC, KREC, and RPP30 plasmid DNA.

Each Eonis run includes all four controls in duplicate in each 96-well DNA extraction plate. The newborn samples and DBS controls are punched into the plate, while the well reserved for the Blank control does not contain any punches. The Blank control undergoes the whole process, from DNA extraction to amplification, and the NTC sample consists of elution solution derived from the DNA extraction process.

The test result validation for each 96-well DNA extraction plate is based on the three kit DBS controls and the Blank reaction well. The Eonis Analysis software will identify each control by their designated locations and check their validity. The possible outcome for each control replicate is 'Valid' or 'Invalid'. If one or more controls are invalid, the results of the clinical specimens processed in the same 96-well DNA extraction plate are reported as 'Invalid'. For the Blank reaction well, both replicates must meet the acceptance criteria for the run to be valid. For C1-C3 DBS controls, the results of both replicates must meet the acceptance criteria for the run to be valid. The thresholds for valid controls are shown in Table 4.

**Table 4. Acceptance limits for kit controls for TREC and KREC**

<b>Kit Control Name</b>	<b>Parameter Name</b>	<b>Lower Limit (≥)</b>	<b>Upper Limit (≤ or No Ct)</b>
NTC (Blank)	TREC Ct	32.65	No Ct
	KREC Ct	32.55	No Ct
	RPP30 Ct	28.35	No Ct
C1	TREC Ct	32.65	No Ct
	KREC Ct	32.55	No Ct
	RPP30 Ct	15.00	28.00
C2	TREC Ct	18.00	No Ct
	KREC Ct	18.00	No Ct
	RPP30 Ct	15.00	28.00
	TREC ΔCt	6.10	13.90
	KREC ΔCt	3.55	16.45

Kit Control Name	Parameter Name	Lower Limit (≥)	Upper Limit (≤ or No Ct)
C3	TREC Ct	18.00	33.00
	KREC Ct	18.00	34.00
	RPP30 Ct	15.00	28.00
	TREC ΔCt	1.10	8.90
	KREC ΔCt	0.00	12.45

**V Substantial Equivalence Information:**

**A Predicate Device Name(s):**

PerkinElmer EnLite™ Neonatal TREC Kit

**B Predicate 510(k) Number(s):**

DEN140010

**C Comparison with Predicate(s):**

Device & Predicate Device(s):	K203035	DEN140010
Device Trade Name	Eonis SCID-SMA Kit	EnLite Neonatal TREC Kit
<b>General Device Characteristic Similarities</b>		
<b>Intended Use/Indications For Use</b>	<p>The Eonis™ SCID-SMA kit is intended for the semi-quantitative determination of TREC (T-cell receptor excision circle) as an aid in screening newborns for Severe Combined Immunodeficiency (SCID).</p> <p>The test is intended for DNA from blood specimens dried on a filter paper and for use on the QuantStudio™ Dx Real-Time PCR instrument.</p> <p>This test is not intended for screening of SCID-like Syndromes, such as DiGeorge Syndrome, or Omenn Syndrome.</p>	<p>The EnLite™ Neonatal TREC Kit is an in vitro diagnostic device intended for the semi-quantitative determination of TREC (T-cell receptor excision circle) DNA in blood specimens dried on filter paper. The test is for use on the VICTOR™ EnLite instrument. The test is indicated for use as an aid in screening newborns for severe combined immunodeficiency disorder (SCID). This test is not intended for use as a diagnostic test for screening of SCID-like Syndromes, such as DiGeorge Syndrome,</p>

	<p>It is also not intended to screen for less acute SCID syndromes such as leaky-SCID or variant SCID.</p> <p>This test is not intended for use as a diagnostic test and a positive screening result should be followed by confirmatory testing.</p>	<p>or Omenn Syndrome. It is also not intended to screen for less acute SCID syndromes such as leaky-SCID or variant SCID.</p>
<b>Specimen</b>	Dried blood spot	Same
<b>Stability/Shelf Life</b>	The kit should be stored at -30 to -16°C. Once opened, the kit should be used within 14 days.	Same
<b>General Device Characteristic Differences</b>		
<b>Intended Use/Indications For Use</b>	<p>For the semi-quantitative determination of KREC (Kappa-deleting recombination excision circle) as an aid in screening newborns for X-linked agammaglobulinemia (XLA).</p> <p>The test is not indicated for screening B-cell deficiency disorders other than XLA, such as atypical XLA, or for screening of XLA carriers.</p>	Not applicable
<b>Test Methodology</b>	Semi-quantitative, multiplex real-time fluorescent-based polymerase chain reaction (PCR) based nucleic acid amplification and detection	Semi-quantitative, polymerase chain reaction (PCR) based nucleic acid amplification and time-resolved fluorescence resonance energy transfer (TR-FRET) based detection
<b>Units of Reporting</b>	TREC results are calibrated using the internal control gene and a mathematical algorithm and reports out the results in the unit of copies/10 <sup>5</sup> cells.	TREC results are calibrated against a calibrator curve on each plate and reports out the results in the unit of copies/μL.
<b>Instrument/Software</b>	QuantStudio Dx Real-Time PCR instrument and Eonis software	VICTOR EnLite instrument and EnLite™ workstation software
<b>Reportable Range</b>	TREC 242 – 4320 copies/10 <sup>5</sup> cells	TREC 29-473 copies/μL blood



	KREC 459 - 24300 copies/10 <sup>5</sup> cells	
<b>Lower Limits of Measure</b>	TREC LoB=0 copies/10 <sup>5</sup> cells, LoD/LoQ= 242 copies/10 <sup>5</sup> cells KREC LoB=0 copies/10 <sup>5</sup> cells, LoD/LoQ =459 copies/10 <sup>5</sup> cells	TREC LoB=3 copies/ μL blood, LoD=20 copies/ μL blood, LoQ=29 copies/ μL blood
<b>Calibrators / Standards</b>	Calibration is based on internal reference (RPP30) in each well and manufacturer calibration for each kit lot.	3 levels of DBS calibrators prepared from porcine whole blood spiked with TREC and beta-actin (reference) plasmids, and a no-template blank
<b>Controls</b>	3 levels of DBS controls prepared from leucocyte-depleted human red blood cells and TREC, KREC and RPP30 plasmids spiked in, plus No template control (NTC)	3 levels of DBS controls prepared from porcine whole blood, TREC and beta-actin (reference) plasmids spiked in.

## VI Standards/Guidance Documents Referenced:

- CLSI EP05-A3: Evaluation of precision of quantitative measurement procedures; approved guideline-third edition
- CLSI EP6-A: Evaluation of the linearity of quantitative measurement procedures: a statistical approach; approved guideline
- CLSI EP07-A2: Interference testing in clinical chemistry: approved guideline-second edition
- CLSI EP25-A: Evaluation of Stability of In Vitro Diagnostic Reagents; approved guideline
- CLSI NBS01-A6: Blood Collection on Filter Paper for Newborn Screening Programs-sixth edition
- CLSI MM13-A Collection Transport Preparation and Storage of Specimens for Molecular Methods; approved guideline
- CLSI MM17 Verification and validation of multiplex nucleic acid assays-second edition
- IEC 62304 Edition 1.1 2015-06 Medical device software-Software life cycle processes

## VII Performance Characteristics (if/when applicable):

### A Analytical Performance:

#### 1. Precision/Reproducibility:

##### *Study 1: Site-to-Site Reproducibility*

The objective of this study was to characterize the within-laboratory precision and reproducibility of the Eonis test system across 3 study sites using a panel of samples with different TREC and KREC levels. At each study site, 2 operators performed 5 runs each, over 5 operating days, 1 run per day, and 5 specimen replicates tested per run. A total of 13

independent samples were tested at each of the 3 study sites using 1 kit lot for a total of 150 replicates per sample (3 sites x 2 operators x 5 runs/days x 5 replicates per run = 150). The study used 1 lot of reagents. The samples were generated to represent a range of concentrations and were spotted onto filter paper, dried overnight, and stored at -30°C to -16°C in a sealed bag with desiccant until use. A description of the samples is provided in Table 5.

**Table 5. Descriptive summary samples used in the precision study**

Sample	TREC	KREC	Preparation
<b>Sample 1</b>	High	Normal	TREC plasmid spiked into cord blood at 1500 copies/ $\mu$ L with the hematocrit adjusted to 48.2 (approx. Ct value 28.63)
<b>Sample 2</b>	Normal	High	KREC plasmid spiked into cord blood at 5000 copies/ $\mu$ L with the hematocrit adjusted to 48.2 (approx. Ct value 28.11)
<b>Sample 3</b>	Normal	Normal	Cord blood with hematocrit adjusted to 48.2
<b>Sample 4</b>	Low (Above cut-off)	Normal	Diluting 3mL of cord blood with hematocrit adjusted to 48.2 into 2.25 mL of Sample 6 targeting near but below 98% quantile (32.71) of TREC Ct distribution of a population of presumptive normal NBS samples
<b>Sample 5</b>	Low (Below cut-off)	Normal	Diluting 857.1 $\mu$ L of cord blood with hematocrit of 48.32 in 2142.9 $\mu$ L of Sample 6 (approx. Ct value 33.07)
<b>Sample 6</b>	Very low (Near 0)	Normal	KREC plasmid spiked into adult whole blood at 2000 copies/ $\mu$ L with hematocrit adjusted to 50.7 (approx. TREC Ct value 38.08; approx. KREC Ct value 30.56)
<b>Sample 7</b>	Normal	Low (Above cut-off)	Diluting 0.75 mL of cord blood with hematocrit adjusted to 48.2 in 2.25 mL of Sample 8 (approx. KREC Ct value 31.96)
<b>Sample 8</b>	Normal	Low (Below cut-off)	TREC plasmid spiked into adult whole blood at 700 copies/ $\mu$ L with hematocrit adjusted to 50.7 (approx. TREC value 30.01 and approx. KREC value 33.30)
<b>Sample 9</b>	Normal	Very low (Near 0)	TREC plasmid spiked at 2000 copies/ $\mu$ L with hematocrit adjusted to 49.5 (approx. TREC Ct value 30.08 and approx. KREC Ct value 0)
<b>Sample 10</b>	Normal	Normal	Pooled adult whole blood from 2 SMA carriers with hematocrit adjusted to 46.1
<b>Sample 11</b>	Normal	Normal	TREC plasmid at 700 copies/ $\mu$ L, KREC plasmid at 1000 copies/ $\mu$ L and SMA Coriell cells at 15000 cells/ $\mu$ L spiked into leukocyte depleted blood with hematocrit adjusted to 49.5

Sample	TREC	KREC	Preparation
<b>Sample 12</b>	Low (Above cut-off)	Low (Above cut-off)	Diluting 4.4 mL of cord blood with hematocrit adjusted to 48.2 into 6.6 mL of adult whole blood with hematocrit adjusted to 50.7 (approx. TREC Ct value 32.48 and approx. KREC Ct value 31.30)
<b>Sample 13</b>	Low (Below cut-off)	Low (Below cut-off)	Adult blood hematocrit adjusted to 50.7 (approx. TREC Ct value 36.34 and approx. KREC Ct value 33.89)

The values in the mean TREC and KREC copies/10<sup>5</sup> cells column are transformed from logarithmic (Ln) mean values, and therefore, they represent geometric means in the copies/10<sup>5</sup> cells scale. The analysis of variance approach was used to calculate the results presented as standard deviations (SDs) in the logarithmic (Ln) scale complemented with total coefficient of variation (%CV) in lognormal scale (Table 6 and Table 8). Summary mean, min and max copies/10<sup>5</sup> cells and SD and %CV results without log transformation for total imprecision are shown in Table 7 and Table 9.

Acceptance Criteria:

- For TREC, the total assay precision within the measuring range must have a SD ≤ 0.9Ln (copies/10<sup>5</sup> cells) where the measuring range starts at 250 copies/10<sup>5</sup> cells.
- For KREC, the total assay precision within the measuring range must have a SD ≤ 1.49 Ln (copies/10<sup>5</sup> cells).

**Table 6. TREC reproducibility data pooled across three laboratories**

Sample	Geometric Mean (Copies 10 <sup>5</sup> /cells)	Mean Ln (Copies 10 <sup>5</sup> /cells)	Repeatability Ln SD	Between Run Ln SD	Between Operator Ln SD	Between Site Ln SD	Total Ln SD	Log-normal CV%
6	23	3.14	0.73	0.60	0.01	0.53	1.09	151
13	34	3.53	1.12	0.31	0.00	0.32	1.20	180
10	180	5.19	0.77	0.01	0.00	0.03	0.77	90
5	299	5.70	0.55	0.16	0.01	0.15	0.59	65
3	912	6.82	0.29	0.06	0.00	0.10	0.31	32
12	936	6.84	0.34	0.00	0.00	0.18	0.38	40
4	1026	6.93	0.28	0.10	0.01	0.13	0.33	34
2	1478	7.30	0.26	0.01	0.04	0.04	0.27	27
8	3013	8.01	0.48	0.20	0.01	0.16	0.55	59
11	3290	8.10	0.25	0.05	0.00	0.09	0.27	28
7	3867	8.26	0.29	0.08	0.00	0.11	0.32	33
9	8225	9.01	0.24	0.08	0.02	0.18	0.31	32
1	9000	9.11	0.21	0.08	0.00	0.14	0.26	27

**Table 7. TREC total variation results without logarithmic transformation**

		Calculations without logarithmic transformation				
Sample	N	Mean (Copies/ 10 <sup>5</sup> cells)	Min (Copies/ 10 <sup>5</sup> cells)	Max (Copies/ 10 <sup>5</sup> cells)	SD (Copies/ 10 <sup>5</sup> cells)	CV%
6	54	40	4	308	55.3	138
13	71	58	3	239	56	97
10	148	223	5	796	126	56
5	150	346	41	971	174	50
3	150	952	220	1792	272	29
12	150	997	239	2563	349	35
4	150	1076	497	2490	335	31
2	150	1526	654	2411	370	24
11	150	3388	513	6401	767	23
7	150	4076	2043	16934	1590	39
8	150	4601	797	218936	17700	384*
9	150	8571	3233	14974	2420	28
1	150	9277	4236	16270	2320	25

\*Dataset for sample 8 has one high outlier (max value) affecting the variability estimate. Without the outlier, the estimated CV% is 41% and similar to other samples within the measuring range.

**Table 8. KREC reproducibility data pooled across three laboratories**

Sample	Geometric Mean (Copies 10 <sup>5</sup> /cells)	Mean Ln (Copies 10 <sup>5</sup> /cells)	Repeatability Ln SD	Between Run Ln SD	Between Operator Ln SD	Between Site Ln SD	Total Ln SD	Log-normal CV%
10	348	5.85	0.72	0.01	0.13	0.22	0.76	89
8	763	6.64	0.46	0.05	0.00	0.13	0.48	50
13	792	6.67	0.43	0.09	0.06	0.09	0.45	47
12	2770	7.93	0.24	0.00	0.00	0.17	0.30	31
3	3849	8.26	0.28	0.07	0.00	0.09	0.31	31
7	4037	8.30	0.23	0.02	0.05	0.11	0.26	27
1	4359	8.38	0.19	0.08	0.03	0.06	0.22	22
6	5007	8.52	0.38	0.01	0.00	0.06	0.38	40
4	9313	9.14	0.20	0.05	0.00	0.02	0.21	21
11	9648	9.17	0.25	0.05	0.00	0.05	0.26	26
5	16531	9.71	0.19	0.07	0.06	0.13	0.24	25
2	34982	10.5	0.27	0.00	0.01	0.01	0.27	28
10	348	5.85	0.72	0.01	0.13	0.22	0.76	89

**Table 9. KREC total variation results without logarithmic transformation**

		Calculations without logarithmic transformation				
Sample	N	Mean (Copies/ 10 <sup>5</sup> cells)	Min (Copies/ 10 <sup>5</sup> cells)	Max (Copies/ 10 <sup>5</sup> cells)	SD (Copies/ 10 <sup>5</sup> cells)	CV%
10	148	431	40	1071	242	56
8	150	841	186	1925	349	41
13	150	862	183	2003	330	38
12	150	2877	1097	5232	792	28
3	150	4025	1075	13447	1300	32
7	150	4162	2241	6788	1030	25
1	150	4457	2104	7746	960	22
6	150	5457	1567	36630	3230	59
4	150	9508	4735	15856	1910	20
11	150	9915	1564	16967	2130	21
5	150	16955	8813	29034	3880	23
2	150	36370	16841	118413	11200	31

*Study 2: Lot-to-Lot Reproducibility*

A study was performed internally (single site) using 3 sets of Eonis test systems, 2 operators, and 3 kit lots. The 13 samples (Table 5) were tested over 23 calendar days with 2 replicates per sample per plate in a total of 54 runs with a randomized plate map. The total number of measurements was 108 per sample.

The values in the mean TREC and KREC copies/10<sup>5</sup> cells column are transformed from logarithmic (Ln) mean values, and therefore, they represent geometric means in the copies/10<sup>5</sup> cells scale. The analysis of variance approach was used to calculate the results presented as SDs in the logarithmic (Ln) scale complemented with total %CVs in lognormal scale (Table 10 and Table 12). Summary mean, min and max copies/10<sup>5</sup> cells, and SD and %CV results without log transformation for total imprecision are shown in Table 11 and Table 13.

**Table 10. TREC precision data across three kit lots and three instruments**

Sample	Geometric Mean (Copies 10 <sup>5</sup> /cells)	Mean Ln (Copies 10 <sup>5</sup> /cells)	Repeatability Ln SD	Between Run Ln SD	Between Instrument Ln SD	Between Lot Ln SD	Total Ln SD	Log-normal CV%
6	20	3.00	0.92	0.00	0.00	0.33	0.97	125
13	34	3.54	0.89	0.00	0.00	0.23	0.92	116
10	159	5.07	0.76	0.17	0.00	0.14	0.79	92
5	464	6.14	0.42	0.06	0.14	0.13	0.47	49
12	1022	6.93	0.34	0.16	0.00	0.08	0.39	40
2	1130	7.03	0.35	0.20	0.00	0.02	0.40	42

Sample	Geometric Mean (Copies 10 <sup>5</sup> /cells)	Mean Ln (Copies 10 <sup>5</sup> /cells)	Repeatability Ln SD	Between Run Ln SD	Between Instrument Ln SD	Between Lot Ln SD	Total Ln SD	Log-normal CV%
4	1176	7.07	0.33	0.13	0.13	0.15	0.41	42
3	2165	7.68	0.37	0.00	0.00	0.05	0.37	38
11	4105	8.32	0.24	0.08	0.16	0.04	0.30	30
7	4146	8.33	0.32	0.15	0.10	0.01	0.37	38
8	4866	8.49	0.35	0.31	0.00	0.07	0.47	50
9	8604	9.06	0.30	0.35	0.19	0.18	0.54	58
1	11048	9.31	0.21	0.14	0.16	0.03	0.30	31

**Table 11. TREC total variation results without logarithmic transformation**

		Calculations without logarithmic transformation				
Sample	N	Mean (Copies/ 10 <sup>5</sup> cells)	Min (Copies/ 10 <sup>5</sup> cells)	Max (Copies/ 10 <sup>5</sup> cells)	SD (Copies/ 10 <sup>5</sup> cells)	CV%
6	29	32	5	150	34.9	109
13	43	49	4	173	41.4	84
10	105	200	11	479	114	57
5	107	508	99	1183	216	43
12	106	1094	189	2247	383	35
2	107	1206	230	2576	414	34
4	108	1259	422	3011	479	38
3	107	2263	107	3084	510	23
11	106	4263	2031	7418	1220	29
7	107	4407	1014	8206	1530	35
8	107	5754	2011	74368	7120	124
9	107	9540	471	19530	3760	39
1	107	11501	3218	19982	3170	28

**Table 12. KREC precision data across three kit lots and three instruments**

Sample	Geometric Mean (Copies 10 <sup>5</sup> /cells)	Mean Ln (Copies 10 <sup>5</sup> /cells)	Repeatability Ln SD	Between Run Ln SD	Between Instrument Ln SD	Between Lot Ln SD	Total Ln SD	Log-normal CV%
9	85	4.44	1.14	0.00	0.00	0.37	1.20	178
10	478	6.17	0.53	0.36	0.00	0.16	0.66	74
13	1033	6.94	0.40	0.13	0.15	0.09	0.45	48
8	1075	6.98	0.47	0.00	0.00	0.13	0.49	52
7	2416	7.79	0.43	0.00	0.08	0.01	0.44	46
12	3361	8.12	0.36	0.10	0.00	0.07	0.38	39

Sample	Geometric Mean (Copies 10 <sup>5</sup> /cells)	Mean Ln (Copies 10 <sup>5</sup> /cells)	Repeatability Ln SD	Between Run Ln SD	Between Instrument Ln SD	Between Lot Ln SD	Total Ln SD	Log-normal CV%
1	5271	8.57	0.35	0.00	0.09	0.04	0.37	38
3	5541	8.62	0.30	0.05	0.00	0.03	0.31	31
6	8691	9.07	0.26	0.26	0.08	0.10	0.39	40
4	11499	9.35	0.29	0.18	0.12	0.06	0.36	38
11	13494	9.51	0.24	0.08	0.17	0.09	0.32	32
5	13767	9.53	0.22	0.21	0.13	0.07	0.34	35
2	39735	10.6	0.20	0.21	0.12	0.07	0.32	33

**Table 13. KREC total variation results without logarithmic transformation**

		Calculations without logarithmic transformation				
Sample	N	Mean (Copies/ 10 <sup>5</sup> cells)	Min (Copies/ 10 <sup>5</sup> cells)	Max (Copies/ 10 <sup>5</sup> cells)	SD (Copies/ 10 <sup>5</sup> cells)	CV%
9	10	146	17	435	139	95
10	104	570	43	1396	306	54
13	106	1139	212	4188	546	48
8	107	1195	356	3180	538	45
7	107	2605	203	5513	898	34
12	107	3568	686	7096	1150	32
1	107	5560	474	10461	1620	29
3	106	5774	863	10031	1500	26
6	107	9308	3333	19079	3420	37
4	108	12212	2707	23846	4180	34
11	107	14103	5330	23372	4030	29
5	106	14492	5673	35070	4960	34
2	108	41717	15231	74089	12400	30

All 84 plates from both studies had results of the Eonis SCID-SMA Kit Controls within the acceptance limits. The data validity-check for each well resulted in a total of 15 out of 3354 well exclusions from the analysis (sample invalid rate is 0.4%). The precision of the device control kits is shown in Table 14a and the internal control RPP30 gene in Table 14b below.

**Table 14a. Descriptive statistics of device kit control results (Ct values) at one site**

Analyte	Sample	N	N Ct N/A	Ct Mean	Ct SD	Ct CV%	Ct Min	Ct Max
KREC	C1	20	20	N/A	N/A	N/A	NaN	NaN
	C2	20	0	33.5	0.56	1.68	32.3	34.6
	C3	20	0	29.5	0.55	1.87	28.8	31.1
	NTC	20	20	N/A	N/A	N/A	NaN	NaN
KREC ΔCt	C1	20	20	N/A	N/A	N/A	NaN	NaN

Analyte	Sample	N	N Ct N/A	Ct Mean	Ct SD	Ct CV%	Ct Min	Ct Max
	C2	20	0	10.1	0.53	5.22	8.98	11.2
	C3	20	0	6.07	0.43	7.10	5.48	7.20
	NTC	20	20	N/A	N/A	N/A	NaN	NaN
RPP30	C1	20	0	22.2	0.49	2.23	21.5	23.1
	C2	20	0	21.7	0.37	1.69	21.1	22.4
	C3	20	0	21.9	0.52	2.35	21.1	23.4
	NTC	20	19	37.5	NA	N/A	37.5	37.5
TREC	C1	20	20	N/A	N/A	N/A	NaN	NaN
	C2	20	0	33.4	0.77	2.31	32.3	35.4
	C3	20	0	28.6	0.64	2.25	27.7	30.5
	NTC	20	20	N/A	N/A	N/A	NaN	NaN
TREC ΔCt	C1	20	20	N/A	N/A	N/A	NaN	NaN
	C2	20	0	33.4	0.77	2.31	32.3	35.4
	C3	20	0	28.6	0.64	2.25	27.7	30.5
	NTC	20	20	N/A	N/A	N/A	NaN	NaN

N/A not applicable; NaN = Not a number (i.e., no Ct value generated)

**Table 14b. RPP30 precision analysis on Ct values**

Sample	N	Mean	Within Run		Between Run		Between Instrument		Within Kit lot		Between Kit lot		Total	
			SD	CV%	SD	CV%	SD	CV%	SD	CV%	SD	CV%	SD	CV%
1	107	23.4	0.70	3.0	0.39	1.7	0.30	1.3	0.42	1.8	0.86	3.7	0.96	4.1
2	108	23.2	0.54	2.3	0.44	1.9	0.23	1.0	0.34	1.5	0.73	3.2	0.81	3.5
3	106	23.2	0.69	3.0	0.25	1.1	0.09	0.4	0.34	1.5	0.74	3.2	0.81	3.5
4	108	24.2	0.99	4.1	0.33	1.4	0.36	1.5	0.45	1.9	1.11	4.6	1.20	4.9
5	106	24.0	0.79	3.3	0.00	0.0	0.22	0.9	0.54	2.2	0.82	3.4	0.98	4.1
6	107	24.4	0.59	2.4	0.24	1.0	0.00	0.0	0.70	2.9	0.64	2.6	0.95	3.9
7	107	24.2	1.17	4.8	0.00	0.0	0.19	0.8	0.64	2.7	1.19	4.9	1.35	5.6
8	107	24.4	0.81	3.3	0.00	0.0	0.20	0.8	0.70	2.9	0.84	3.4	1.09	4.5
9	106	25.6	0.69	2.7	0.00	0.0	0.00	0.0	0.37	1.4	0.69	2.7	0.78	3.0
10	107	23.9	0.62	2.6	0.14	0.6	0.34	1.4	0.56	2.4	0.72	3.0	0.92	3.8
11	107	23.0	0.45	2.0	0.34	1.5	0.00	0.0	0.40	1.8	0.56	2.4	0.69	3.0
12	107	23.7	0.81	3.4	0.28	1.2	0.25	1.1	0.55	2.3	0.89	3.8	1.05	4.4
13	106	24.5	0.81	3.3	0.50	2.1	0.2	0.8	0.72	2.9	0.97	4.0	1.21	5.0

Qualitative imprecision was determined in accordance with CLSI document EP12-A2 to classify the results, using the screening performance study cut-offs (262 copies/10<sup>5</sup> cells for TREC, 484 copies/10<sup>5</sup> cells for KREC). The C5-C95 interval was determined to be 79–626 copies/10<sup>5</sup> cells for TREC, and 189–1064 copies/10<sup>5</sup> cells for KREC. Concentrations outside of these intervals were consistently analyte negative (concentrations <C5, equals to screen positive) or consistently analyte positive (concentrations >C95, equals to screen negative).



2. Linearity:

Linearity was determined with 1 kit lot and 1 operator. Three (3) sets of contrived samples were used for this evaluation and are described below. The samples were diluted to 9 levels with 4 replicates for each level. Samples were punched at 4 replicates using a DBS puncher and then transferred into a new 96-well plate in a randomized order and were extracted and run using JANUS automated liquid handlers. Acceptance criteria were that all analytes, TREC, KREC and RPP30, must have PCR efficiency between 80% and 120% with a coefficient of determination, denoted  $R^2 \geq 95\%$  within their measuring ranges. The acceptance criteria for linearity for TREC and KREC concentrations were, on their measuring ranges, the maximal difference of the best fitting model from the linear model must be  $\pm 25\%$  (calculated in copies/ $10^5$  cells).

*Sample Set 1: Estimation of PCR Efficiency for TREC, KREC and RPP30*

The first set of samples were prepared by spiking leucocyte-depleted human red blood cells, which had the hematocrit adjusted to 51.9, with linearized TREC and KREC plasmids and cell lines covering the required measuring ranges (TREC: 52 to 7800 copies/ $\mu\text{L}$  blood; KREC: 67 to 10,100 copies/ $\mu\text{L}$  blood; RPP30: 2670 to 400,000 copies/ $\mu\text{L}$  blood, per spiked cell numbers), a semi-log linear regression line plot of Ct value vs. log of input nucleic acid in the unit of copies/ $\mu\text{L}$  blood was first analyzed for each individual analyte (TREC, KREC and RPP30) and the corresponding coefficient of determination ( $R^2$ ) of the linear regression of the mean Ct values and PCR efficiency of each analyte were calculated. As shown in Table 15, all 3 analytes had  $R^2 > 99\%$ , and their estimates of PCR efficiency were between 80% and 120%.

**Table 15. Estimate of  $R^2$  and PCR efficiency of TREC, KREC and RPP30**

Analyte Name	$R^2$	PCR efficiency estimate	Lower 95% CL	Upper 95% CL
TREC	0.993	83.9%	76.7%	92.6%
KREC	0.995	87.6%	81.2%	95.0%
RPP30	0.992	99.9%	89.9%	112.3%

*Sample Set 2 and 3: Estimation of Linearity of the TREC and KREC Assay*


Linearity was determined in accordance with CLSI document EP06-A2 and was designed to evaluate the linearity of TREC and KREC in the unit of copies/ $10^5$  cells calculated using  $\Delta\text{Ct}$  TREC and  $\Delta\text{Ct}$  KREC for each dilution. A high concentration sample was prepared by resuspending WBC (White Blood Cells) pellets concentrated from the 60 mL of human umbilical cord bloods into 1mL of human umbilical cord blood to achieve TREC and KREC Ct close to the lowest values (i.e., high concentrations) observed in a population of newborn specimens (28.34. for TREC and 26.91 for KREC). Samples were then diluted with leukocyte depleted whole blood into nine levels of serial dilutions down to a level close to their LoQs (242 copies/ $10^5$  cells for TREC and 459 copies/ $10^5$  cells for KREC). In addition, a third set of samples were created to provide additional information on the KREC linearity at lower KREC concentrations because the second set did not cover the measuring range. This low concentration sample set 3 was prepared by diluting umbilical cord blood with

leucocyte-depleted human red blood. HL-60 cells were spiked into the leucocyte-depleted blood to adjust the RPP30 Ct values to match the RPP30 Ct value in cord blood. Using the observed value for undiluted sample and the consecutive dilution factors, the expected values at copies/10<sup>5</sup> cells scale were calculated. Linearity of TREC and KREC was then determined by fitting regression of Ln Measured vs. Ln Expected value. Weights for regression fit were calculated from the observed variability of each level. Using the obtained equation, the predicted values at each dilution level and the deviations between measured and predicted values were calculated. The descriptive statistics of each linearity sample level are presented in Table 16 and Table 17. The data for KREC represents the combined datasets. The measured results and fitted models are shown in Figure 1 and 3 the relative deviations in Figures 2 and 4. For TREC, all points were within the study acceptance criteria of ±25%, with highest relative deviation of -14.3%. For KREC, all points were within the study acceptance criteria of ±25%, with highest relative deviation of 17.3%.


**Table 16. Descriptive statistics and calculated weights - TREC**

RC*	N	Copies/10 <sup>5</sup> Mean	Ln Copies Mean	Copies/10 <sup>5</sup> SD	Ln Copies SD	Weight	Copies/10 <sup>5</sup> Min	Copies/10 <sup>5</sup> Max
1.000	4	4316	8.37	461	0.11	331	3790	4890
0.690	4	2441	7.80	834	0.37	29.2	1430	3330
0.476	3	1588	7.37	459	0.31	31.2	1110	1910
0.328	4	1353	7.21	257	0.19	111	1100	1710
0.226	4	1287	7.16	260	0.20	100	999	1630
0.156	4	508	6.23	148	0.34	34.6	309	630
0.108	4	185	5.22	113	0.71	7.93	66.5	340
0.074	3	114	4.74	64.5	0.48	13.0	82.8	198
0.051	4	93.7	4.54	77.8	0.78	6.57	40.8	190

\*Relative concentration



**Figure 1. Fitted linear model TREC**



**Figure 2. Relative deviations and the study limits – TREC**

**Table 17. Descriptive statistics and calculated weights - KREC**

RC*	Sample	N	Hit rate	Copies/10 <sup>5</sup> cells Mean	Ln Copies Mean	Ln Copies SD	Weight	Copies/10 <sup>5</sup> cells Min	Copies/10 <sup>5</sup> cells Max
1.000	Lin	4	1	24343	10.1	0.14	204	19700	26800
0.690	Lin	4	1	15835	9.67	0.19	111	12500	19500
0.476	Lin	3*	1	10405	9.25	0.03	3330	10100	10800
0.328	Lin	4	1	8519	9.05	0.05	1600	7990	9020
0.226	Lin	4	1	6768	8.82	0.08	625	6290	7470
1.000	L2	20	1	6438	8.77	0.26	296	4170	11500
	L3	20	1	5064	8.53	0.29	238	2380	8880
	L1	20	1	3905	8.27	0.36	154	2500	8160
0.156	Lin	4	1	3533	8.17	0.45	19.8	2170	5800
0.700	L1	20	1	3072	8.03	0.33	184	1440	5130
0.108	Lin	4	1	2836	7.95	0.34	34.6	1720	3750
0.600	L1	20	1	2644	7.88	0.32	195	1500	4650
0.074	Lin	3*	1	2490	7.82	0.16	117	2090	2790
0.500	L1	20	1	2080	7.64	0.43	108	567	4320
0.300	L3	20	1	1588	7.37	0.50	80.0	375	3540
0.051	Lin	4	1	1249	7.13	0.37	29.2	780	1910
0.200	L1	19	0.95	1188	7.08	0.36	147	619	2200
0.150	L2	20	1	871	6.77	0.53	71.2	356	2850
	L3	20	1	804	6.69	0.64	48.8	82.9	1570
0.100	L2	19	0.95	757	6.63	0.47	86.0	344	2970
0.080	L2	20	1	478	6.17	0.92	23.6	16.2	1580
0.040		14	0.7	305	5.72	1.03	13.2	54.9	1030
0.080	L3	18	0.9	296	5.69	0.68	38.9	74.1	936
	L1	16	0.8	255	5.54	0.63	40.3	75.2	620

RC*	Sample	N	Hit rate	Copies/10 <sup>5</sup> cells Mean	Ln Copies Mean	Ln Copies SD	Weight	Copies/10 <sup>5</sup> cells Min	Copies/10 <sup>5</sup> cells Max
0.030	L2	13	0.65	226	5.42	0.68	28.1	52.9	480
0.020	L3	4	1	200	5.30	0.48	17.4	152	411
0.060	L3	18	0.9	187	5.23	0.97	19.1	29.8	804
0.030	L3	14	0.7	178	5.18	0.57	43.1	75.7	407
0.010	L3	9	0.45	174	5.16	0.64	22.0	80.1	608
0.040	L3	17	0.85	166	5.11	1.03	16.0	10.5	392
0.020	L2	13	0.65	151	5.02	0.90	16.0	22.2	502
0.010	L2	9	0.45	147	4.99	0.34	77.9	61.7	190
0.030	L1	13	0.65	129	4.86	0.65	30.8	44.3	466
0.010	L1	3	0.15	117	4.76	0.56	9.57	61.8	166

\*Relative Concentration




Figure 3. Fitted linear model – KREC




Figure 4. Relative deviations and the study limits – KREC

The TREC analyte is demonstrated to be linear from 94 copies/10<sup>5</sup> cells to 4316 copies/10<sup>5</sup> cells with observed maximum deviation of -14.3%. The KREC analyte is demonstrated to be linear from 117 to 24343 copies/10<sup>5</sup> cells with observed maximum deviation of 17.3%.

### 3. Analytical Specificity

#### (i) *Interference*

The Eonis SCID-SMA kit was evaluated for interference from the following potential endogenous and exogenous sources:

- Hemoglobin up to 200g/L
- Conjugated bilirubin up to 16.6 mg/dL (166mg/L)
- Unconjugated bilirubin up to 10mg/dL (100mg/L)
- Intralipid up to 1500 mg/dL (15mg/mL)
- Li-hep up to 7500 USP/dL
- Na Citrate up to 0.0645 mol/L
- EDTA up to 9.8 mg/mL

The interference study was performed by 3 operators using 1 Eonis SCID-SMA Kit and 3 QuantStudio Dx Real-time PCR instruments. Three levels of TREC and KREC target DNA (normal, near cut-off and endogenous were tested in conjunction with each level of interferent (low and high). Blood samples were collected from healthy individuals and the hematocrit was adjusted to 40-55% corresponding to the hematocrit of neonates (except blood used in the hemoglobin sample series). Blood samples were then spiked with TREC and KREC plasmid DNA to obtain the desired concentration levels. Blood samples with the added potential interferent were used to prepare dried blood spots. The blood spots were dried overnight and stored at -80°C. The specimens used are shown in Table 18.

**Table 18. Samples used in the interference study**

<b>Sample Name (Abbreviation)</b>	<b>TREC Level</b>	<b>KREC Level</b>	<b>Interferent Substance Level</b>
Normal Targets+ Low Interferent (NL)	Normal	Normal	Low (0%)
Normal Targets+ High Interferent (NH)	Normal	Normal	High (100%)
Cut-off Targets + Low Interferent (CL)	Near Cut-off	Near Cut-off	Low (0%)
Cut-off Targets + High Interferent (CH)	Near Cut-off	Near Cut-off	High (100%)
Endogenous Targets + Low Interferent (EL)	Endogenous	Endogenous	Low (0%)
Endogenous Targets + High Interferent (EH)	Endogenous	Endogenous	High (100%)

The study generated a total of 544 sample results not including controls (7 interfering substances x 2 interferent levels x 3 target DNA levels x 13 repeats). Acceptance criteria

were based on observed imprecision in the precision studies. The results indicated there was no interference observed for any of the tested interferents. Hemoglobin at up to 200 g/L, conjugated bilirubin at up to 16.6 mg/dL (=166 mg/L), unconjugated bilirubin at up to 10 mg/dL (=100 mg/L), Intralipid at up to 1500 mg/dL (15 mg/mL), Li-heparin at up to 7500 USP/dL, Na-citrate at up to 0.0645 mol/L, and EDTA at up to 9.8 mg/mL did not significantly interfere with the Eonis SCID-SMA kit.

(ii) *Primer and Probe Specificity*

An *in silico* study was performed to assess the specificity of the primers and probes used in the Eonis assay. Primer and probe specificity was assessed by using the online Standard Nucleotide BLAST tool (Optimize for highly similar sequences, megablast) against the Human genomic plus transcript (Human G + T) database published by the National Center for Biotechnology Information (NCBI).

The BLAST analysis confirmed that there were only hits with 100% homology to the T receptor alpha delta locus (TREC) on chromosome 14 along, the full length of the primer and probe sequences, and 100% homology to the immunoglobulin kappa locus (IGK) on chromosome 2, and to the ribonuclease P subunit p30 locus (RPP30) on chromosome 10. There was no significant homology to other genome sequences.

4. Assay Reportable Range:

Refer to the Results interpretation section.

5. Traceability, Stability, Expected Values (Controls, Calibrators, or Methods):





(i) *Traceability*

There are no recognized standards or reference material for TREC or KREC DNA. The kit contains controls employed in the determination of TREC and KREC values.




(ii) *Specimen Stability*

The influence of storage time, temperature, and humidity was studied using DBS samples prepared with two brands of filter paper. Samples were prepared by mixing cell lines and red blood cell concentrate (mimicking a true positive sample for all 3 conditions), using cord blood (mimicking a true negative sample) or by mixing cell lines, red blood cell concentrate and cord blood for the near cut-off level samples.


After collection and preparation, the blood spots were dried overnight at ambient temperature (+19–+25°C). The dried blood spots were then stored at different storage conditions. Storage of specimens in an environment with elevated temperatures and humidity increases the risk of inaccurate TREC/KREC results in normal or near cut-off samples. The TREC/KREC results are shown in Figure 5a and Figure 5b.



**Figure 5a. TREC copies across different storage conditions in samples mimicking true positive (Sample 1), true negative (Sample 4) and in-between (Sample 2 and Sample 3)**







**Figure 5b. KREC copies across different storage conditions in samples mimicking true positive (Sample 1), true negative (Sample 4) and in-between (Sample 2 and Sample 3)**

*(iii) Reagent Stability:*

Stability studies were conducted to demonstrate the real time, in-use, and on-board stability as well as transport shipping stress simulation of the Eonis DNA Extraction Kit components and the Eonis SCID-SMA kit. The tests were performed after a transport simulation treatment that simulates the worst-case conditions that can occur during packaging and shipping. When a transport simulation treatment was performed, all the kit components were used in the study. The stability testing was conducted in accordance with guideline EN ISO 23640 In vitro diagnostic medical devices – Evaluation of stability of in vitro diagnostic reagents (2015) and CLSI Guideline EP25-A Evaluation of Stability of In Vitro Diagnostic Reagents, Approved Guideline- Volume 29, No. 20 (2009).

The testing was performed at 7 months of the components' shelf life. All Mean Ct values and variation (standard deviation) results fulfilled the acceptance criteria at all tested time points. At all tested time points, the results passed the data validity check. (Refer to DEN200044 for details of the study design.) The following stability conditions apply based on the results:

- Eonis DNA Extraction Kit components (Wash Solution and Elution Solution) are stable for 14 days at +19 - +25°C after first use.
- Eonis DNA extraction Kit components are stable in troughs for 30 minutes at +19 - +25°C.
- Note: The Eonis DNA extraction Kit components are stable on board the JANUS Extraction instrument in troughs for 2 hours.
- PCR Reagents 1 and 2 once thawed from -30°C to -16°C, are stable for 14 days at +2°C to +8°C.
- PCR Reagents 1 and 2 are stable in their tubes uncapped on-board the JANUS PCR Mastermix instrument for 2 hours at +19 to +25°C.
- The mixture of PCR Reagents 1 and 2 is stable for 15 minutes in a trough at +19 to +25°C.
- The mixture of PCR Reagents 1 and 2 is stable for 2 hours in a sealed plate protected from light.
- The SCID-SMA Kit controls are stable for 14 days after first use when stored at -30°C to -16°C in a sealed bag with desiccant.

- The SCID-SMA Kit controls are stable on-board the JANUS Extraction instrument at +19 to +25°C for 4 hours after they are punched into 96-well plates.
- Based on the results of the real-time and transport simulation stability study a shelf life 180 days (6 months) is claimed for all Eonis SCID-SMA kit components, when stored at the labeled storage temperature of -30°C - -16°C.

## 6. Detection Limit:

### (i) *Limit of Blank (LoB)*

The LoB is the highest measurement result that is likely to be observed for a blank sample that contains no analyte. The LoB study was performed using contrived analyte negative samples created by spiking a Coriell cell line into leukocyte depleted blood which was TREC and KREC negative. Five (5) individual contrived analyte-negative samples were prepared by spiking the Coriell cell line into 5 different lots of leukocyte-depleted human blood with the hematocrit adjusted to 40-55%. The study was performed by 1 operator using 2 sets of Instrument Systems and 2 kit lots run over 5 days for a total of 10 runs. Six replicates were tested in each run totaling 300 results (150 results/kit lot). For the RPP30 LoB, blank samples were used (no punch/no template DNA), 20 replicates were tested in each run totaling 200 results (100 results/kit lot). Each plate contained kit controls NTC, C1, C2 and C3, assayed in duplicate. Acceptance criteria were that the LoB of RPP30 must be zero; or any samples that have RPP30 Ct  $\geq$  32.00 will be reported as 'Invalid' for TREC and KREC. The LoB of TREC and KREC must be no greater than the LoD.

For TREC and KREC, all of the replicates had no Ct value reported. For RPP30, there were two replicates having Ct values reported at 34.49 and 36.23; however, as this would lead to an invalid result for the sample, the percentage of false-positive results for RPP30, defined as the percentage of replicates of the blank sample that had a valid RPP30 Ct value  $\leq$  32.00 was also zero for both kit lots. The LoB of TREC is 0 copies/ $\mu$ L blood, and the LoB of KREC is 0 copies/  $\mu$ L blood.

### (ii) *Limit of Detection (LoD) and Limit of Quantitation (LoQ)*

LoD and LoQ were evaluated separately for TREC and KREC. The LoD and LoQ studies were performed using contrived samples prepared by diluting cord blood with adult leucocyte-depleted human red blood that had TREC and/or KREC levels that were undetectable by the Eonis SCID-SMA Kit. HL-60 cells were spiked into the adult leucocyte-depleted blood to adjust the RPP30 Ct values to match the RPP30 Ct value in cord blood (within 1 Ct difference). Each cord blood sample was diluted into different levels in which at least the lowest three levels had yield hit rates within the range of 0.10 to 0.90 using the corresponding diluents. The samples were then spotted onto filter paper and dried overnight. For each analyte, TREC or KREC, the study was performed using 2 sets of instruments and 2 kit lots. One operator performed 2 runs per day, for 4 days, totaling 8 runs.

For each sample and each dilution, including negative samples (diluents), 5 replicates were tested in each run, totaling 960 results with 20 replicates per dilution per lot. The final dilutions used in the study and the corresponding number of positive results observed, total

number of measurements and hit rates observed for TREC and KREC are summarized in Tables 19 and 20, respectively.

**Table 19. TREC calculated hit rate percentages with sample concentrations in copies/10<sup>5</sup> cells**

Kit Lot	Cord Blood Sample	Sample Level	Dilution	Concentration in copies/10 <sup>5</sup> cells	Observed Positives	Observed Negatives	Hit Rate %	
1	L1	Cord Blood	1.00	973	0	20	100%	
		Level 1	2.00	528	0	20	100%	
		Level 2	3.33	428	0	20	100%	
		<b>Level 3</b>	<b>5.00</b>	<b>137</b>	<b>0</b>	<b>20</b>	<b>100%</b>	
		Level 4	10.0	77	5	15	75%	
		Level 5	12.5	44	6	14	70%	
		Level 6	16.7	44	9	11	55%	
		Level 7	25.0	36	12	8	40%	
		Diluent	N/A	21	18	2	10%	
	L2	Cord Blood	1.00	2670	0	20	100%	
		<b>Level 1</b>	<b>6.67</b>	<b>237</b>	<b>0</b>	<b>20</b>	<b>100%</b>	
		Level 2	10.0	166	2	18	90%	
		Level 3	12.5	122	2	18	90%	
		Level 4	25.0	50	11	9	45%	
		Level 5	33.3	52	10	10	50%	
		Level 6	50.0	38	13	7	35%	
		Level 7	100	15	17	3	15%	
		Diluent	N/A	0	20	0	0%	
	L3	Cord Blood	1.00	1525	0	20	100%	
		<b>Level 1</b>	<b>3.33</b>	<b>365</b>	<b>0</b>	<b>20</b>	<b>100%</b>	
		Level 2	6.67	194	2	18	90%	
		Level 3	12.5	83	3	17	85%	
		Level 4	16.7	44	7	13	65%	
		Level 5	25.0	39	5	15	75%	
		Level 6	33.3	55	9	11	55%	
		Level 7	50.0	17	14	6	30%	
		Level 8	100	22	16	4	20%	
		Diluent	N/A	0	20	0	0%	
	2	L1	Cord Blood	1.00	1141	0	20	100%
			Level 1	2.00	395	0	20	100%

Kit Lot	Cord Blood Sample	Sample Level	Dilution	Concentration in copies/10 <sup>5</sup> cells	Observed Positives	Observed Negatives	Hit Rate %
		Level 2	3.33	596	0	20	100%
		<b>Level 3</b>	<b>5.00</b>	<b>106</b>	<b>1</b>	<b>19</b>	<b>95%</b>
		Level 4	10.0	63	7	13	65%
		Level 5	12.5	84	16	4	20%
		Level 6	16.7	32	9	11	55%
		Level 7	25.0	17	17	3	15%
		Diluent	N/A	45	19	1	5%
	<b>L2</b>	Cord Blood	1.00	2490	0	20	100%
		<b>Level 1</b>	<b>6.67</b>	<b>240</b>	<b>0</b>	<b>20</b>	100%
		Level 2	10.0	118	3	17	85%
		Level 3	12.5	109	2	18	90%
		Level 4	25.0	54	8	12	60%
		Level 5	33.3	30	6	14	70%
		Level 6	50.0	25	11	9	45%
		Level 7	100	29	17	3	15%
		Diluent	N/A	0	20	0	0%
	<b>L3</b>	<b>Cord Blood</b>	<b>1.00</b>	<b>1510</b>	<b>0</b>	<b>20</b>	<b>100%</b>
		Level 1	3.33	321	2	18	90%
		Level 2	6.67	136	0	20	100%
		<b>Level 3</b>	<b>12.5</b>	<b>72</b>	<b>1</b>	<b>19</b>	<b>95%</b>
		Level 4	16.7	59	6	14	70%
		Level 5	25.0	34	11	9	45%
		Level 6	33.3	19	13	7	35%
		Level 7	50.0	44	15	5	25%
		Level 8	100	34	17	3	15%
		Diluent	N/A	0	20	0	0%

**Table 20. KREC calculated hit rate percentages with sample concentrations in copies/10<sup>5</sup> cells**

Kit Lot	Cord Blood Sample	Sample Level	Dilution	Concentration in copies/10 <sup>5</sup> cells	Observed Positives	Observed Negatives	Hit Rate %
<b>1</b>	<b>L1</b>	Cord Blood	1.00	3905	0	20	100%
		Level 1	1.43	3072	0	20	100%
		Level 2	1.67	2644	0	20	100%

Kit Lot	Cord Blood Sample	Sample Level	Dilution	Concentration in copies/10 <sup>5</sup> cells	Observed Positives	Observed Negatives	Hit Rate %	
		Level 3	2.00	2080	0	20	100%	
		<b>Level 4</b>	<b>5.00</b>	<b>1188</b>	<b>0</b>	<b>19</b>	<b>100%</b>	
		Level 5	12.5	255	4	16	80%	
		Level 6	33.3	129	7	13	65%	
		Level 7	100	117	17	3	15%	
		Diluent	N/A	0	20	0	0%	
	<b>L2</b>	Cord Blood	1.00	6438	0	20	100%	
		Level 1	6.67	871	0	20	100%	
		Level 2	10.0	757	0	19	100%	
		<b>Level 3</b>	<b>12.5</b>	<b>478</b>	<b>0</b>	<b>20</b>	<b>100%</b>	
		Level 4	25.0	305	6	14	70%	
		Level 5	33.3	226	7	13	65%	
		Level 6	50.0	151	7	13	65%	
		Level 7	100	147	11	9	45%	
		Diluent	N/A	0	20	0	0%	
	<b>L3</b>	Cord Blood	1.00	5064	0	20	100%	
		Level 1	3.33	1588	0	20	100%	
		<b>Level 2</b>	<b>6.67</b>	<b>804</b>	<b>0</b>	<b>20</b>	<b>100%</b>	
		Level 3	12.5	296	2	18	90%	
		Level 4	16.7	187	2	18	90%	
		Level 5	25.0	166	3	17	85%	
		Level 6	33.3	178	5	14	74%	
		Level 7	50.0	200	16	4	20%	
		Level 8	100	174	11	9	45%	
		Diluent	N/A	0	20	0	0%	
	<b>2</b>	<b>L1</b>	Cord Blood	1.00	4064	0	20	100%
			Level 1	1.43	4146	0	20	100%
			Level 2	1.67	3828	0	20	100%
			Level 3	2.00	2039	0	20	100%
			Level 4	<b>5.00</b>	<b>846</b>	<b>0</b>	<b>20</b>	<b>100%</b>
Level 5			12.5	209	2	18	90%	
Level 6			33.3	191	6	14	70%	
Level 7			100	128	15	5	25%	
Diluent			N/A	0	20	0	0%	
<b>L2</b>		Cord Blood	1.00	5825	0	20	100%	

Kit Lot	Cord Blood Sample	Sample Level	Dilution	Concentration in copies/10 <sup>5</sup> cells	Observed Positives	Observed Negatives	Hit Rate %
		Level 1	6.67	1064	0	20	100%
		Level 2	10.0	556	0	20	100%
		Level 3	12.5	412	0	20	100%
		<b>Level 4</b>	<b>25.0</b>	<b>235</b>	<b>1</b>	<b>19</b>	<b>95%</b>
		Level 5	33.3	245	3	17	85%
		Level 6	50.0	116	7	13	65%
		Level 7	100	134	10	10	50%
		Diluent	N/A	0	20	0	0%
	<b>L3</b>	Cord Blood	1.00	4915	0	20	100%
		Level 1	3.33	1510	0	20	100%
		<b>Level 2</b>	<b>6.67</b>	<b>659</b>	<b>0</b>	<b>20</b>	<b>100%</b>
		Level 3	12.5	311	2	18	90%
		Level 4	16.7	287	3	17	85%
		Level 5	25.0	296	3	17	85%
		Level 6	33.3	167	5	15	75%
		Level 7	50.0	158	9	11	55%
		Level 8	100	143	15	5	25%
		Diluent	N/A	0	20	0	0%

The LoD was calculated using Probit analysis at 95% probability for each reagent kit lot. Probit analysis in copies/10<sup>5</sup> cells unit was performed. The data collected from all 3 cord blood samples and their dilutions were pooled for the calculation. The results for both kit lots with 95% confidence intervals for TREC and KREC are summarized in Tables 21 and 22, respectively. The LoD estimate is the concentration value that Probit analysis gives an estimated hit rate of 95%.

**Table 21. TREC LoD estimates based on Probit analysis with hit rate percentages and sample concentrations in copies/10<sup>5</sup> cells unit**

Kit lot	LoD estimate (copies/10 <sup>5</sup> cells)	95% confidence interval
1	174	137 - 242
2	242	184 - 360

**Table 22. KREC LoD estimates based on Probit analysis with hit rate percentages and sample concentrations in copies/10<sup>5</sup> cells unit**

Kit lot	LoD estimate (copies/10 <sup>5</sup> cells)	95% confidence interval
1	459	371 - 669
2	340	292 - 439

The LoQ was evaluated as the functional sensitivity that represents the measurand concentration associated with a desired within-laboratory precision, as exemplified by CLSI guideline EP17-A2 Appendix D. Only the dilutions that yielded 100% hit rates were included in the LoQ calculation. For each qualified dilution, the mean and the SD of the concentrations in Ln (copies/10<sup>5</sup> cells) were calculated. A power function model (SD vs. the mean concentration) was used to fit the datasets for each reagent kit lot. The LoQ estimate for each reagent kit lot was determined as the predicted lowest concentration that has within-laboratory precision equal to an SD ≤ 0.90 Ln (copies/10<sup>5</sup> cells) for TREC and ≤ 1.49 Ln (copies/10<sup>5</sup> cells) for KREC. The results are shown in Tables 23 and 24.

**Table 23. TREC LoQ estimates based on conservative approach in copies/10<sup>5</sup> cells**

Kit lot	LoD (copies/10 <sup>5</sup> cells)	LoQ estimate (copies/10 <sup>5</sup> cells)	Final LoQ (=LoD) (copies/10 <sup>5</sup> cells)
1	174	30	242
2	242	123	

**Table 24. KREC LoQ estimates based on conservative approach in copies/10<sup>5</sup> cells**

Kit lot	LoD (copies/10 <sup>5</sup> cells)	LoQ estimate (copies/10 <sup>5</sup> cells)	Final LoQ (=LoD) (copies/10 <sup>5</sup> cells)
1	459	7	459
2	340	0	

The following results were determined for LoD and LoQ of TREC and KREC:

LoD = LoQ  
TREC: 242 copies/10<sup>5</sup> cells  
KREC: 459 copies/10<sup>5</sup> cells

7. Assay Cut-Off:

Refer to the Clinical Cut-off in Section VII.D.

8. Accuracy (Instrument):

Not applicable

9. Carry-Over:

A study was conducted to evaluate the potential for sample carryover to demonstrate control of contamination during testing of samples with the Eonis reagent kits. The carryover test was performed using 1 kit lot. The high-analyte positive sample was prepared by spiking plasmid DNAs of each analyte into leukocyte-depleted human blood at levels lower than 29.71 for TREC Ct, 29.22 for KREC Ct, and 22.31 for RPP30 Ct, corresponding to the 10% quantiles of a population of newborn specimens. The high analyte-positive sample and the analyte-negative sample were placed in alternating rows and columns to produce a checkerboard configuration. The observed carryover percentage was 0% for the three analytes.

**B Comparison Studies:**

1. Method Comparison with Predicate Device:

Not applicable

2. Matrix Comparison:

*Filter Paper Reproducibility*

A filter paper reproducibility study was conducted using filter paper from 2 manufacturers in order to evaluate manufacturer difference, lot differences, as well as to assess homogeneity across dried blood spots (DBS). Between paper manufacturer and paper lots were evaluated using samples spotted on 3 lots of each filter paper brand.

The study was conducted using a subset of the samples that were used in the precision and 3-site reproducibility studies (although some at different concentrations) covering high, normal and low levels of each analyte (Table 5).

- Sample 1 High TREC: Prepared by spiking TREC plasmid into cord blood at 2000 copies/ $\mu$ L with hematocrit adjusted to 51.8 (approx. Ct value 28.19)
- Sample 2 High KREC: Prepared by spiking TREC plasmid into cord blood at 2000 copies/ $\mu$ L with hematocrit adjusted to 51.8 (approx. Ct value 28.65)
- Sample 3: Prepared by spotting cord blood with hematocrit adjusted to 51.8 onto 3 lots of each filter paper brand
- Sample 11: Prepared by spiking TREC plasmid at 700 copies/ $\mu$ L, KREC plasmid at 1000 copies/ $\mu$ L and Coriell cells at 15000 cells/ $\mu$ L into leukocyte depleted blood with hematocrit adjusted to 49.
- Sample 12 Low levels (above cut-off) TREC and KREC: Prepared by diluting 400  $\mu$ L of cord blood with hematocrit adjusted to 50.9 in 600  $\mu$ L of adult whole blood with hematocrit adjusted to 52.9 (approx. TREC Ct value 33.29 and approx. KREC Ct value 32.20)
- Sample 13 Low levels (below cut-off) TREC and KREC: Prepared by spotting adult whole blood with hematocrit adjusted to 52.9 onto 3 lots of each paper brand (approx. TREC Ct value no Ct and approx. KREC Ct value 3.44)



Each analyte was manually spotted on 3 filter paper lots at 85 µL per spot.

One operator performed 1 run per day for 5 days for a total of 5 runs and using 1 kit lot and 1 Eonis SCID-SMA test system. Six samples were prepared on 3 lots of each filter paper brand for a total of 36 conditions (6 samples x 3 lots x 2 filter paper brands = 36). Each run consisted of 3 96-well extraction plates (2 randomized sample sets per plate plus 2 sets of kit controls) and 1 384-well PCR plate in each testing day.

To assess the homogeneity of the analytes across the DBS, 5 replicates of each condition were punched from 5 different locations of each spot in the following order: top, left, center, right, bottom. The study generated 900 results (25 replicates x 6 samples x 3 filter paper lots x 2 filter paper brands) with 150 results per sample and 450 results per filter paper brand. The total variation lnSD results are compared to the product specification for total variation (“Acceptance criteria”) and are shown in Table 25. The total variation lognormal CV% results are presented in Table 26.

**Table 25. TREC and KREC precision analysis with between variation estimates across paper manufacturers and across paper lots**

Analyte	Sample	N	Mean copies/ 10 <sup>5</sup> cells	Mean Ln copies/10 <sup>5</sup> cells	Within Paper Lot lnSD	Between Paper Lot lnSD	Between Manufacturer lnSD	Total lnSD	Specifi- cation
TREC	2	147	584	6.37	0.61	0.00	0.02	0.61	0.90
	3	149	871	6.77	0.49	0.04	0.09	0.50	0.90
	12	150	1300	7.17	0.28	0.03	0.01	0.28	0.90
	11	150	4105	8.32	0.23	0.06	0.05	0.24	0.90
	1	150	14328	9.57	0.25	0.00	0.03	0.25	0.90
KREC	13	150	1064	6.97	0.43	0.00	0.02	0.43	1.49
	12	150	4403	8.39	0.24	0.03	0.03	0.25	1.49
	1	150	4817	8.48	0.26	0.03	0.01	0.26	1.49
	3	150	5486	8.61	0.21	0.01	0.04	0.22	1.49
	11	150	11968	9.39	0.23	0.06	0.03	0.24	1.49
	2	147	22026	10.0	0.27	0.00	0.08	0.28	1.49

**Table 26. TREC and KREC precision analysis with between variation estimates across paper manufacturers**

Analyte	Sample	N	Mean copies/10 <sup>5</sup> cells	Within Paper Lot Lognormal CV%	Between Paper Lognormal CV%	Between Manufacturer Lognormal CV%	Total Lognormal CV%
TREC	2	147	584	67	0	6	67
	3	149	871	52	4	2	53
	12	150	1300	28	3	9	29

Analyte	Sample	N	Mean copies/10 <sup>5</sup> cells	Within Paper Lot Lognormal CV%	Between Paper Lognormal CV%	Between Manufacturer Lognormal CV%	Total Lognormal CV%
	11	150	4105	24	6	1	25
	1	150	14328	25	0	5	25
KREC	13	150	1064	45	0	3	45
	12	150	4403	25	3	2	25
	1	150	4817	26	3	3	26
	3	150	5486	21	1	1	22
	11	150	11968	23	6	4	24
	2	147	22026	27	0	3	28

ANOVA results in Table 27 below demonstrate that there were statistically significant differences between paper manufacturers (Sample #2/KREC) and between paper lots (Sample #11/TREC and KREC); however, the highest between lot variation was 0.06 (Sample #11/TREC and KREC), and the highest between manufacturer variation was 0.09 (Sample #2/KREC). Per study acceptance criteria, the variation between these samples would not be clinically significant. The data demonstrated the similarity in performance between two major brands of filter paper.

**Table 27. TREC and KREC between-paper lot variation and between-paper manufacturer variation, and the corresponding p-values from ANOVA**

Analyte	Sample	N	Mean Copies (copies/10 <sup>5</sup> cells)	Mean Ln Copies	Between Paper Lot SD	P-value for testing level differences between paper lots	Between Manufacturer SD	P-value for testing level differences between manufacturers
TREC	2	147	584	6.37	0.00	0.781	0.02	0.759
	3	149	871	6.77	0.04	0.272	0.09	0.133
	12	150	1300	7.17	0.03	0.226	0.01	0.725
	11	150	4105	8.32	0.06	0.011	0.05	0.091
	1	150	14328	9.57	0.00	0.776	0.03	0.291
KREC	13	150	1064	6.97	0.00	0.206	0.02	0.652
	12	150	4403	8.39	0.03	0.736	0.03	0.365
	1	150	4817	8.48	0.03	0.063	0.01	0.622
	3	150	5486	8.61	0.01	0.262	0.04	0.143
	11	150	11968	9.39	0.06	0.009	0.03	0.291
	2	147	22026	10.0	0.00	0.726	0.08	0.017

## C Clinical Studies:

### (i) *Specimens*

A clinical study was conducted at one foreign testing site, the Statens Serum Institute (SSI) in Denmark. SSI is the same site that conducted the study for the PerkinElmer EnLite Neonatal TREC Assay (DEN140010). The majority of the DBS specimens tested for the clinical study (retrospective, anonymized residual [leftover] newborn specimens) were procured by SSI from the screened Danish population. Specimen collection was prospective-retrospective and enriched with confirmed positive specimens for SCID and XLA which came from two biobanks: Danish Neonatal Screening Biobank and the California Biobank. Additionally, 3 of the 4 confirmed positive samples for XLA were procured from the SSI. Confirmed positive specimens were also procured from the California Department of Public Health Genetic Disease Screening Program Biobank (CA), including SCID (all 17 specimens) and XLA (1 of 4 confirmed specimens). All samples were tested at SSI as part of the evaluation. Age matched control specimens were used and presumed to be normal and were archived just before and after the time period during which the confirmed positive specimens for each of the two disorders (SCID, XLA) were archived. These specimens were used as controls for storage conditions to ensure specimen integrity remained intact within the timeframe during which they were stored. Therefore, in addition to consecutive specimens, archived specimens with storage time matched to the storage time of positive specimens were included. Specimen inclusion/exclusion criteria were provided and included specimen age, and storage conditions.

### (ii) *Inclusion/Exclusion criteria for assessment of specimens*

The age of the newborns at the time of specimen collection was less than or equal to 7 days to mimic the expected intended use population and a typical sample collection age for first tier screening to enable early diagnosis. Routine newborn screening specimens originating from the study site were collected at least 5.5 years ago to be able to assess, based on medical records, that the subject had not been identified to present clinical signs attributable to SCID or XLA. The confirmed positive SCID or XLA specimens had been collected  $\leq 17$  years before the start date of the study to avoid inclusion of sample stored for longer than that time period. Confirmed positive XLA samples were difficult to procure; therefore, samples violating the upper limit of the storage time or the limit for the age of the newborn at the time of sample collection were included.

Inclusion criteria included the following:

- Age of the newborn at the time of specimen collection  $\leq 7$  days
- Specimen storage conditions are known
- Specimen collected at least 5.5 years ago (from study start)

For specimens confirmed positive for SCID:

- Specimen confirmed positive for typical SCID (documentation of confirmed diagnosis and/or documentation of flow cytometry result)
- Date of specimen collection is  $\leq 17$  years from the start date of the study

For specimens confirmed positive for XLA:

- Specimen confirmed positive for XLA (e.g., documentation of BTK gene mutation and/or flow cytometry result)
- Date of specimen collection is  $\leq 17$  years from the start date of the study

In addition, for storage age-matched control specimens:

- Date of specimen collection is known
- Specimen collected just before or after the time period during which the confirmed positive specimens for SCID and/or XLA was collected

Exclusion criteria included the following:

- The DBS is not of appropriate quality (i.e., DBS is scratched, abraded, show signs of incomplete saturation or supersaturation, dilution, discoloration, or contamination)
- Parents/guardians have opted out for use of the specimen beyond routine screening
- Repeat specimen from a newborn whose specimen is already included in the study
- Less than one disk with a diameter of 3.2 mm (1/8 inch) of DBS available in the cut-off phase, and less than two disks available in the pivotal phase for routine specimens

In addition, for specimens confirmed positive for SCID:

- Other than typical SCID, such as less acute SCID (e.g., leaky SCID, variant SCID) or SCID like syndromes (e.g., DiGeorge syndrome or Omenn syndrome)

In addition, for specimens confirmed positive for XLA:

- Atypical XLA or female carriers

(iii) *Study Design:*

A total of 6924 unique specimens were tested in the study: 3432 specimens were tested in the cut-off phase and 3484 specimens were tested in the pivotal phase (this number includes 17 SCID cases, 6 XLA cases and 52 SMA cases; SMA cases are deducted from the table – refer to DEN200044). The main study consisted of two phases. Specimens in the cut-off study were not retested in the pivotal study.

Cut-off Phase: Newborn specimens representative of samples of the target population at the study site were tested with the device to produce newborn population distribution data for TREC and KREC and to select clinical cut-off values for the pivotal study. The number of false positives for TREC and KREC were determined from the data using lower percentiles 0.3%, 0.5, and 1% and are shown in Table 28.

**Table 28. Descriptive statistics for the copies (per  $10^5$  cells) for KREC and TREC and the percentile cut-offs used in the study**

Analyte	Copies / $10^5$ cells					
	Median	Min	Max	0.3 percentile	0.5 percentile	1.0 percentile
TREC	2520	117	9990	262	413	563
KREC	3590	12	20100	129	261	484

**Pivotal Phase:** Presumed unaffected newborn specimens representing the target population at the study site were supplemented with confirmed positive specimens collected from SCID and 6 XLA affected newborns. A total of 3432 unique samples as shown in Table 29 (after subtracting 52 SMA cases – refer to DEN200044), were assessed in the pivotal phase. After exclusion of the samples not fulfilling the inclusion/exclusion criteria, a total of 3383 samples were included in the final analysis. Demographic characteristics of the routine screening specimens and confirmed positive samples were collected. These included: age at the time of sample collection (mean 2.5 days); weight distribution of the newborns at the time of specimen collection (95% ≥ 2500g); gestational age (98% ≥ 34 weeks); and gender (50.6% male and 49.4% female).

**Table 29. Summary of Pivotal Phase Specimens Assayed and Included in the Data Analysis**

Sample	Assayed Samples	Excluded due to failure to meet inclusion/exclusion criteria	Excluded due to Invalid Result	Samples included in final analysis
Routine screening specimens	3044	26	-	3018
SCID positives	17		-	17
XLA positives	6		-	6
Storage time matched controls	365	23	-	342
<b>Total</b>	<b>3432</b>	<b>49</b>	<b>0</b>	<b>3383</b>

Confirmatory test results were used as the comparator for the confirmed positive SCID or XLA cases including the presence of BTK mutation for XLA. Determination of the true clinical status with regard to SCID has been previously described (DEN140010). The true clinical status with regard to XLA was assessed using data available in patient registries. The clinical assessment of the study subjects was obtained from their medical records to confirm that the subject, at 5.5 years of age (defined as 66 months) or younger, has not been identified with XLA (has no clinical signs attributable to XLA) and is apparently healthy. Danish hospital registry data were reviewed by filtering out subjects that had at least one of the 275 ICD10 diagnosis codes that were listed as potentially linked to XLA. By the age of 5.5 years, a subject with a XLA would have either:

- Been diagnosed with the XLA; or
- Presented with serious medical complications suggestive of XLA as listed below:

**Table 30. Symptoms and signs associated with XLA**

Age of onset	Symptoms, signs
0-2 years	A severe life-threatening bacterial infection such as sepsis, meningitis, cellulitis, or empyema

Age of onset	Symptoms, signs
	Paucity of lymphoid tissue (small adenoids, tonsils, and lymph nodes on physical examination)
< 5 years	Recurrent otitis, pneumonitis, sinusitis, and conjunctivitis

To establish the screening performance, the cut-offs for TREC and KREC were determined by calculating the TREC and KREC concentrations corresponding to the 0.3th and 1.0th population percentiles (262 copies/10<sup>5</sup> cells for TREC, 484 copies/10<sup>5</sup> cells for KREC) established in the cut-off study with an independent dataset and were applied in the pivotal study

Testing was in singlicate unless repeat testing was necessary due to an invalid result. The specimens having TREC and KREC levels below the cut-off values (i.e., potential positives) in the initial round of testing were re-tested in duplicate. Any invalids were also retested in singlicate per protocol. If either replicate was positive, the sample was considered positive. All 3 replicates were required to be negative for the sample to be negative. The final results (presumptive normal, presumptive positive, invalid result) were classified after the second round of testing according to the testing interpretation algorithm.

The Eonis SCID-SMA kits' performance with each subjects' sample was compared to the clinical data/outcome captured for that subject. Status of SCID patients were determined previously (DEN140010). A total of 18 subjects had screen positive results with initial KREC results below the cut-off (suggestive for XLA). After replicate testing, 9 subjects with repeat test KREC results below cut-off (one or both replicates below cut-off) remained. For all 18 subjects, the status was alive at 5.5 years after birth, and 6 subjects had at least one ICD10 code listed in the registry. Of these 6 subjects, none had multiple diagnosis that were indicative of XLA, nor did they have diagnosis that mandated further investigation. Thus, all were concluded to be unaffected for XLA (Table 31).

**Table 31. Summary of the routine screening samples tested in the pivotal study**

Routine screening samples	TREC	KREC
Screened samples	3018	3018
Initial screen positive	10	18
Initial screen negative	3008	3000
Retest rate	0.3%	0.6%
Final screen positive	9	9
Final screen negative	3008	3004
False-positive rate	0.3%	0.3%

The retest rate for the routine screening population was 0.9 % with the study cut-offs for TREC and KREC. The results of the screening performance after retesting per protocol for SCID (TREC) below in Tables 32, along with the false positive rate and 95% confidence intervals (CI) in Table 33.

**Table 32. Screening performance of TREC with the Eonis SCID-SMA kit**

TREC		Clinical Status		Total (%)
		SCID affected (%)	Normal (%)	
Screening result	Presumptive positive (%)	17 (100 %)	9 (0.3 %)	26 (0.9 %)
	Presumptive normal (%)	0 (0.0 %)	3008 (99.7 %)	3008 (99.1 %)
	<b>Total (%)</b>	17 (100 %)	3017 (100 %)	3034 (100 %)

**Table 33. SCID false positive rate observed in the pivotal trial after retest**

	Positive %-agreement	Negative %-agreement	Overall %-agreement	False positive rate	False negative rate	Positive predictive value	Negative predictive value
<b>Percent</b>	100 %	99.7 %	99.7 %	0.3 %	0.0 %	0.6 %	100 %
<b>Confidence interval</b>	80.5 % - 100 %	99.4 % - 99.9 %	99.4 % - 99.9 %	0.1 % - 0.6 %	0 % - 19.5 %	NA	NA

The results of the screening performance after retesting per protocol for XLA (KREC) is shown in Tables 34 and 35, along with the false positive rate after retest and 95% confidence intervals (CI) in Table 36.

**Table 34. Screening performance of the Eonis SCID-SMA device for XLA after retesting. The screening performance is presented with 95% confidence intervals.**

KREC		Clinical Status		Total (%)
		XLA affected (%)	Normal (%)	
Screening result	Presumptive positive (%)	6 (100 %)	9* (0.3 %)	15 (0.5 %)
	Presumptive normal (%)	0 (0.0 %)	3004 (99.7 %)	3004 (99.5 %)
	<b>Total (%)</b>	6 (100 %)	3013 (100 %)	3019 (100 %)

\*Includes one false positive result from a female subject. In female population the false positive rate is estimated to be 0.1 %.

Because XLA is a disease that affects only males, the results for KREC are shown for males in Table 35.

**Table 35. Screening performance of the Eonis SCID-SMA device for XLA after retesting in males alone**

KREC (Male only)		Clinical Status		Total (%)
		XLA affected (%)	Normal (%)	
Screening result	Presumptive positive (%)	6 (100%)	8 (0.5%)	14 (0.9%)
	Presumptive normal (%)	0 (0.0%)	1515 (99.5%)	1515 (99.1%)
	Total (%)	6 (100%)	1523 (100%)	1529 (100%)

**Table 36. XLA False positive rate observed in the trial after retest**

	Positive %-agreement	Negative %-agreement	Overall %-agreement	False positive rate	False negative rate
Percent	100 %	99.4 %	99.4 %	0.6 %	0.0 %
Confidence interval	39.8 % - 100 %	99.1 % - 99.6 %	99.1 % - 99.6 %	0.4 % - 0.9 %	0 % - 60.2 %

**D Clinical Cut-Off:**

TREC: 262 copies/10<sup>5</sup> cells  
 KREC: 484 copies/10<sup>5</sup> cells

**E Expected Values/Reference Range:**

The incidence of SCID is approximately 1 in 50,000 live births.<sup>1</sup> The incidence of XLA is approximately 1 in 190,000 live male births.<sup>2</sup> The distribution of a newborn population was determined by analyzing 3341 newborn DBS from the Danish Newborn Screening Biobank (NBS-Biobank) with the Eonis SCID-SMA kit. The specimens were archived routine screening specimens. The TREC and KREC descriptive statistics (median and range of results observed) and the value at which 0.3, 0.5, or 1.0 percent of the specimens were below of the specimens are displayed in Table 28. Data separated by gender for KREC is shown in Table 37.

**Table 37. Descriptive statistics and percentiles for KREC copies/10<sup>5</sup> cells by gender**


Analyte	Gender	N	Copies/10 <sup>5</sup> cells					
			Median	Min	Max	0.3 percentile	0.5 percentile	1.0 percentile
KREC	Male	1699	3460	73	16200	78	180	433
	Female	1642	3740	12	20100	232	363	576

<sup>1</sup> [6] Puck JM. (2011): Neonatal screening for severe combined immunodeficiency. Curr Opin Pediatr. 2011 Dec;23(6):667-73.


<sup>2</sup> [18] Winkelstein et al., X-Linked Agammaglobulinemia, Report on a United States Registry of 201 Patients, Medicine, 2006



The frequency distributions (N=3341) of TREC and KREC analyte are shown in Figures 6 and 7.



**Figure 6. TREC (copies/10<sup>5</sup> cells) distribution of presumed normal newborn samples**



**Figure 7. KREC (copies/10<sup>5</sup> cells) distribution of presumed normal newborn samples.**

**F Other Supportive Instrument Performance Characteristics Data:**

DNA Extraction Equivalency

This study was conducted to demonstrate that Eonis test system has equivalent performance when the DNA extraction is performed either manually or by a different automated liquid handler and is designed based on the recommendation of CLSI Guideline EP05-A.

Two operators performed 3 Eonis extractions and PCR consolidations/day for 5 days using 1 kit lot, totaling 30 runs using 7 samples from the precision study (Table 5). Each day, an operator would run 1 extraction/PCR consolidation on the same JANUS Extraction/JANUS PCR Mastermix, 1 extraction/PCR consolidation on a second commercial liquid handler, and 1 manual extraction manual PCR setup/consolidation. The study generated a total of 1050 results (5 replicates x 7 samples x 3 extraction/PCR methods x 2 operators x 5 days) with 150 results per sample and 350 results per extraction/PCR method.

There were 45 replicates performed with the JANUS handler, 50 replicates for the second commercial liquid handler, and 50 replicates for the manual extraction process. The results for TREC and KREC, evaluated using variance analysis method, are presented in Table 38 and Table 39.

**Table 38. DNA extraction equivalency for TREC**

Sample	N	Mean (Copies/10 <sup>5</sup> cells)	Mean Ln Copies	Repeatability Ln SD	Between Extraction Method Ln SD	Total Ln SD
13	56	36	3.57	0.90	0.02	0.95
4	145	781	6.66	0.35	0.15	0.39
12	145	992	6.90	0.30	0.09	0.32
3	145	1188	7.08	0.30	0.07	0.31
2	145	1920	7.56	0.26	0.19	0.34
11	145	4964	8.51	0.19	0.13	0.27
1	145	13494	9.51	0.15	0.23	0.30

**Table 39. DNA extraction equivalency for KREC**

Sample	N	Mean (Copies/10 <sup>5</sup> cells)	Mean Ln Copies	Repeatability Ln SD	Between Extraction Method Ln SD	Total Ln SD
13	145	1097	7.00	0.41	0.04	0.42
12	145	3429	8.14	0.23	0.11	0.26
3	145	5768	8.66	0.26	0.14	0.29
1	145	6438	8.77	0.18	0.15	0.24
11	145	14618	9.59	0.19	0.14	0.25
4	145	17327	9.76	0.19	0.21	0.29
2	145	54176	10.9	0.17	0.26	0.34

**VIII Proposed Labeling:**

The labeling supports the finding of substantial equivalence for this device.

**IX Conclusion:**

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.