

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION  
DECISION SUMMARY**

**A. 510(k) Number:**

K183597

**B. Purpose for Submission:**

To establish substantial equivalence to a predicate device and to obtain clearance for a new assay: the QIAstat-Dx Respiratory Panel.

**C. Measurands:**

Influenza A Matrix gene (M)  
Influenza A H1 Hemagglutinin gene (HA)  
Influenza A H3 Hemagglutinin gene (HA)  
Influenza A H1 pdm09 Neuraminidase gene (NA)  
Influenza B Nucleoprotein gene (NP)  
Coronavirus 229E Membrane protein gene (M)  
Coronavirus OC43 Nucleocapsid protein gene (N)  
Coronavirus NL63 Nucleocapsid protein gene (N)  
Coronavirus HKU1 Nucleocapsid protein gene (N)  
Parainfluenza virus 1 Hemagglutinin-neuraminidase glycoprotein gene (HN)  
Parainfluenza virus 2 Hemagglutinin-neuraminidase glycoprotein gene (HN)  
Parainfluenza virus 3 Phosphoprotein gene (P)  
Parainfluenza virus 4 Nucleocapsid protein gene (N)  
Rhinovirus/Enterovirus 5'-UTR region  
Adenovirus Hexon gene  
Respiratory Syncytial Virus A/B Matrix protein gene (M)  
Human metapneumovirus A/B Nucleoprotein gene (N)  
*Legionella pneumophila* Macrophage infectivity potentiator gene (MIP)  
*Chlamydomphila pneumoniae* Major outer membrane protein gene (ompA)  
*Mycoplasma pneumoniae* Cytadhesin gene (P1)  
*Bordetella pertussis* Transposase Insertion element (IS481)

**D. Type of Test:**

Multiplexed Real-Time reverse transcription nucleic acid amplification assay for the amplification and detection of multiple respiratory pathogen nucleic acids.

**E. Applicant:**

QIAGEN GmbH

## **F. Proprietary and Established Names:**

QIAstat-Dx<sup>®</sup> Respiratory Panel  
QIAstat-Dx<sup>®</sup> Analyzer

## **G. Regulatory Information:**

1. Regulation section:

21 CFR 866.3980, Respiratory Viral Panel Multiplex Nucleic Acid Assay

2. Classification:

Class II

3. Product code:

OCC, OEM, OOU, OEP, OOI, OTG, OZX, OZY, OQW, OZZ

4. Panel:

Microbiology (83)

## **H. Intended Use:**

1. Intended use(s):

The QIAstat-Dx Respiratory Panel is a multiplexed nucleic acid test intended for use with QIAstat-Dx system for the simultaneous qualitative detection and identification of multiple respiratory viral and bacterial nucleic acids in nasopharyngeal swabs (NPS) eluted in Universal Transport Media (UTM) obtained from individuals suspected of respiratory tract infections. The following organism types and subtypes are identified using the QIAstat-Dx Respiratory Panel: Adenovirus, Coronavirus 229E, Coronavirus HKU1, Coronavirus NL63, Coronavirus OC43, Human Metapneumovirus A+B, Influenza A, Influenza A H1, Influenza A H3, Influenza A H1N1/pdm09, Influenza B, Parainfluenza Virus 1, Parainfluenza Virus 2, Parainfluenza Virus 3, Parainfluenza Virus 4, Rhinovirus/Enterovirus, Respiratory Syncytial Virus A+B, *Bordetella pertussis*, *Chlamydomphila pneumoniae* and *Mycoplasma pneumoniae*.

The detection and identification of specific viral and bacterial nucleic acids from individuals presenting with signs and symptoms of a respiratory infection aids in the diagnosis of respiratory infection if used in conjunction with other clinical and epidemiological information. The results of this test should not be used as the sole basis for diagnosis, treatment or other management decisions. Negative results in the setting of a respiratory illness may be due to infection with pathogens that are not detected by the test or lower respiratory tract infection that is not detected by a nasopharyngeal swab

specimen. Positive results do not rule out co-infection with other organisms: the agent(s) detected by the QIAstat-Dx Respiratory Panel may not be the definite cause of disease. Additional laboratory testing (e.g., bacterial and viral culture, immunofluorescence and radiography) may be necessary when evaluating a patient with possible respiratory tract infection.

Due to the small number of positive specimens collected for certain organisms during the prospective clinical study, performance characteristics for *Bordetella pertussis* and Parainfluenza Virus 1 were established primarily with retrospective clinical specimens. Performance characteristics for *Chlamydophila pneumoniae*, Parainfluenza Virus 2, Parainfluenza Virus 4, Influenza A subtype H1 and Coronavirus 229E were established primarily using contrived clinical specimens.

Due to the genetic similarity between Human Rhinovirus and Enterovirus, the QIAstat-Dx Respiratory Panel cannot reliably differentiate them. A positive QIAstat-Dx Respiratory Panel Rhinovirus/Enterovirus result should be followed-up using an alternate method (e.g., cell culture or sequence analysis).

Performance characteristics for Influenza A were established when Influenza A H1N1-2009 and A H3 were the predominant Influenza A viruses in circulation. Performance of detecting Influenza A may vary if other Influenza A strains are circulating or a novel Influenza A virus emerges. If infection with a novel Influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent Influenza viruses and sent to state or local health departments for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.

2. Indication(s) for use:

Same as Intended Use

3. Special conditions for use statement(s):

For Prescription Use

4. Special instrument requirements:

To be used only with the QIAstat-Dx Analyzer

**I. Device Description:**

**QIAstat-Dx Respiratory Panel Reagent Kit**

Each QIAstat-Dx Respiratory Panel reagent kit contains six (6) QIAstat-Dx Cartridges in individually wrapped foil pouches and six (6) individually packaged transfer pipettes for dispensing liquid sample into the QIAstat-Dx Cartridge.

### **QIAstat-Dx Respiratory Panel Cartridge**

The QIAstat-Dx Respiratory Panel cartridge is a disposable plastic device that allows the RP assay to be performed on the QIAstat-Dx Analyzer. The cartridge includes all the necessary reagents to perform extraction, amplification, and detection of target nucleic acids from the eluted NPS specimen. All sample preparation and assay steps are performed within the cartridge.

### **QIAstat-Dx Analyzer**

The QIAstat-Dx Analyzer 1.0 is the unit that hosts a cartridge and, on command from the user, is able to run predefined assay protocols. Within the cartridge, multiple steps are automatically performed in sequence by using pneumatic pressure and a multiport valve to transfer sample and fluids via the Transfer Chamber to their intended destinations. Following the introduction of the sample from a disposable transfer pipette, the following assay steps occur automatically and sequentially:

- Resuspension of internal control and Proteinase K enzyme;
- Cell lysis using mechanical (rotation) and chemical (chaotropic and isotonic) means;
- Mixing of the purified nucleic acid with lyophilized “Master Mix” reagents;
- Sequential transfer of mixed eluate/Master Mix from the Transfer Chamber to each Reaction Chamber containing the specified, air dried primers and probes;
- Within each Reaction Chamber, real-time, multiplex PCR (“rt-PCR”) testing is performed. Increase in fluorescence (indicative of detection of each target analyte) is detected directly within each Reaction Chamber;
- The detected signal per fluorescent marker per Reaction Chamber is then used by the system software to generate the assay result.

### **Materials Required but not Provided**

- QIAstat-Dx Analyzer Instrument
- QIAstat-Dx Analyzer User Manual

### **Interpretation of Results**

The QIAstat-Dx Analyzer automatically interprets and saves test results. After ejecting the cartridge, the results summary screen is automatically displayed. Detected analytes (i.e.; positive results) are displayed at the top of the list under the category ‘Detected’ in red font with a plus sign (+) next to the name. Equivocal results are next listed in yellow font with a question mark (?) next to the name. The last section of the results screen shows all targets tested with either a plus sign if it was detected, a question mark if the result was equivocal, and a minus sign (-) with the name in green colored font if the analyte was tested but not detected.

### **Quality Control**

The QIAstat-Dx Respiratory Panel Cartridge includes Internal Control (IC), a titered lyophilized MS2 bacteriophage that provides verification that all steps of the analysis process including sample resuspension, lysis, nucleic acid purification, reverse transcription, and PCR were successful. The results screen displays a message indicating that the internal controls “Passed” when the test was run successfully. An IC message of “Failed” indicates

that the internal control was not amplified; ‘Positive’ test results are still reported as positive, but all ‘Negative’ results are invalid. Positive and negative external controls are recommended by the manufacturer but are not provided.

**J. Substantial Equivalence Information:**

1. Predicate device name(s):

FilmArray Respiratory Panel (RP)

2. Predicate 510(k) number(s):

K123620

3. Comparison with predicate:

**Table 1: QIAstat-Dx Respiratory Panel – Comparison With Predicate**

	<b>QIAstat-Dx Respiratory Panel</b>	<b>FilmArray Respiratory Panel</b>
510(k) Number	K183597	K123620
Assay Targets	17 Respiratory virus targets plus 3 atypical bacteria	Same
Product Code	OCC	Same
Device Technology	Multiplex real time PCR	Same
Results Interpretation	Automated	Same
Time to Result	approximately 1 hour	Same
Specimen Types	Nasopharyngeal swab (NPS) eluted in UTM	NPS
Instrument	QIAstat-Dx Analyzer	FilmArray Instrument
Intended Use	The QIAstat-Dx Respiratory Panel is a multiplexed nucleic acid test intended for use with QIAstat-Dx system for the simultaneous qualitative detection and identification of multiple respiratory viral and bacterial nucleic acids in nasopharyngeal swabs (NPS) eluted in universal transport media (UTM) obtained from individuals suspected of respiratory tract infections. The following organism types and subtypes are identified using the QIAstat-Dx Respiratory Panel:	FilmArray Respiratory Panel (RP) is a multiplexed nucleic acid test intended for use with the FilmArray Instrument for the simultaneous qualitative detection and identification of multiple respiratory viral and bacterial nucleic acids in nasopharyngeal swabs (NPS) obtained from individuals suspected of respiratory tract infections. The following organism types and subtypes are identified using the FilmArray RP: Adenovirus, Coronavirus 229E, Coronavirus HKU1, Coronavirus

	<p>Adenovirus, Coronavirus 229E, Coronavirus HKU1, Coronavirus NL63, Coronavirus OC43, Human Metapneumovirus A+B, Influenza A, Influenza A H1, Influenza A H3, Influenza A H1N1/pdm09, Influenza B, Parainfluenza Virus 1, Parainfluenza Virus 2, Parainfluenza Virus 3, Parainfluenza Virus 4, Rhinovirus/Enterovirus, Respiratory Syncytial Virus A+B, <i>Bordetella pertussis</i>, <i>Chlamydophila pneumoniae</i> and <i>Mycoplasma pneumoniae</i>.</p> <p>The detection and identification of specific viral and bacterial nucleic acids from individuals presenting with signs and symptoms of a respiratory infection aids in the diagnosis of respiratory infection if used in conjunction with other clinical and epidemiological information. The results of this test should not be used as the sole basis for diagnosis, treatment or other management decisions. Negative results in the setting of a respiratory illness may be due to infection with pathogens that are not detected by the test or lower respiratory tract infection that is not detected by a nasopharyngeal swab specimen. Positive results do not rule out co-infection with other organisms: the agent(s) detected by the QIAstat-Dx Respiratory Panel may not be the definite cause of disease. Additional laboratory testing (e.g., bacterial and viral culture, immunofluorescence and radiography) may be necessary when evaluating a patient with possible respiratory tract infection.</p> <p>Due to the small number of positive specimens collected for certain organisms during the prospective clinical study, performance characteristics for <i>Bordetella pertussis</i> and Parainfluenza Virus 1 were</p>	<p>NL63, Coronavirus OC43, Human Metapneumovirus, Influenza A, Influenza A subtype H1, Influenza A subtype H3, Influenza A subtype H1-2009, Influenza B, Parainfluenza Virus 1, Parainfluenza Virus 2, Parainfluenza Virus 3, Parainfluenza Virus 4, Human Rhinovirus/Enterovirus, Respiratory Syncytial Virus, <i>Bordetella pertussis</i>, <i>Chlamydophila pneumoniae</i>, and <i>Mycoplasma pneumoniae</i>. The detection and identification of specific viral and bacterial nucleic acids from individuals exhibiting signs and symptoms of a respiratory infection aids in the diagnosis of respiratory infection if used in conjunction with other clinical and epidemiological information. The results of this test should not be used as the sole basis for diagnosis, treatment, or other management decisions. Negative results in the setting of a respiratory illness may be due to infection with pathogens that are not detected by this test or, lower respiratory tract infection that is not detected by a nasopharyngeal swab specimen. Positive results do not rule out coinfection with other organisms: the agent(s) detected by the Film Array RP may not be the definite cause of disease. Additional laboratory testing (e.g. bacterial and viral culture, immunofluorescence, and radiography) may be necessary when evaluating a patient with possible respiratory tract infection.</p> <p>Due to the small number of positive specimens collected for certain organisms during the prospective clinical study, performance characteristics for <i>Bordetella pertussis</i>, Coronavirus 229E, Coronavirus OC43, Influenza A H1, Influenza A H3, Influenza A H1-2009, Influenza B,</p>
--	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

	<p>established primarily with retrospective clinical specimens. Performance characteristics for <i>Chlamydophila pneumoniae</i>, Parainfluenza Virus 2, Parainfluenza Virus 4, Influenza A subtype H1 and Coronavirus 229E were established primarily using contrived clinical specimens.</p> <p>Due to the genetic similarity between Human Rhinovirus and Enterovirus, the QIAstat-Dx Respiratory Panel cannot reliably differentiate them. A positive QIAstat-Dx Respiratory Panel Rhinovirus/Enterovirus result should be followed-up using an alternate method (e.g., cell culture or sequence analysis).</p> <p>Performance characteristics for Influenza A were established when Influenza A H1N1-2009 and A H3 were the predominant Influenza A viruses in circulation. Performance of detecting Influenza A may vary if other Influenza A strains are circulating or a novel Influenza A virus emerges. If infection with a novel Influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent Influenza viruses and sent to state or local health departments for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.</p>	<p><i>Mycoplasma pneumoniae</i>, Parainfluenza Virus 1, Parainfluenza Virus 2, and Parainfluenza Virus 4 were established primarily with retrospective clinical specimens. Performance characteristics for <i>Chlamydophila pneumoniae</i> were established primarily using contrived clinical specimens.</p> <p>Due to the genetic similarity between Human Rhinovirus and Enterovirus, the FilmArray RP cannot reliably differentiate them. A positive FilmArray RP Rhinovirus/Enterovirus result should be followed-up using an alternate method (e.g., cell culture or sequence analysis).</p> <p>The FilmArray RP assay for Coronavirus OC43 may cross-react with some isolates of Coronavirus HKU1. A dual positive result may be due to cross-reactivity or may indicate a co-infection.</p> <p>Performance characteristics for Influenza A were established when Influenza A H1-2009, A H1, and A H3 were the predominant Influenza A viruses in circulation. Performance of detecting Influenza A may vary if other Influenza A strains are circulating or a novel Influenza A virus emerges. If infection with a novel Influenza A virus is suspected based on current clinical and epidemiological screening criteria recommended by public health authorities, specimens should be collected with appropriate infection control precautions for novel virulent Influenza viruses and sent to state or local health departments for testing. Viral culture should not be attempted in these cases unless a BSL 3+ facility is available to receive and culture specimens.</p>
--	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**K. Standard/Guidance Documents Referenced (if applicable):**

Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 1: General requirements. IEC 61010-1:2010.

Safety requirements for electrical equipment for measurement, control, and laboratory use – Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment. IEC 61010-2-2015.

Electrical equipment for measurement, control and laboratory use – EMC requirements – Part 1: General requirements. IEC 61326-1:2013

Electrical equipment for measurement, control and laboratory use. EMC requirements. Particular requirements. In vitro diagnostic (IVD) medical equipment. IEC 61326-2-6:2012

**L. Test Principle:**

Multiplexed reverse transcription nucleic acid amplification

**M. Performance Characteristics (if/when applicable):**

1. Analytical performance:

*a. Precision/Reproducibility:*

Reproducibility Study

A reproducibility study of QIAstat-Dx Respiratory Panel was conducted by operators at three sites using panels of blinded coded specimens containing high negative, low positive, and moderate positive mixed analyte samples. A total of twelve sample mixes were prepared for the study. Nine operators from three sites (five operators from site 1, two operators each from sites 2 and 3) participated in the study. The study was conducted over five days testing four replicates per day. Samples were prepared by spiking individual pathogens into HeLa cells in UTM to final concentrations of 0.1x, 1x, or 3x LOD.

The percent agreement with expected results for all analytes was  $\geq 95\%$  for samples tested at 1x and 3x LOD. All of the sample mixes generated negative test results for analytes not included in the specific mix tested. There were no significant differences observed within run (replicates tested by one operator), between runs (five different days), between sites (three sites), or between operators (nine operators).

The Reproducibility Study site-to-site qualitative results (agreements with expected results) are presented in Tables 2, 3, and 4 below.



**Table 2: QIAstat-Dx Reproducibility with Samples at 0.1x LOD**

Analyte (QIAstat-Dx Target)	Site	Positive Detected	Percent Agreement with Expected Results	95% CI
Adenovirus (Adenovirus)	Site 1	10/20	50.0%	29.9-70.1%
	Site 2	9/19	47.4%	27.3-68.3%
	Site 3	10/19	52.6%	31.7-72.7%
	All Sites	29/58	50.0%	37.5-62.5%
<i>B. pertussis</i> ( <i>B. pertussis</i> )	Site 1	9/20	45.0%	25.8-65.8%
	Site 2	7/19	36.8%	19.1-59.0%
	Site 3	9/20	45.0%	25.8-65.8%
	All Sites	25/59	42.4%	30.6-55.1%
<i>C. pneumoniae</i> ( <i>C. pneumoniae</i> )	Site 1	11/20	55.0%	34.2-74.2%
	Site 2	11/19	57.9%	36.3-76.9%
	Site 3	14/20	70.0%	48.1-85.5%
	All Sites	36/59	61.0%	48.3-72.4%
Coronavirus 229E (Coronavirus 229E)	Site 1	9/20	45.0%	25.8-65.8%
	Site 2	12/19	63.2%	41.0-80.9%
	Site 3	5/20	25.0%	11.2-46.9%
	All Sites	26/59	44.1%	32.2-56.7%
Coronavirus HKU1 (Coronavirus HKU1)	Site 1	17/20	85.0%	64.0-94.8%
	Site 2	10/19	52.6%	31.7-72.7%
	Site 3	9/20	45.0%	25.8-65.8%
	All Sites	36/59	61.0%	48.3-72.4%
Coronavirus NL63 (Coronavirus NL63)	Site 1	13/20	65.0%	43.3-81.9%
	Site 2	12/19	63.2%	41.0-80.9%
	Site 3	14/19	73.7%	51.2-88.2%
	All Sites	39/58	67.2%	54.4-77.9%
Coronavirus OC43 (Coronavirus OC43)	Site 1	13/20	65.0%	43.3-81.9%
	Site 2	15/20	75.0%	53.1-88.8%
	Site 3	15/20	75.0%	53.1-88.8%
	All Sites	43/60	71.7%	59.2-81.5%
Enterovirus (Rhinovirus / Enterovirus)	Site 1	8/20	40.0%	21.9-61.3%
	Site 2	6/19	31.6%	15.4-54.0%
	Site 3	7/20	35.0%	18.1-56.7%
	All Sites	21/59	35.6%	24.6-48.3%

Human Metapneumovirus  (hMPV)	Site 1	6/20	30.0%	14.5-51.9%
	Site 2	9/19	47.4%	27.3-68.3%
	Site 3	9/20	45.0%	25.8-65.8%
	All Sites	24/59	40.7%	29.1-53.4%
Influenza A/SwineNY/03/ 2009  (Influenza A)	Site 1	19/20	95.0%	76.4-99.1%
	Site 2	18/20	90.0%	69.9-97.2%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	57/60	95.0%	86.3-98.3%
Influenza A/Port Chalmers/1/73  (Influenza A)	Site 1	10/20	50.0%	29.9-70.1%
	Site 2	9/19	47.4%	27.3-68.3%
	Site 3	16/19	84.2%	62.4-94.5%
	All Sites	35/58	60.3%	47.5-71.9%
Influenza A/NJ/8/76  (Influenza A)	Site 1	14/20	70.0%	48.1-85.5%
	Site 2	9/19	47.4%	27.3-68.3%
	Site 3	12/20	60.0%	38.7-78.1%
	All Sites	35/59	59.3%	46.6-70.9%
Influenza A/NJ/8/76  (Influenza A H1)	Site 1	13/20	65.0%	43.3-81.9%
	Site 2	13/19	68.4%	46.0-84.6%
	Site 3	15/20	75.0%	53.1-88.8%
	All Sites	41/59	69.5%	56.9-79.7%
Influenza B/Taiwan/2/62  (Influenza B)	Site 1	7/20	35.0%	18.1-56.7%
	Site 2	9/19	47.4%	27.3-68.3%
	Site 3	8/20	40.0%	21.9-61.3%
	All Sites	24/59	40.7%	29.1-53.4%
Influenza A/SwineNY/03/ 2009  (Influenza A H1N1 pdm09)	Site 1	14/20	70.0%	48.1-85.5%
	Site 2	16/20	80.0%	58.4-91.9%
	Site 3	15/20	75.0%	53.1-88.8%
	All Sites	45/60	75.0%	62.8-84.2%
Influenza A/Port Chalmers/1/73  (Influenza A H3)	Site 1	13/20	65.0%	43.3-81.9%
	Site 2	16/19	84.2%	62.4-94.5%
	Site 3	17/19	89.5%	68.6-97.1%
	All Sites	46/58	79.3%	67.2-87.7%
<i>Mycoplasma pneumoniae</i>  ( <i>M. pneumoniae</i> )	Site 1	13/20	65.0%	43.3-81.9%
	Site 2	14/20	70.0%	48.1-85.5%
	Site 3	14/20	70.0%	48.1-85.5%
	All Sites	41/60	68.3%	55.8-78.7%

Parainfluenza virus 1 (PIV 1)	Site 1	14/20	70.0%	48.1-85.5%
	Site 2	12/19	63.2%	41.0-80.9%
	Site 3	9/19	47.4%	27.3-68.3%
	All Sites	35/58	60.3%	47.5-71.9%
Parainfluenza virus 2 (PIV 2)	Site 1	9/20	45.0%	25.8-65.8%
	Site 2	11/19	57.9%	36.3-76.9%
	Site 3	12/20	60.0%	38.7-78.1%
	All Sites	32/59	54.2%	41.7-66.3%
Parainfluenza virus 3 (PIV 3)	Site 1	13/20	65.0%	43.3-81.9%
	Site 2	17/20	85.0%	64.0-94.8%
	Site 3	17/20	85.0%	64.0-94.8%
	All Sites	47/60	78.3%	66.4-86.9%
Parainfluenza virus 4 (PIV 4)	Site 1	10/20	50.0%	29.9-70.1%
	Site 2	11/19	57.9%	36.3-76.9%
	Site 3	9/20	45.0%	25.8-65.8%
	All Sites	30/59	50.9%	38.4-63.2%
RSV A (RSV)	Site 1	6/20	30.0%	14.5-51.9%
	Site 2	7/20	35.0%	18.1-56.7%
	Site 3	9/20	45.0%	25.8-65.8%
	All Sites	22/60	36.7%	25.6-49.3%
RSV B (RSV)	Site 1	14/20	70.0%	48.1-85.5%
	Site 2	15/19	79.0%	56.7-91.5%
	Site 3	10/20	50.0%	29.9-70.1%
	All Sites	39/59	66.1%	53.4-76.9%
Rhinovirus (Rhinovirus / Enterovirus)	Site 1	15/20	75.0%	53.1-88.8%
	Site 2	15/20	75.0%	53.1-88.8%
	Site 3	18/20	90.0%	69.9-97.2%
	All Sites	48/60	80.0%	68.2-88.2%

**Table 3: QIAstat-Dx Reproducibility with Samples at 1x LOD**

Analyte (QIAstat-Dx Target)	Site	Positive Detected	Percent Agreement with Expected Results	95% CI
Adenovirus (Adenovirus)	Site 1	20/20	100%	83.9-100%
	Site 2	18/18	100%	82.4-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	58/58	100%	93.8-100%
<i>B. pertussis</i> ( <i>B. pertussis</i> )	Site 1	18/20	90%	69.9-97.2%
	Site 2	20/20	100%	83.9-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	58/60	96.7%	88.6-99.1%
<i>C. pneumoniae</i> ( <i>C. pneumoniae</i> )	Site 1	20/20	100%	83.9-100%
	Site 2	20/20	100%	83.9-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	60/60	100%	94.0-100%
Coronavirus 229E (Coronavirus 229E)	Site 1	18/20	90%	69.9-97.2%
	Site 2	20/20	100%	83.9-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	58/60	96.7%	88.6-99.1%
Coronavirus HKU1 (Coronavirus HKU1)	Site 1	20/20	100%	83.9-100%
	Site 2	20/20	100%	83.9-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	60/60	100%	94.0-100%
Coronavirus NL63 (Coronavirus NL63)	Site 1	20/20	100%	83.9-100%
	Site 2	18/18	100%	82.4-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	58/58	100%	93.8-100%
Coronavirus OC43 (Coronavirus OC43)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	59/59	100%	93.9-100%
Enterovirus (Rhinovirus / Enterovirus)	Site 1	19/20	95.0%	76.4-99.1%
	Site 2	20/20	100%	83.9-100%
	Site 3	19/20	95.0%	76.4-99.1%
	All Sites	58/60	97.5%	88.6-99.1%

Human Metapneumovirus  (hMPV)	Site 1	19/20	95.0%	76.4-99.1%
	Site 2	20/20	100.0%	83.9-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	59/60	98.3%	91.1-99.7%
Influenza A/SwineNY/03/ 2009  (Influenza A)	Site 1	20/20	100.0%	83.9-100%
	Site 2	19/19	100.0%	83.2-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	59/59	100.0%	93.9-100%
Influenza A/Port Chalmers/1/73  (Influenza A)	Site 1	19/20	95.0%	76.4-99.1%
	Site 2	18/18	100.0%	82.4-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	57/58	98.3%	90.9-99.7%
Influenza A/NJ/8/76  (Influenza A)	Site 1	19/20	95.0%	76.4-99.1%
	Site 2	20/20	100.0%	83.9-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	59/60	98.3%	91.1-99.7%
Influenza A/NJ/8/76  (Influenza A H1)	Site 1	20/20	100.0%	83.9-100%
	Site 2	20/20	100.0%	83.9-100%
	Site 3	19/20	95.0%	76.4-99.1%
	All Sites	59/60	98.3%	91.1-99.7%
Influenza B/Taiwan/2/62  (Influenza B)	Site 1	19/20	95.0%	76.4-99.1%
	Site 2	20/20	100.0%	83.9-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	59/60	98.3%	91.1-99.7%
Influenza A/SwineNY/03/ 2009  (Influenza A H1N1 pdm09)	Site 1	20/20	100.0%	83.9-100%
	Site 2	19/19	100.0%	83.2-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	59/59	100.0%	93.9-100%
Influenza A/Port Chalmers/1/73  (Influenza A H3)	Site 1	20/20	100.0%	83.9-100%
	Site 2	18/18	100.0%	82.4-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	58/58	100.0%	93.8-100%
<i>Mycoplasma pneumoniae</i>  ( <i>M. pneumoniae</i> )	Site 1	20/20	100.0%	83.9-100%
	Site 2	19/19	100.0%	83.2-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	59/59	100.0%	93.9-100%

Parainfluenza virus 1 (PIV 1)	Site 1	20/20	100.0%	83.9-100%
	Site 2	18/18	100.0%	82.4-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	58/58	100.0%	93.8-100%
Parainfluenza virus 2 (PIV 2)	Site 1	19/20	95.0%	76.4-99.1%
	Site 2	20/20	100.0%	83.9-100%
	Site 3	19/20	95.0%	76.4-99.1%
	All Sites	58/60	96.7%	88.6-99.1%
Parainfluenza virus 3 (PIV 3)	Site 1	20/20	100.0%	83.9-100%
	Site 2	19/19	100.0%	83.2-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	59/59	100.0%	93.9-100%
Parainfluenza virus 4 (PIV 4)	Site 1	20/20	100.0%	83.9-100%
	Site 2	20/20	100.0%	83.9-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	60/60	100.0%	94.0-100%
RSV A (RSV)	Site 1	20/20	100.0%	83.9-100%
	Site 2	19/19	100.0%	83.2-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	59/59	100.0%	93.9-100%
RSV B (RSV)	Site 1	20/20	100.0%	83.9-100%
	Site 2	20/20	100.0%	83.9-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	60/60	100.0%	94.0-100%
Rhinovirus (Rhinovirus / Enterovirus)	Site 1	20/20	100.0%	83.9-100%
	Site 2	19/19	100.0%	83.2-100%
	Site 3	20/20	100.0%	83.9-100%
	All Sites	59/59	100.0%	93.9-100%

**Table 4: QIAstat-Dx Reproducibility with Samples at 3x LOD**

Analyte (QIAstat-Dx Target)	Site	Positive Detected	Percent Agreement with Expected Results	95% CI
Adenovirus (Adenovirus)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	59/59	100%	93.9-100%
<i>B. pertussis</i> ( <i>B. pertussis</i> )	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	59/59	100%	93.9-100%
<i>C. pneumoniae</i> ( <i>C. pneumoniae</i> )	Site 1	20/20	100%	83.9-100%
	Site 2	19/20	95.0%	76.4-99.1%
	Site 3	20/20	100%	83.9-100%
	All Sites	59/60	98.3%	91.1-99.7%
Coronavirus 229E (Coronavirus 229E)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	59/59	100%	93.9-100%
Coronavirus HKU1 (Coronavirus HKU1)	Site 1	20/20	100%	83.9-100%
	Site 2	20/20	100%	83.9-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	60/60	100%	94.0-100%
Coronavirus NL63 (Coronavirus NL63)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	59/59	100%	93.9-100%
Coronavirus OC43 (Coronavirus OC43)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	19/19	100%	83.2-100%
	All Sites	58/58	100%	93.8-100%
Enterovirus (Rhinovirus / Enterovirus)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	59/59	100%	93.9-100%

Human Metapneumovirus  (hMPV)	Site 1	20/20	100%	83.9-100%
	Site 2	20/20	100%	83.2-100%
	Site 3	19/19	100%	83.9-100%
	All Sites	59/59	100%	93.9-100%
Influenza A/SwineNY/03/ 2009  (Influenza A)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	19/19	100%	83.2-100%
	All Sites	58/58	100%	93.8-100%
Influenza A/Port Chalmers/1/73  (Influenza A)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	59/59	100%	93.9-100%
Influenza A/NJ/8/76  (Influenza A)	Site 1	20/20	100%	83.9-100%
	Site 2	20/20	100%	83.9-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	60/60	100%	94.0-100%
Influenza A/NJ/8/76  (Influenza A H1)	Site 1	19/20	95.0%	76.4-99.1%
	Site 2	20/20	100%	83.9-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	59/60	98.3%	91.1-99.7%
Influenza B/Taiwan/2/62  (Influenza B)	Site 1	19/20	95.0%	76.4-99.1%
	Site 2	19/19	100%	83.2-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	58/59	98.3%	91.0-99.7%
Influenza A/SwineNY/03/ 2009  (Influenza A H1N1 pdm09)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	19/19	100%	83.2-100%
	All Sites	58/58	100%	93.8-100%
Influenza A/Port Chalmers/1/73  (Influenza A H3)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	59/59	100%	93.9-100%
<i>Mycoplasma pneumoniae</i>  ( <i>M. pneumoniae</i> )	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	19/19	100%	83.2-100%
	All Sites	58/58	100%	93.8-100%



Parainfluenza virus 1 (PIV 1)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	59/59	100%	93.9-100%
Parainfluenza virus 2 (PIV 2)	Site 1	19/20	95.0%	76.4-99.1%
	Site 2	20/20	100%	83.9-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	59/60	98.3%	91.1-99.7%
Parainfluenza virus 3 (PIV 3)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	19/19	100%	83.2-100%
	All Sites	58/58	100%	93.8-100%
Parainfluenza virus 4 (PIV 4)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	59/59	100%	93.9-100%
RSV A (RSV)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	19/19	100%	83.2-100%
	All Sites	58/58	100%	93.8-100%
RSV B (RSV)	Site 1	20/20	100%	83.9-100%
	Site 2	20/20	100%	83.9-100%
	Site 3	20/20	100%	83.9-100%
	All Sites	60/60	100%	94.0-100%
Rhinovirus (Rhinovirus / Enterovirus)	Site 1	20/20	100%	83.9-100%
	Site 2	19/19	100%	83.2-100%
	Site 3	19/19	100%	83.2-100%
	All Sites	58/58	100%	93.8-100%

**Table 5: Percent Variance Between Site, Day, Operator, and Instrument According to ANOVA for 1x LOD Reproducibility Study**

Analyte	N	Mean Ct	Between Site	Between Instrument	Between Day	Between Operator	Residual	Total
Adenovirus	58	34.2	0.00%	4.30%	1.63%	10.4%	83.7%	100%
<i>B. pertussis</i>	58	35.3	0.00%	0.00%	2.08%	0.00%	97.9%	100%
<i>C. pneumoniae</i>	60	33.7	0.00%	0.00%	0.00%	0.00%	100%	100%
CoV 229E	58	36.1	0.00%	2.58%	0.00%	0.00%	97.4%	100%
CoV HKU1	60	36.3	3.14%	0.00%	8.93%	0.00%	87.9%	100%
CoV NL63	58	35.8	0.00%	0.00%	0.00%	10.5%	89.5%	100%
CoV OC43	59	33.1	0.00%	0.00%	0.00%	1.64%	98.4%	100%

Enterovirus	58	35.8	0.00%	0.00%	0.00%	0.00%	100%	100%
hMPV	59	34.5	0.00%	0.00%	0.00%	0.00%	100%	100%
Influenza A <sup>1</sup>	59	33.5	16.8%	0.00%	0.00%	0.00%	83.2%	100%
Influenza A <sup>2</sup>	57	36.0	0.00%	4.95%	12.0%	0.00%	83.1%	100%
Influenza A <sup>3</sup>	59	36.4	0.00%	0.00%	0.00%	0.00%	100%	100%
Influenza A <sup>4</sup>	59	34.1	0.00%	0.00%	0.00%	4.18%	95.8%	100%
Influenza B	59	33.9	0.00%	2.92%	0.00%	0.00%	97.1%	100%
Influenza H1N1 pdm09 <sup>5</sup>	59	33.2	21.2%	0.00%	0.00%	0.00%	78.8%	100%
Influenza A H3 <sup>6</sup>	58	34.0	2.22%	0.00%	0.00%	2.04%	95.7%	100%
<i>M. pneumoniae</i>	59	33.3	11.3%	0.00%	0.00%	23.2%	65.5%	100%
PIV 1	58	34.9	0.00%	0.00%	0.00%	0.00%	100%	100%
PIV 2	58	33.7	0.00%	0.00%	0.00%	7.88%	92.1%	100%
PIV 3	59	33.3	0.00%	7.62%	0.00%	0.00%	92.4%	100%
PIV 4	60	33.6	0.00%	0.00%	0.00%	0.00%	100%	100%
RSV A	59	34.2	11.3%	0.00%	0.00%	0.00%	88.7%	100%
RSV B	60	33.9	0.00%	0.00%	2.45%	0.00%	97.6%	100%
Rhinovirus	59	33.7	0.00%	0.00%	0.00%	3.02%	97.0%	100%

<sup>1</sup> Strain A/Swine NY/03/2009; Flu A target.

<sup>2</sup> Strain A/Port Chalmers/1/73' Flu A target.

<sup>3</sup> Strain A/NJ/8/1976; Flu A target.

<sup>4</sup> Strain A/NJ/8/1976; H1 target.

<sup>5</sup> Strain A/Swine NY/03/2009; pdm09 target.

<sup>6</sup> Strain A/Port Chalmers/1/73; H3 target.

**Table 6: Percent Variance Between Site, Day, Operator, and Instrument According to ANOVA for 3x LOD Reproducibility Study**

Analyte	N	Mean Ct	Between Site	Between Instrument	Between Day	Between Operator	Residual	Total
Adenovirus	59	33.1	0.00%	32.3%	4.56%	0.00%	63.1%	100%
<i>B. pertussis</i>	59	34.2	0.00%	22.7%	1.23%	0.00%	76.1%	100%
<i>C. pneumoniae</i>	59	33.2	0.00%	0.00%	2.32%	0.00%	97.7%	100%
CoV 229E	59	34.6	12.4%	0.00%	3.77%	0.13%	83.7%	100%
CoV HKU1	60	34.8	28.2%	7.31%	0.00%	0.32%	64.2%	100%
CoV NL63	59	34.0	0.00%	7.24%	0.00%	0.00%	92.8%	100%
CoV OC43	58	31.1	1.36%	0.00%	1.53%	0.00%	97.1%	100%
Enterovirus	59	34.1	0.00%	7.02%	0.80%	0.00%	92.2%	100%
hMPV	59	33.2	1.11%	0.00%	0.00%	0.00%	98.9%	100%
Influenza A <sup>1</sup>	58	31.8	0.00%	3.43%	0.00%	0.15%	96.4%	100%
Influenza A <sup>2</sup>	59	34.4	0.00%	0.00%	0.99%	0.00%	99.0%	100%
Influenza A <sup>3</sup>	60	35.1	7.75%	0.00%	3.86%	0.00%	88.4%	100%
Influenza A <sup>4</sup>	59	32.6	0.00%	0.01%	0.00%	0.28%	99.7%	100%
Influenza B	58	32.2	9.20%	0.00%	0.00%	4.58%	86.2%	100%
Influenza H1N1 pdm09 <sup>5</sup>	58	31.5	5.65%	1.97%	0.00%	0.35%	92.0%	100%
Influenza A H3 <sup>6</sup>	59	32.4	0.00%	1.76%	0.00%	0.00%	98.2%	100%
<i>M. pneumoniae</i>	58	31.5	0.00%	0.00%	0.00%	0.00%	100%	100%
PIV 1	59	33.6	4.70%	0.00%	0.00%	0.00%	95.3%	100%
PIV 2	59	32.5	0.00%	1.96%	0.00%	0.04%	98.0%	100%
PIV 3	58	31.8	0.00%	0.00%	9.47%	0.00%	90.5%	100%
PIV 4	59	31.9	0.00%	0.00%	0.00%	0.00%	100%	100%
RSV A	58	32.8	2.59%	0.00%	10.5%	0.00%	86.9%	100%
RSV B	60	32.7	0.00%	8.47%	0.00%	0.76%	90.8%	100%
Rhinovirus	58	32.1	2.77%	0.00%	0.00%	0.00%	97.2%	100%

- <sup>1</sup> Strain A/Swine NY/03/2009; Flu A target.
- <sup>2</sup> Strain A/Port Chalmers/1/73' Flu A target.
- <sup>3</sup> Strain A/NJ/8/1976; Flu A target.
- <sup>4</sup> Strain A/NJ/8/1976; H1 target.
- <sup>5</sup> Strain A/Swine NY/03/2009; pdm09 target.
- <sup>6</sup> Strain A/Port Chalmers/1/73; H3 target.

*b. Linearity/assay reportable range:*

Not applicable; this is a qualitative assay.

*c. Traceability, Stability, Expected values (controls, calibrators, or methods):*

Specimen Stability

To provide data supporting the specimen storage recommendations stated in the product package insert, an analytical study was carried out to evaluate specimen stability.

Contrived positive NPS samples (eluted NPS samples) were prepared by diluting specific panel targets into a sample medium consisting of UTM combined with HeLa cells. Sample mixes were prepared at both 5x LOD and 1x LOD concentrations in discrete sample mixes. Negative data was obtained from specimen mixes where the analyte was expected to be absent based on the composition of the specific mix. All prepared samples were tested at N=10 at each of the following three time points/conditions: Time 0 (fresh), 4 hours at 15 to 25 °C, 72 hours at 2 to 8 °C, and 30 days at -15 to -25 °C. Positive samples were considered stable as long as they tested positive in the QIAstat-Dx Respiratory Panel with at least a 90% detection rate. The acceptance criteria of ≥90% detection was achieved for all analytes at the 5x LOD concentration tested under the four conditions described above. Details regarding the detection rates are presented in the tables below.

**Table 7a: Specimen Stability Study Results – 5x LOD**

Mix	Analyte	Fresh	RT 4h	2 to 8 °C 72h	-15 to -25 °C 30d
Mix 1	Influenza A	9/10	10/10	9/10	10/10
	Influenza A H1	9/10	10/10	10/10	10/10
	Cor HKU1	9/10	10/10	9/10	10/10
	PIV 2	10/10	10/10	9/10	10/10
	RSV B	10/10	10/10	10/10	10/10
	<i>C. pneumoniae</i>	10/10	10/10	10/10	10/10
Mix 2	Influenza B	10/10	10/10	10/10	10/10
	Cor 229E	10/10	10/10	10/10	10/10
	PIV 4	10/10	10/10	10/10	10/10
	Enterovirus D68	9/10	10/10	9/10	10/10
	hMPV A1	10/10	10/10	10/10	10/10
	<i>B. pertussis</i>	10/10	10/10	10/10	10/10
Mix 3	Influenza A	10/10	10/10	10/10	10/10
	Influenza H1N1pdm09	10/10	10/10	10/10	10/10
	Cor OC43	10/10	10/10	10/10	10/10
	PIV 3	10/10	10/10	9/10	10/10
	Rhinovirus A2	10/10	10/10	10/10	10/10

	RSV A	10/10	10/10	10/10	10/10
	<i>M. pneumoniae</i>	10/10	10/10	10/10	10/10
Mix 4	Influenza A	10/10	9/10	10/10	10/10
	Influenza A H3	10/10	9/10	10/10	10/10
	Cor NL63	9/10	9/10	10/10	10/10
	PIV 1	9/10	9/10	10/10	10/10
	Adenovirus B3	10/10	9/10	9/10	10/10

While some variance is expected due to the fact that 1x LOD is essentially equivalent to a 95% detection rate for the analyte, the acceptance criteria was not met for all analytes and all time points/conditions when tested at the 1x LOD concentration.

**Table 7b: Specimen Stability Study Results – 1x LOD**

Mix	Analyte	Fresh	RT 4h	2 to 8 °C 72h	-15 to -25 °C 30d
Mix 5	Influenza A	6/10	7/10	9/10	9/10
	Influenza A H1	8/10	9/10	9/10	9/10
	Cor HKU1	8/10	9/10	9/10	10/10
	PIV 2	9/10	9/10	9/10	9/10
	RSV B	9/10	10/10	9/10	10/10
	<i>C. pneumoniae</i>	9/10	8/10	7/10	10/10
Mix 6	Influenza B	10/10	9/10	8/10	10/10
	Cor 229E	8/10	9/10	9/10	10/10
	PIV 4	9/10	9/10	9/10	10/10
	Enterovirus D68	7/10	9/10	6/10	9/10
	hMPV A1	10/10	9/10	6/10	10/10
	<i>B. pertussis</i>	7/10	8/10	9/10	10/10
Mix 7	Influenza A	10/10	10/10	10/10	10/10
	Influenza H1N1pdm09	10/10	8/10	9/10	10/10
	Cor OC43	9/10	9/10	10/10	10/10
	PIV 3	10/10	9/10	9/10	10/10
	Rhinovirus A2	10/10	9/10	10/10	10/10
	RSV A	10/10	9/10	9/10	10/10
	<i>M. pneumoniae</i>	10/10	8/10	10/10	10/10
Mix 8	Influenza A	8/10	9/10	10/10	9/10
	Influenza A H3	9/10	10/10	9/10	10/10
	Cor NL63	8/10	10/10	10/10	9/10
	PIV 1	7/10	8/10	8/10	9/10
	Adenovirus B3	4/10	7/10	9/10	9/10

Supplemental testing was performed for select analytes at select conditions at the 1x LOD concentration to further test specimen stability where acceptance criteria had not been previously met.

**Table 7c: Additional Specimen Stability Study Results – 1x LOD**

Mix	Analyte	Fresh	RT 4h	2 to 8 °C 72h	-15 to -25 °C 30d
Mix 9	Influenza A	9/10	10/10	nd	nd
	Influenza A H1	7/10	9/10	nd	nd
	Adenovirus B3	8/10	9/10	nd	nd
Mix 10	Enterovirus D68	9/10	nd	8/10	nd

	hMPV A1	10/10	nd	10/10	nd
	<i>C. pneumoniae</i>	10/10	nd	10/10	nd

nd – Testing not done.

### Shelf Life

The stated shelf life for the QIAstat-Dx Respiratory Panel is 12 months when stored at 15 to 25 °C. Stability data to support the proposed shelf life and shipping conditions was obtained by testing three separate production lots of QIAstat-Dx cartridges using QC material under a real time stability study. The study was performed at room temperature at 30-day intervals spanning a planned 13-month period. Data generation is ongoing and is being maintained under Qiagen’s Quality Systems internal protocols.

#### *d. Detection limit:*

The objective of the Analytical Sensitivity Study was to identify the limit of detection (LOD) of the QIAstat-Dx Respiratory Panel prepared from high-titer stocks obtained from commercial suppliers or clinical isolates for commercially unavailable target analytes. For the purposes of the study, the LOD level was defined as the concentration of analyte that produced positive QIAstat-Dx Respiratory Panel test results approximately 95% of the time when tested in multiple replicates.

The LOD was assessed in a two-step process for every analyte. The first step was setup of a preliminary LOD by testing serial dilutions for every pathogen in a minimum of 3 log serial dilutions around the expected LOD. For this step of the study, samples were prepared in a simulated matrix consisting of UTM plus HeLa cells and four replicates of each analyte were tested. The data obtained from the first step was then used to choose the concentration of analyte likely to provide a minimum of 19 out of 20 positive results for the LOD confirmation. Confirmed LOD concentrations prepared in simulated matrix were then verified by testing the analytes in clinical matrix whereby at least 19 out of 20 replicates were again detected by the assay (see “matrix equivalency” below). If the assay failed to verify the LOD in clinical matrix, a sample of 10x more concentrated titer was prepared and retested for LOD verification in clinical matrix. The concentrations of analytes which provided at least 19/20 positive results in clinical matrix is listed as the claimed “LOD” in table 8.

**Table 8: QIAstat-Dx Respiratory Panel Limits of Detection**

Pathogen	Strain	LOD	Units
Influenza A H1N1	A/NJ/8/76	341	CEID <sub>50</sub> /mL
	A/Brisbane/59/07	4	TCID <sub>50</sub> /mL
	A/New Caldonia/20/99	15	TCID <sub>50</sub> /mL
Influenza A H3N2	A/Virginia/ATCC6/2012	0.1	PFU/mL
	A/Wisconsin/67/2005	3.8	TCID <sub>50</sub> /mL
	A/Port Chalmers/1/73	499	CEID <sub>50</sub> /mL
Influenza A H1N1/pdm09	A/Virginia/ATCC1/2009	6.7	PFU/mL
	A/SwineNY/03/2009	5.6	TCID <sub>50</sub> /mL
Influenza B	B/Virginia/ATCC5/2012	0.03	PFU/mL
	B/FL/04/06	1080	CEID <sub>50</sub> /mL
	B/Taiwan/2/62	5000	CEID <sub>50</sub> /mL
Coronavirus 229E	n/a	0.2	TCID <sub>50</sub> /mL
	n/a	3.6	TCID <sub>50</sub> /mL
Coronavirus OC43	n/a	0.1	TCID <sub>50</sub> /mL
	n/a	0.1	TCID <sub>50</sub> /mL
Coronavirus NL63	n/a	0.01	TCID <sub>50</sub> /mL
Coronavirus HKU1	n/a	40,000	Copies/mL
Parainfluenza Virus 1	C35	0.2	TCID <sub>50</sub> /mL
	n/a	0.2	TCID <sub>50</sub> /mL
Parainfluenza Virus 2	Greer	7.3	TCID <sub>50</sub> /mL
	n/a	1.3	TCID <sub>50</sub> /mL
Parainfluenza Virus 3	C 243	2.3	TCID <sub>50</sub> /mL
	n/a	11.5	TCID <sub>50</sub> /mL
Parainfluenza Virus 4a	M-25	0.5	TCID <sub>50</sub> /mL
Parainfluenza Virus 4b	n/a	9.5	TCID <sub>50</sub> /mL
RSV A	A2	12.0	PFU/mL
	Long	33.0	PFU/mL
RSV B	18537	0.03	PFU/mL
	CH93(18)-18	0.4	TCID <sub>50</sub> /mL
Human metapneumovirus	Peru6-2003	0.01	TCID <sub>50</sub> /mL
	IA10-2003	0.5	TCID <sub>50</sub> /mL
	IA14-2003	0.4	TCID <sub>50</sub> /mL
	Peru2-2002	1480	TCID <sub>50</sub> /mL
Adenovirus	GB	4993	TCID <sub>50</sub> /mL
	RI-67	15.8	TCID <sub>50</sub> /mL
	Adenoid 75	7331	TCID <sub>50</sub> /mL
	Adenoid 71	69.5	TCID <sub>50</sub> /mL
	Adenoid 6	28.1	TCID <sub>50</sub> /mL
	Tonsil 99	88.8	TCID <sub>50</sub> /mL
Enterovirus	US/IL/14-18952	8.9	TCID <sub>50</sub> /mL
	Echovirus 6	0.9	TCID <sub>50</sub> /mL
Rhinovirus	1059	8.9	TCID <sub>50</sub> /mL
	HGP	8.9	TCID <sub>50</sub> /mL
	11757	50.0	TCID <sub>50</sub> /mL
	Type 1A	8.9	TCID <sub>50</sub> /mL
<i>Mycoplasma pneumoniae</i>	M129-B7	0.1	CCU/mL
	PI 1428	1.0	CCU/mL
<i>Chlamydia pneumoniae</i>	TW183	14.2	IFU/mL
	CWL-029	120	IFU/mL
<i>Bordetella pertussis</i>	I028	0.3	CFU/mL
	18323	2.6	CFU/mL

*e. Analytical reactivity:*

Various respiratory panel virus and bacteria strains were tested to examine the ability of the QIAstat-Dx Respiratory Panel to detect a wide variety of analyte strains in a clinical setting. Samples were prepared in a simulated matrix consisting of UTM plus HeLa cells and tested according to the package insert.

Log dilutions of each analyte were prepared as sample mixes containing multiple analytes in each mix. All analyte dilutions were run in triplicate. Acceptance criterion for the study was 3/3 positive results for the analyte being tested. If the criterion was not met, a 10x more concentrated titer of the analyte was tested in triplicate. Results in tables 9 - 20 show the analyte, strain, and concentration at which the acceptance criterion was met. The lowest level of each strain that generated positive results on all three replicates was identified as the lowest level detected by the QIAstat-Dx Respiratory Panel.

**Table 9: QIAstat-Dx RP Analytical Reactivity Results for Influenza A H1**

Analyte	Strain	Concentration	Detected / Tested	QIAstat-Dx RP Result
Influenza A H1N1	A/Brisbane/59/07	0.4 TCID <sub>50</sub> /mL	3/3	Influenza A H1
	A/New Caldonia/20/99	1.5 TCID <sub>50</sub> /mL	3/3	Influenza A H1
	A/NJ/8/76	34.1 CEID <sub>50</sub> /mL	3/3	Influenza A H1
	A/Denver/1/57	340 CEID <sub>50</sub> /mL	3/3	Influenza A H1
	A/Mal/302/54	15.8 CEID <sub>50</sub> /mL	3/3	Influenza A H1
	A/Weiss/43	28117 CEID <sub>50</sub> /mL	3/3	Influenza A H1
	A/PR/8/34	390 PFU/mL	3/3	Influenza A H1
	A/Fort Monmouth/1/1947	28.1 CEID <sub>50</sub> /mL	3/3	Influenza A H1
	A/WS/33	15.8 TCID <sub>50</sub> /mL	3/3	Influenza A H1
A/Swine/Iowa/15/1930	889 CEID <sub>50</sub> /mL	3/3	Influenza A H1	

**Table 10: QIAstat-Dx RP Analytical Reactivity Results for Influenza A H3**

Analyte	Strain	Concentration	Detected / Tested	QIAstat-Dx RP Result
Influenza A H3N2	A/Port Chalmers/1/73	499 CEID <sub>50</sub> /mL	3/3	Influenza A H3
	A/Virginia/ATCC6/2012	0.1 PFU/mL	3/3	Influenza A H3
	A/Wisconsin/67/2005	3.8 TCID <sub>50</sub> /mL	3/3	Influenza A H3
	A/Wisconsin/15/2009	5.8 CEID <sub>50</sub> /mL	3/3	Influenza A H3
	A/Victoria/3/75	16 CEID <sub>50</sub> /mL	3/3	Influenza A H3
	A/Aichi/2/68	31 PFU/mL	3/3	Influenza A H3
	A/Hong Kong/8/68	1581 TCID <sub>50</sub> /mL	3/3	Influenza A H3
	A/Alice <sup>1</sup>	500 TCID <sub>50</sub> /mL	3/3	Influenza A H3
	MRC-2 <sup>2</sup>	8891 CEID <sub>50</sub> /mL	3/3	Influenza A H3
A/Switzerland/9715293/2013	1000 CEID <sub>50</sub> /mL	3/3	Influenza A H3	

<sup>1</sup> Recombinant strain; carries A/England/42/72 genes

<sup>2</sup> Recombinant strain; carries A/England/42/72 and A/PR/8/34 genes

**Table 11: QIAstat-Dx RP Analytical Reactivity Results for Influenza A H1N1 pdm09**

Analyte	Strain	Concentration	Detected / Tested	QIAstat-Dx RP Result
Influenza A H1N1 pdm09	A/Virginia/ATCC1/2009	6.7 PFU/mL	3/3	Influenza A H1N1/pdm09
	A/SwineNY/03/2009	5.6 TCID <sub>50</sub> /mL	3/3	Influenza A H1N1/pdm09
	A/Virginia/ATCC2/2009	61 PFU/mL	3/3	Influenza A H1N1/pdm09
	A/Virginia/ATCC3/2009	1800 PFU/mL	3/3	Influenza A H1N1/pdm09
	Swine NY/01/2009	138 TCID <sub>50</sub> /mL	3/3	Influenza A H1N1/pdm09
	Swine NY/02/2009	1.4 TCID <sub>50</sub> /mL	3/3	Influenza A H1N1/pdm09
	A/California/07/2009	1400 CEID <sub>50</sub> /mL	3/3	Influenza A H1N1/pdm09
	Canada/6294/09	1.7 TCID <sub>50</sub> /mL	3/3	Influenza A H1N1/pdm09
	Mexico/4108/09	14.1 TCID <sub>50</sub> /mL	3/3	Influenza A H1N1/pdm09
	Netherlands/2629/2009	16 TCID <sub>50</sub> /mL	3/3	Influenza A H1N1/pdm09

**Table 12: QIAstat-Dx RP Analytical Reactivity Results for Influenza A**

Analyte	Strain	Concentration	Detected / Tested	QIAstat-Dx RP Result
Influenza A H2N2	Japan/305/1957 <sup>1</sup>	3.26 x 10 <sup>-3</sup> ng/μL	3/3	Influenza A (no subtype)
	Korea/426/1968 <sup>2</sup>	6.25 x 10 <sup>-5</sup> ng/μL	3/3	Influenza A (no subtype)
Influenza A H5N3	A/Duck/Singapore/645/1997 <sup>3</sup>	2.48 x 10 <sup>-3</sup> ng/μL	3/3	Influenza A (no subtype)
Influenza A H10N7	Chicken/Germany/N/49 <sup>3</sup>	6.80 x 10 <sup>-2</sup> ng/μL	3/3	Influenza A (no subtype)
Influenza A H1N2	Recombinant Kilbourne F63 <sup>4</sup>	1.48 x 10 <sup>-2</sup> ng/μL	3/3	Influenza A H1

<sup>1</sup> Nucleic acid

<sup>2</sup> Nucleic acid; recombinant cross with A/PR/8/34

<sup>3</sup> Nucleic acid; avian source

<sup>4</sup> Nucleic acid; recombinant cross of A/NWS/1934 x A/Rockefeller Institute/5/1957



**Table 13: QIAstat-Dx RP Analytical Reactivity Results for Influenza B**

Analyte	Strain	Concentration	Detected / Tested	QIAstat-Dx RP Result
Influenza B	B/Virginia/ATCC5/2012	0.03 PFU/mL	3/3	Influenza B
	B/FL/04/06	108 CEID <sub>50</sub> /mL	3/3	Influenza B
	B/Taiwan/2/62	49.9 CEID <sub>50</sub> /mL	3/3	Influenza B
	B/Allen/45	n/a	0/3	Negative
	B/Hong Kong/5/72	n/a	0/3	Negative
	B/Maryland/1/59	338 CEID <sub>50</sub> /mL	3/3	Influenza B
	B/GL/1739/54	50.0 CEID <sub>50</sub> /mL	3/3	Influenza B
	B/Wisconsin/1/2010	0.3 CEID <sub>50</sub> /mL	3/3	Influenza B
	B/Massachusetts/2/2012	2300 CEID <sub>50</sub> /mL	3/3	Influenza B
	B/Florida/02/06	n/a	1/3	Influenza B and Negative <sup>1</sup>
	B/Brisbane/60/2008	1.8 CEID <sub>50</sub> /mL	3/3	Influenza B
B/Malaysia/2506/2004	1.58 CEID <sub>50</sub> /mL	3/3	Influenza B	

<sup>1</sup> Multiple dilutions and repeats failed to achieve acceptance criterion of 3/3

**Table 14: QIAstat-Dx RP Analytical Reactivity Results for Coronavirus**

Analyte	Strain	Concentration	Detected / Tested	QIAstat-Dx RP Result
Coronavirus 229E	n/a	3.6 TCID <sub>50</sub> /mL	3/3	Coronavirus 229E
	n/a	0.2 TCID <sub>50</sub> /mL	3/3	Coronavirus 229E
Coronavirus OC43	n/a	0.1 TCID <sub>50</sub> /mL	3/3	Coronavirus OC43
	n/a	0.1 TCID <sub>50</sub> /mL	3/3	Coronavirus OC43
Coronavirus NL63	n/a	0.01 TCID <sub>50</sub> /mL	3/3	Coronavirus NL63
	n/a	1.6 TCID <sub>50</sub> /mL	3/3	Coronavirus NL63
Coronavirus HKU1	n/a	3.0 x 10 <sup>4</sup> copies/mL	3/3	Coronavirus HKU1
	Clinical isolate	4.0 x 10 <sup>8</sup> copies/mL	3/3	Coronavirus HKU1
	Clinical isolate	7.0 x 10 <sup>7</sup> copies/mL	3/3	Coronavirus HKU1
	Clinical isolate	7.0 x 10 <sup>7</sup> copies/mL	3/3	Coronavirus HKU1

**Table 15: QIAstat-Dx RP Analytical Reactivity Results for Parainfluenza Virus**

Analyte	Strain	Concentration	Detected / Tested	QIAstat-Dx RP Result
Parainfluenza Virus 1	n/a	0.02 TCID <sub>50</sub> /mL	3/3	PIV 1
	C35	0.2 TCID <sub>50</sub> /mL	3/3	PIV 1
	n/a	n/a <sup>1</sup>	3/3	PIV 1
Parainfluenza Virus 2	Greer	2.3 TCID <sub>50</sub> /mL	3/3	PIV 2
	n/a	1.3 TCID <sub>50</sub> /mL	3/3	PIV 2
	n/a	1.3 TCID <sub>50</sub> /mL	3/3	PIV 2
Parainfluenza Virus 3	n/a	11.5 TCID <sub>50</sub> /mL	3/3	PIV 3
	C 243	2.3 TCID <sub>50</sub> /mL	3/3	PIV 3
	n/a	n/a <sup>1</sup>	3/3	PIV 3
Parainfluenza Virus 4	M-25	0.5 TCID <sub>50</sub> /mL	3/3	PIV 4
	n/a	9.6 TCID <sub>50</sub> /mL	3/3	PIV 4
	n/a	28.2 TCID <sub>50</sub> /mL	3/3	PIV 4
	CH 19503	1 TCID <sub>50</sub> /mL	3/3	PIV 4

<sup>1</sup> Stock titer not available from supplier.

**Table 16: QIAstat-Dx RP Analytical Reactivity Results for RSV A+B**

Analyte	Strain	Concentration	Detected / Tested	QIAstat-Dx RP Result
Respiratory Syncytial Virus A+B	18537	0.03 PFU/mL	3/3	RSV A+B
	A2	12 PFU/mL	3/3	RSV A+B
	Long	33 PFU/mL	3/3	RSV A+B
	CH93(18)-18	0.4 TCID <sub>50</sub> /mL	3/3	RSV A+B
	n/a	0.3 TCID <sub>50</sub> /mL	3/3	RSV A+B
	B WV/14617/85	15.8 TCID <sub>50</sub> /mL	3/3	RSV A+B

**Table 17: QIAstat-Dx RP Analytical Reactivity Results for hMPV A+B**

Analyte	Strain	Concentration	Detected / Tested	QIAstat-Dx RP Result
Human Metapneumovirus	IA10-2003	0.5 TCID <sub>50</sub> /mL	3/3	hMPV A+B
	IA14-2003	0.4 TCID <sub>50</sub> /mL	3/3	hMPV A+B
	Peru2-2002	1479 TCID <sub>50</sub> /mL	3/3	hMPV A+B
	Peru6-2003	0.01 TCID <sub>50</sub> /mL	3/3	hMPV A+B
	IA3-2002	66 TCID <sub>50</sub> /mL	3/3	hMPV A+B
	IA27-2004	1.3 TCID <sub>50</sub> /mL	3/3	hMPV A+B
	Peru3-2003	31.6 TCID <sub>50</sub> /mL	3/3	hMPV A+B
	IA8-2003	0.4 TCID <sub>50</sub> /mL	3/3	hMPV A+B
	Peru1-2002	2188 TCID <sub>50</sub> /mL	3/3	hMPV A+B

**Table 18: QIAstat-Dx RP Analytical Reactivity Results for Adenovirus**

Analyte	Strain	Concentration	Detected / Tested	QIAstat-Dx RP Result
Adenovirus	Tonsil 99	88.8 TCID <sub>50</sub> /mL	3/3	Adenovirus
	GB	4993 TCID <sub>50</sub> /mL	3/3	Adenovirus
	Adenoid 71	69.5 TCID <sub>50</sub> /mL	3/3	Adenovirus
	Adenoid 6	28.1 TCID <sub>50</sub> /mL	3/3	Adenovirus
	Adenoid 75	7331 TCID <sub>50</sub> /mL	3/3	Adenovirus
	RI-67	15.8 TCID <sub>50</sub> /mL	3/3	Adenovirus
	Huie	88.9 TCID <sub>50</sub> /mL	3/3	Adenovirus
	Gomen	0.3 TCID <sub>50</sub> /mL	3/3	Adenovirus
	Slobitski	16 TCID <sub>50</sub> /mL	3/3	Adenovirus
	AV-1645 [128]	2.8 TCID <sub>50</sub> /mL	3/3	Adenovirus
	Compton	0.28 TCID <sub>50</sub> /mL	3/3	Adenovirus
	Holden	8.9 TCID <sub>50</sub> /mL	3/3	Adenovirus
	Trim	160 TCID <sub>50</sub> /mL	3/3	Adenovirus
Dugan	0.2 TCID <sub>50</sub> /mL	3/3	Adenovirus	
Tak (73-3544)	28117 TCID <sub>50</sub> /mL	3/3	Adenovirus	

**Table 19: QIAstat-Dx RP Analytical Reactivity Results for Rhinovirus/Enterovirus**

Analyte	Strain	Concentration	Detected / Tested	QIAstat-Dx RP Result
Enterovirus	US/IL/14-18952	8.9 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	D-1 (Cox)	0.9 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	H	8.9 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	M.K. (Kowalik)	n/a	3/3	Rhinovirus/Enterovirus
	Gregory	889 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	Bastianni	281 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	Griggs	1.6 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	Conn-5	158 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	Ohio-1	2812 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	Nancy	0.9 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	CHHE-29	0.03 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	Kuykendall	28.1 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
Rhinovirus	1059	8.9 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	2060	8.9 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	HGP	8.9 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	11757	49.9 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	FEB	281 TCID <sub>50</sub> /mL	3/3	Rhinovirus/Enterovirus
	33342	200 PFU/mL	3/3	Rhinovirus/Enterovirus

**Table 20: QIAstat-Dx RP Analytical Reactivity Results for *Mycoplasma pneumoniae*, *Bordetella pertussis*, and *Chlamydomphila pneumoniae***

Analyte	Strain	Concentration	Detected / Tested	QIAstat-Dx RP Result
<i>M. pneumoniae</i>	PI 1428	1 CCU/mL	3/3	<i>M. pneumoniae</i>
	M129-B7	0.1 CCU/mL	3/3	<i>M. pneumoniae</i>
	FH	0.2 CCU/mL	3/3	<i>M. pneumoniae</i>
<i>B. pertussis</i>	I028	0.3 CFU/mL	3/3	<i>B. pertussis</i>
	19323	2.6 CFU/mL	3/3	<i>B. pertussis</i>
	10-536	n/a	3/3	<i>B. pertussis</i>
<i>C. pneumoniae</i>	TW183	14.2 IFU/mL	3/3	<i>C. pneumoniae</i>
	CWL-029	120 IFU/mL	3/3	<i>C. pneumoniae</i>
	AR-39	29 IFU/mL	3/3	<i>C. pneumoniae</i>

*f. Analytical specificity:*

To determine the analytical specificity of the QIAstat-Dx Respiratory Panel, 21 on-panel pathogens and 52 off-panel pathogens were tested for any potential to cross-react with primers and probes specific for other analytes in the assay. Viral targets were tested at 10<sup>5</sup> units/mL and bacteria/fungal targets were tested at 10<sup>6</sup> units/mL wherever possible. Two off-panel bacterial targets were tested at lower concentrations due to limits of availability from the supplier: *Bordetella hinzii* and *Legionella feeleii* were tested at 5.0 x 10<sup>3</sup> CFU/mL and 1.0 x 10<sup>4</sup> CFU/mL, respectively. The on-panel and off-panel testing samples were prepared by single spiking organisms into simulated NPS sample matrix (UTM + HeLa cells). All organisms were tested in triplicate using three different lots of QIAstat-Dx RP cartridges and up to 22 different analyzers. Acceptance criteria for the on-panel pathogens required all replicates to

provide a positive result for the specific target present in the sample and a negative result for all targets absent from the sample. Tables of on-panel and off-panel organisms used in this study are presented below.

**Table 21: On-Panel Targets for QIAstat-Dx RP Analytical Specificity**

Pathogen	Strain	Concentration Tested
Influenza A H1N1	A/NJ/8/76	1.0 x10 <sup>5</sup> CEID <sub>50</sub> /mL
Influenza A H3N2	A/Virginia/ATCC6/2012	1.0 x10 <sup>5</sup> PFU/mL
Influenza A/2009/H1N1	A/Virginia/ATCC1/2009	1.0 x10 <sup>5</sup> PFU/mL
Influenza B	B/FL/04/06	1.0 x10 <sup>5</sup> CEID <sub>50</sub> /mL
Coronavirus 229E	n/a	1.0 x10 <sup>5</sup> TCID <sub>50</sub> /mL
Coronavirus OC43	n/a	1.0 x10 <sup>5</sup> TCID <sub>50</sub> /mL
Coronavirus NL63	n/a	1.0 x10 <sup>5</sup> TCID <sub>50</sub> /mL
Coronavirus HKU1	n/a	1.0 x10 <sup>5</sup> Copies/mL
Parainfluenza Virus 1	C35	1.0 x10 <sup>5</sup> TCID <sub>50</sub> /mL
Parainfluenza Virus 2	Greer	1.0 x10 <sup>5</sup> TCID <sub>50</sub> /mL
Parainfluenza Virus 3	C 243	1.0 x10 <sup>5</sup> TCID <sub>50</sub> /mL
Parainfluenza Virus 4	PIV 4a	1.0 x10 <sup>5</sup> TCID <sub>50</sub> /mL
RSV A	A2	1.0 x10 <sup>5</sup> TCID <sub>50</sub> /mL
hMPV	IA10-2003	1.0 x10 <sup>5</sup> TCID <sub>50</sub> /mL
Adenovirus	Adenoid 71	1.0 x10 <sup>5</sup> TCID <sub>50</sub> /mL
Adenovirus	Gomen	1.0 x10 <sup>5</sup> TCID <sub>50</sub> /mL
Enterovirus	US/IL/14-18952	1.0 x10 <sup>5</sup> TCID <sub>50</sub> /mL
Rhinovirus	Type 1A	1.0 x10 <sup>5</sup> TCID <sub>50</sub> /mL
<i>Mycoplasma pneumoniae</i>	M129	1.0 x10 <sup>6</sup> CCU/mL
<i>Bordetella pertussis</i>	E431	1.0 x10 <sup>6</sup> CFU/mL
<i>Chlamydomphila pneumoniae</i>	AR-39	1.0 x10 <sup>6</sup> IFU/mL

All results from the on-panel target list met the acceptance criteria at the concentrations tested.

**Table 22: Off-Panel Targets for QIAstat-Dx RP Analytical Specificity**

<b>Bacteria</b>	<b>Bacteria (continued)</b>	<b>Viruses</b>	<b>Fungi</b>
<i>Acinetobacter calcoaceticus</i>	<i>Moraxella catarrhalis</i>	Bocavirus <sup>2</sup>	<i>Aspergillus flavis</i>
<i>Bordetella avium</i>	<i>Mycobacterium tuberculosis</i> <sup>1</sup>	Cytomegalovirus	<i>Aspergillus fumigatus</i>
<i>Bordetella bronchioseptica</i>	<i>Mycoplasma hominis</i>	Epstein-Barr virus	<i>Candida albicans</i>
<i>Bordetella hinzii</i>	<i>Mycoplasma orale</i>	HSV-1	<i>Cryptococcus neoformans</i>
<i>Bordetella holmesii</i>	<i>Neisseria elongata</i>	HSV-2	
<i>Bordetella parapertussis</i>	<i>Neisseria gonorrhoeae</i>	Measles virus	
<i>Chlamydia trachomatis</i>	<i>Neisseria meningitidis</i>	MERS CoV <sup>3</sup>	
<i>Corynebacterium diphtheriae</i>	<i>Proteus mirabilis</i>	Mumps virus	
<i>Enterobacter aerogenes</i>	<i>Pseudomonas aeruginosa</i>		
<i>Escherichia coli</i>	<i>Serratia marcescens</i>		
<i>Haemophilus influenzae</i>	<i>Staphylococcus aureus</i>		
<i>Klebsiella pneumoniae</i>	<i>Staphylococcus epidermidis</i>		
<i>Klebsiella oxytoca</i>	<i>Stenotrophomonas maltophilia</i>		
<i>Lactobacillus acidophilus</i>	<i>Streptococcus agalactiae</i>		
<i>Lactobacillus plantarum</i>	<i>Streptococcus pneumonia</i>		
<i>Legionella bozemanii</i>	<i>Streptococcus salivarius</i>		
<i>Legionella dumofii</i>	<i>Streptococcus pyogenes</i>		
<i>Legionella feeleii</i>	<i>Ureaplasma urealyticum</i>		
<i>Legionella longbeacheae</i>			
<i>Legionella micdadei</i>			
<i>Legionella pneumophila</i>			

<sup>1</sup> Genomic DNA

<sup>2</sup> Clinical isolate

<sup>3</sup> Synthetic RNA

The following false positive results were observed in the off-panel testing: Positive *M. pneumoniae* results were observed for *Enterobacter aerogenes* (1/9), *Streptococcus pyogenes* (2/9), and *Aspergillus fumigatus* (1/3); a positive Rhinovirus/Enterovirus result was observed for *Legionella micdadei* (1/6); and positive *Bordetella pertussis* results were observed for *Bordetella bronchioseptica* (3/3) and *Bordetella holmesii* (3/3). The unexpected positive results for *M. pneumoniae* may be due to contamination of the off-panel pathogen sources with *M. pneumoniae* analyte. Cross-reactivity observed with other *Bordetella* species is likely due to the insertion sequence transposon which is the molecular target for *Bordetella pertussis* (IS481). IS481 is found in other *Bordetella* species, albeit typically in fewer copy numbers. A precaution has been added to the package insert warning of the possibility of cross-reactivity of non-pertussis species of *Bordetella* with the QIAstat-Dx Respiratory Panel test.

**g. Potentially Interfering Substances:**

An analytical study was performed to assess the potential interference effects of 30 substances naturally present in respiratory specimens or that may be artificially introduced into the nasal cavity/nasopharynx. Samples were tested in triplicate with and without addition of the potentially inhibitory substance for direct sample-to-

sample comparison. Sample mixes were prepared in a simulated matrix consisting of UTM plus HeLa cells to an initial concentration of 10x LOD. Following addition of an equal volume of UTM (control) or interferent the final test concentration was 5x LOD for all panel targets. All pathogen-containing samples without spiked substances generated positive signals for all pathogens present in the sample mixes. Positive influenza signals were generated when influenza vaccine (Fluenz Tetra and FluMist) was tested as a potential interferent. In addition to causing false positive results due the presence of live influenza virus in the vaccine ‘interferent’, it was observed that Fluenz Tetra and FluMist are capable of causing false negative results for QIAstat-Dx RP targets. A precaution has been added to the package insert warning that influenza vaccine present in the patient specimen may cause erroneous results. No other potentially interfering substances tested in this study were found to affect the accuracy of target detection for the QIAstat-Dx Respiratory Panel.

**Table 23: QIAstat-Dx Respiratory Panel Interfering Substances Tested**

<b>Substance</b>	<b>Concentration</b>
Human genomic DNA	20 ng/μL
Whole Blood	1% (v/v)
Mucin	1% (v/v)
Tobramycin	0.6 mg/mL
Mupirocin	2% (w/v)
Saline nasal spray	1% (v/v)
Afrin nasal spray	1% (v/v)
Petroleum jelly	1% (w/v)
Vicks <sup>®</sup> Analgesic ointment	1% (w/v)
Fluenz Tetra nasal vaccine	0.00001% v/v
FluMist live influenza vaccine	0.00001% v/v
Bleach	5% (v/v)

Additional interference was tested using the same methods for specimen collection materials. For swabs, the specific type of swab was dipped in UTM prior to mixing with the pathogen mix. For transport media (VTM), the pathogen mixes were prepared in the specified media at 100% (i.e.; complete replacement of the UTM used as part of the standard setup). The materials and conditions tested are presented in the table below.

**Table 24: QIAstat-Dx Respiratory Panel Collection Material Interference Tested**

Collection Material	Condition
Swab – Copan 168C	Swab +1mL UTM
Swab – Copan FloQ	Swab +1mL UTM
Swab – Copan 175KS01	Swab +1mL UTM
Swab - Puritan	Swab +1mL UTM
VTM Sigma Virocult	100%
VTM Remel M4-RT	100%
VTM Remel M4	100%
VTM Remel M5	100%
VTM Remel M6	100%
BD Universal Viral Transport	100%

No interference was observed using the collection materials listed above.

*h. Microbial Interference:*

QIAstat-Dx Respiratory Panel testing was performed in the presence of non-panel respiratory pathogens to determine whether the pathogens are capable of interfering with the detection of the panel targets. Samples were tested in triplicate with and without addition of the potentially inhibitory organism for direct sample-to-sample comparison. Sample mixes were prepared in a simulated matrix consisting of UTM plus HeLa cells to an initial concentration of 10x LOD. Following addition of an equal volume of UTM (control) or interferent the final test concentration was 5x LOD for all panel targets. All pathogen-containing samples without spiked interferent generated positive signals for all pathogens present in the sample mixes. The potentially interfering organisms did not inhibit detection of any panel targets at the concentrations tested. No false positive results occurred because of the presence of the interfering organism during testing. A list of the potential interfering organisms and their final test concentrations are listed in the table below.

**Table 25: List of Potentially Interfering Organisms and Concentrations Tested**

Interferent Tested	Interferent Concentration
<i>Staphylococcus aureus</i>	1.00 x 10 <sup>6</sup> CFU/mL
<i>Neisseria meningitidis</i>	5.0 x 10 <sup>4</sup> CFU/mL
<i>Corynebacterium diphtheriae</i>	5.0 x 10 <sup>3</sup> CFU/mL
Cytomegalovirus	1.00 x 10 <sup>5</sup> TCID <sub>50</sub> /mL

i. *Carry-over:*

Not applicable. The QIAstat-Dx Respiratory Panel Assay consists of single-use disposable cassettes containing all the reagents, reservoirs, and reaction chambers necessary to perform the test.

j. *Assay cut-off:*

Assay cut-off for the QIAstat-Dx Respiratory Panel was set according to data obtained during LOD studies. The assay software Result Call Algorithm (RCA) is the software processing fluorescence measurements during PCR cycling to produce a qualitative result. The RCA parameters were initially set during system development and later adjusted using empirical data obtained during the LOD studies. The refined RCA was then confirmed during clinical evaluation with the results of the method comparison testing.

2. Comparison studies:

a. *Method comparison with predicate device:*

Not applicable. Performance of the QIAstat-Dx Respiratory Panel was evaluated against the comparator method in prospective and retrospective clinical studies and with contrived specimens where necessary.

b. *Matrix comparison:*

A matrix equivalency study was performed to verify that the initial LOD determinations performed in simulated matrix were comparable to the LOD values obtained in clinical matrix. In order to assess the performance in clinical matrix, a concentration of 1x LOD (as determined in simulated matrix) for at least one strain per respiratory panel (RP) pathogen was prepared in true negative clinical NPS sample matrix and tested in 20 replicates. Results of the matrix equivalency study (as determined empirically by an LOD confirmation experiment) are shown in the table below.

**Table 26: Results of Matrix Equivalency Study**

<b>Pathogen</b>	<b>Strain</b>	<b>Detection Rate<sup>1</sup></b>	<b>Equivalent?</b>
Influenza B	B/Florida/4/2006	18/20	-
Influenza A H3N2	A/Port Chalmers/1/73	Flu A: 20/20 H3: 19/20	+ +
Coronavirus 229E	n/a	20/20	+
Coronavirus OC43	n/a	20/20	+
Coronavirus NL63	n/a	20/20	+
Coronavirus HKU1	n/a	20/20	+
Parainfluenza Virus 1	n/a	18/20	-



Parainfluenza Virus 2	Greer	17/20	-
Parainfluenza Virus 3	C 243	20/20	+
Parainfluenza Virus 4	M-25	20/20	+
Rhinovirus	A2	19/20	+
Enterovirus	US/IL/14-18952	19/20	+
Adenovirus	GB	19/20	+
RSV B	CH93(18)-18	19/20	+
hMPV	hMPV-16	20/20	+
<i>Mycoplasma pneumoniae</i>	PI 1428	20/20	+
<i>Bordetella pertussis</i>	I028	19/20	+
<i>Chlamydomphila pneumoniae</i>	TW183	20/20	+

<sup>1</sup>Detection Rate in clinical matrix at LOD concentration determined in simulated matrix.

The results of the matrix equivalency study demonstrated that the two matrices are not equivalent for all analytes detected by the QIAstat-Dx Respiratory Panel and that analytical studies performed in simulated matrix (at the LOD determined in simulated matrix) represent a more challenging analyte concentration for those studies. Claimed LOD concentrations listed in the LOD section of this summary and the package insert represent either the concentration confirmed in this study or the higher (more concentrated) analyte level as verified in clinical matrix according to minimum empirical acceptance criteria ( $\geq 95\%$  detection rate with a minimum of 20 replicates).

### 3. Clinical Studies:

#### Prospective Specimens

A multi-center study was conducted at five study sites located throughout the U.S. plus one international site between December 2017 and April 2019. The QIAstat-Dx Respiratory Panel was used to evaluate fresh, prospectively collected nasopharyngeal swab specimens eluted in UTM from children and adults of all ages presenting with flu-like symptoms and meeting inclusion/exclusion criteria. Retrospective (archived) samples were also included as part of the performance testing. Each study location was representative of the intended use setting for the QIAstat-Dx Respiratory Panel assay and testing was performed by trained clinical laboratory personnel. A residual NPS specimen in UTM was tested for each subject with the QIAstat-Dx Respiratory Panel and an FDA cleared multiplex respiratory pathogen panel comparator.

A total of 2341 nasopharyngeal swab specimens in UTM were enrolled in the study. Of those, 37 specimens did not meet eligibility criteria or produced invalid results upon repeat testing. A total of 2304 nasopharyngeal swab specimens were considered evaluable. Of the 2304 specimens that met eligibility criteria, 310 were retrospective (archived) specimens and 1994 were fresh prospective or frozen prospective samples. Patient age and gender distribution for the evaluable specimens is presented in tables 27 and 28 below.

**Table 27: Prospective Clinical Study Participant Demographics by Age**

Age Group (Years)	Count	Percent
<5	627	31.4%
6-21	239	12.0%
22-49	330	16.5%
>50	798	40.0%
Total	1994	100%

**Table 28: Clinical Study Participant by Gender and Site**

Site	Female	Male
1	232	186
2	0	0
3	230	196
4	271	177
5	133	170
6	204	195
Total	1070 (53.7%)	924 (46.3%)

Of the 1994 evaluable prospective specimens, 95.88% (1912/1994) yielded valid results on the first attempt (i.e., first loaded cartridge). Invalid or no result were obtained for the remaining 82 specimens (4.11 %). Forty-two (42) specimens were invalid due to cartridge internal control failure (2.11%). Of these, 20 (1.00%) provided a result for positively detected targets and 22 (1.10%) provided a negative result. For 40 (2.00%) specimens no results were obtained due to incomplete runs. Of these, 1 specimen was aborted by users (0.05%), 21 were due to instrument errors (1.05%), and 18 were due to cartridge related errors (0.90%). Seventy-two (72) of the 82 initially failed (no result or invalid) specimens yielded valid results after a single retesting using a new cartridge/sample. The remaining 10 specimens failed on the second attempt (2 due to cartridge failures, 1 due to instrument errors, and 7 due to internal control failures). Of these internal control failures, detected pathogens were reported for 4 specimens.

Compared to an FDA-cleared molecular assay, the performance of the QIAstat-Dx Respiratory Panel for NPS swabs eluted in UTM is presented below.

**Table 29: QIAstat-Dx Respiratory Panel Eluted Nasopharyngeal Swab Performance Compared to FDA-cleared Molecular Comparator – Prospective Specimens**

Analyte	Group	TP/(TP+FN)	PPA	95% CI	TN/(TN+FP)	NPA	95% CI
Adenovirus <sup>1</sup>	Fresh	55/58	94.8%	85.9-98.2	833/839	99.3%	98.4-99.7
	Frozen	31/32	96.9%	84.3-99.4	1047/1057	99.1%	98.3-99.5
	Overall	86/90	95.6%	89.1-98.3	1880/1896	99.2%	98.6-99.5
Coronavirus 229E	Fresh	8/9	88.9%	56.5-98.0	886/886	100%	99.6-100
	Frozen	0/0	n/a	n/a	1089/1089	100%	99.6-100
	Overall	8/9	88.9%	56.5-98.0	1975/1975	100%	99.8-100
Coronavirus HKU1 <sup>2</sup>	Fresh	3/3	100%	43.8-100	890/892	99.8%	99.2-99.9
	Frozen	48/49	98.0%	89.3-99.6	1035/1040	99.5%	98.9-99.8
	Overall	51/52	98.1%	89.9-99.7	1925/1932	99.6%	99.3-99.8
Coronavirus NL63 <sup>3</sup>	Fresh	4/5	80.0%	37.6-96.4	890/890	100%	99.6-100
	Frozen	36/42	85.7%	72.2-93.3	1046/1048	99.8%	99.3-99.9
	Overall	40/47	85.1%	72.3-92.6	1936/1938	99.9%	99.6-100
Coronavirus OC43 <sup>4</sup>	Fresh	3/3	100%	43.8-100	892/892	100%	99.6-100
	Frozen	23/26	88.5%	71.0-96.0	1059/1063	99.6%	99.0-99.9
	Overall	26/29	89.7%	73.6-96.4	1951/1955	99.8%	99.5-99.9
Human <sup>5</sup> metapneumovirus	Fresh	62/67	92.5%	83.7-96.8	828/829	99.9%	99.3-100
	Frozen	53/55	96.4%	87.7-99.0	1030/1034	99.6%	99.0-99.8
	Overall	115/122	94.3%	88.6-97.2	1858/1863	99.7%	99.4-99.9
Rhinovirus/ Enterovirus <sup>6</sup>	Fresh	144/157	91.7%	86.3-95.1	715/739	96.8%	95.2-97.8
	Frozen	124/137	90.5%	84.4-94.4	941/953	98.7%	97.8-99.3
	Overall	268/294	91.2%	87.4-93.9	1656/1692	97.9%	97.1-98.5
Influenza A <sup>7</sup>	Fresh	132/133	99.2%	95.8-99.9	753/757	99.5%	98.6-99.8
	Frozen	110/111	99.1%	95.1-99.8	972/977	99.5%	98.8-99.8
	Overall	242/244	99.2%	97.0-99.8	1725/1734	99.5%	99.0-99.7
Influenza A H1 <sup>8</sup>	Fresh	0/1	0.0%	0.0-79.3	894/894	100%	99.6-100
	Frozen	0/0	n/a	n/a	1089/1089	100%	99.6-100
	Overall	0/1	0.0%	0.0-79.3	1983/1983	100%	99.8-100
Influenza A H1N1/pdm09 <sup>9</sup>	Fresh	62/63	98.4%	91.5-99.7	826/831	99.4%	98.6-99.7
	Frozen	18/18	100%	82.4-100	1071/1071	100%	99.6-100
	Overall	80/81	98.8%	93.3-99.8	1897/1902	99.7%	99.4-99.9
Influenza A H3 <sup>10</sup>	Fresh	67/67	100%	94.5-100	825/826	99.9%	99.3-100
	Frozen	89/90	98.9%	82.4-100	992/998	99.4%	98.7-99.7
	Overall	156/157	99.4%	93.3-99.8	1817/1824	99.6%	99.2-99.8
Influenza B <sup>11</sup>	Fresh	64/67	95.5%	87.6-98.5	827/892	99.9%	99.3-100
	Frozen	58/62	93.5%	84.6-97.5	1026/1026	100%	99.6-100
	Overall	122/129	94.6%	89.2-97.3	1853/1854	99.9%	99.7-100
Parainfluenza virus 1 <sup>12</sup>	Fresh	3/3	100%	43.8-100	892/892	100%	99.6-100
	Frozen	13/14	92.9%	68.5-98.7	1072/1075	99.7%	99.2-99.9
	Overall	16/17	94.1%	73.0-99.0	1964/1967	99.8%	99.6-99.9
Parainfluenza virus 2	Fresh	2/2	100%	34.2-100	893/893	100%	99.6-100
	Frozen	0/0	n/a	n/a	1089/1089	100%	99.6-100
	Overall	2/2	100%	34.2-100	1982/1982	100%	99.8-100
Parainfluenza virus 3 <sup>13</sup>	Fresh	102/104	98.1%	93.3-99.5	788/793	99.4%	98.5-99.7
	Frozen	9/9	100%	70.1-100	1081/1081	100%	99.6-100
	Overall	111/113	98.2%	93.8-99.5	1869/1874	99.7%	99.4-99.9
Parainfluenza virus 4 <sup>14</sup>	Fresh	3/3	100%	43.8-100	892/892	100%	99.6-100
	Frozen	0/0	n/a	n/a	1087/1089	99.8%	99.3-99.9
	Overall	3/3	100%	43.8-100	1979/1981	99.9%	99.6-100

Respiratory Syncytial Virus <sup>15</sup>	Fresh	73/76	96.1%	88.9-98.6	819/820	99.9%	99.3-100
	Frozen	139/144	96.5%	92.1-98.5	941/945	99.6%	98.9-99.8
	Overall	212/220	96.4%	93.0-98.1	1760/1765	99.7%	99.3-99.9
<i>Bordetella pertussis</i> <sup>16</sup>	Fresh	2/2	100%	34.2-100	893/893	100%	99.6-100
	Frozen	1/1	100%	20.7-100	1082/1088	99.4%	98.8-99.7
	Overall	3/3	100%	43.8-100	1975/1981	99.7%	99.3-99.9
<i>Chlamydomphila pneumoniae</i> <sup>17</sup>	Fresh	4/4	100%	51.0-100	891/891	100%	99.6-100
	Frozen	1/1	100%	20.7-100	1087/1088	99.9%	99.5-100
	Overall	5/5	100%	56.6-100	1978/1979	99.9%	99.7-100
<i>Mycoplasma pneumoniae</i> <sup>18</sup>	Fresh	18/18	100%	82.4-100	875/877	99.8%	99.2-100
	Frozen	1/1	100%	20.7-100	1085/1088	99.7%	99.2-99.9
	Overall	19/19	100%	83.2-100	1960/1965	99.7%	99.4-99.9

<sup>1</sup> Adenovirus was detected in 3/4 FN specimens using an independent molecular method. Adenovirus was detected in 6/16 FP specimens using an independent molecular method.

<sup>2</sup> The single FN specimen was negative for Coronavirus HKU1 when tested using an independent molecular method. Coronavirus HKU1 was detected 0/7 FP specimens using an independent molecular method.

<sup>3</sup> Coronavirus NL63 was detected in 7/7 FN specimens using an independent molecular method. Coronavirus NL63 was detected in 1/2 FP specimens using an independent molecular method.

<sup>4</sup> The 3 FN specimens were negative for Coronavirus OC43 when tested using an independent molecular method. Coronavirus OC43 was detected in 3/4 FP specimens using an independent molecular method.

<sup>5</sup> Human metapneumovirus (hMPV) was detected in 4/7 FN specimens using an independent molecular method. hMPV was detected in 3/5 FP specimens using an independent molecular method.

<sup>6</sup> Rhinovirus was detected in 18/26 FN specimens using an independent molecular method. Rhinovirus was detected in 14/36 FP specimens using an independent molecular method.

<sup>7</sup> Influenza A was detected in 1/2 FN specimens by an independent molecular method. Three (3) FP samples were not available for testing. Influenza A was detected in the 3/6 remaining FP samples by an independent molecular method.

<sup>8</sup> Influenza A H1 was detected in 1/1 FN specimen by an independent molecular method. Note: Pre-2009 H1 has not been in circulation since being replaced by 2009 H1 and the discrepant results are likely a false positive result.

<sup>9</sup> Influenza A H1N1/pdm09 was detected in 1/1 FN by an independent molecular method. Influenza A H1 was detected in 3/5 FP specimens by an independent molecular method.

<sup>10</sup> Influenza A H3 was detected in 1/1 FN by an independent molecular method. Influenza H3 was detected in 7/7 FP specimens by an independent molecular method.

<sup>11</sup> Influenza B was detected in 6/6 FN specimens available for testing by an independent molecular method; one discordant sample was not tested by an independent molecular method. Influenza B was detected in 1/1 FP specimens available for testing by an independent molecular method.

<sup>12</sup> 1/1 FN specimens tested negative for Parainfluenza virus 1 by an independent molecular method. Parainfluenza Virus 1 was detected in 3/3 FP specimens by an independent molecular method.

<sup>13</sup> Parainfluenza Virus 3 was detected in 1/2 FN specimens by an independent molecular method. Parainfluenza 3 was detected in 3/5 FP specimens by an independent molecular method.

<sup>14</sup> Parainfluenza Virus 4 was detected in 2/2 FP specimens by an independent molecular method.

<sup>15</sup> Respiratory Syncytial Virus was detected in 2/8 FN specimens by an independent molecular method. Respiratory Syncytial Virus was detected in 3/5 FP specimens by an independent molecular method.

<sup>16</sup> *Bordetella pertussis* was detected in 1/6 FP specimens by an independent molecular method.

<sup>17</sup> *Chlamydomphila pneumoniae* was detected in 1/1 FP specimens by an independent molecular method.

<sup>18</sup> *Mycoplasma pneumoniae* was detected in 1/4 specimens by an independent molecular method.

The QIAstat-Dx Respiratory Panel detected a total of 191 specimens with distinctive multiple organism detections (9.6% of all specimens) in the prospective study. The three organisms most prevalent in multiple detections by the QIAstat-Dx Respiratory Panel in the prospective study were Rhinovirus/Enterovirus (108/191, 56.5%), Respiratory Syncytial Virus (77/191, 40.8%), and Adenovirus (53/191, 27.7%). A summary of the coinfections is described in Table 30.

**Table 30: Coinfections Detected in Prospective Samples During the Clinical Study**

	<b>Detected Coinfections</b>	<b>Discrepant Coinfections</b>
Total Coinfections	191	51
Total Double Coinfections	166	42
Total Triple Coinfections	22	7
Total Quadruple Coinfections	3	2

Archived Specimens

Some of the analytes on the QIAstat-Dx Respiratory Panel were of low prevalence and were not encountered in sufficiently large numbers during the prospective study to adequately demonstrate clinical performance. To supplement the results of the prospective clinical study, an evaluation of preselected frozen archived retrospective specimens was performed. The specimens selected for testing had previously tested positive for one of the following targets at the clinical laboratory by their standard of care method: *Bordetella pertussis*, Coronavirus 229E, Coronavirus HKU1, Coronavirus NL63, Coronavirus OC43, Influenza A H1N1 2009, *Mycoplasma pneumoniae*, *Chlamydomphila pneumoniae*, Parainfluenza Virus 1, Parainfluenza Virus 2, and Parainfluenza Virus 4. Testing was performed by operators who were blinded as to the expected test result. A total of 310 clinical frozen archived retrospective sample were tested by both the comparator method and QIAstat-Dx Respiratory Panel. If the comparator method did not confirm the preselected target as positive, it was excluded from the data analysis for that target. Demographic information for the archived specimens is presented below. Archived samples consisted of 158 male and 152 female specimens.

**Table 31: Participant Age for Archived Clinical Samples**

<b>Age Group (Years)</b>	<b>Count</b>	<b>Percent</b>
<5	139	44.9%
6-21	85	27.4%
22-49	53	17.1%
>50	33	10.7%
Total	310	100%

**Table 32: QIAstat-Dx Respiratory Panel Eluted Nasopharyngeal Swab Performance Compared to FDA-cleared Molecular Comparator – Archived Specimens**

Analyte	TP/(TP+FN)	PPA	95% CI	TN/(TN+FP)	NPA	95% CI
Adenovirus <sup>1</sup>	9/9	100%	70.1-100	297/304	97/8%	95.4-98.9
Coronavirus 229E	26/27	96.3%	81.7-99.3	286/286	100%	98.7-100
Coronavirus HKU1 <sup>2</sup>	14/14	100%	78.5-100	298/299	99.7%	98.1-99.9
Coronavirus NL63 <sup>3</sup>	24/24	100%	86.2-100	286/288	99.3%	97.5-99.8
Coronavirus OC43	28/28	100%	87.9-100	282/282	100%	98.6-100
Human <sup>4</sup> metapneumovirus	2/2	100%	34.2-100	311/311	100%	98.7-100
Rhinovirus/ Enterovirus	44/49	89.8%	78.2-95.5	254/264	96.2%	93.2-97.9
Influenza A	17/17	100%	81.5-100	296/296	100%	98.7-100
Influenza A H1	0/0	n/a	n/a	313/313	100%	98.8-100
Influenza A H1N1/pdm09 <sup>5</sup>	7/8	87.5%	52.9-97.8	304/304	100%	98.8-100
Influenza A H3	8/8	100%	67.5-100	305/305	100%	98.8-100
Influenza B	1/1	100%	20.7-100	312/312	100%	98.8-100
Parainfluenza virus 1	40/40	100%	91.2-100	267/267	100%	98.6-100
Parainfluenza virus 2	3/3	100%	43.8-100	309/309	100%	98.8-100
Parainfluenza virus 3 <sup>6</sup>	1/4	25.0%	4.6-69.9	309/309	100%	98.8-100
Parainfluenza virus 4 <sup>7</sup>	22/24	91.7%	74.2-97.7	278/278	100%	98.6-100
Respiratory Syncytial Virus <sup>8</sup>	11/12	91.7%	64.6-98.5	300/301	99.7%	98.4-99.9
<i>Bordetella pertussis</i>	33/33	100%	89.6-100	261/261	100%	98.5-100
<i>Chlamydomphila pneumoniae</i> <sup>9</sup>	54/61	88.5%	78.2-94.3	250/250	100%	98.5-100
<i>Mycoplasma pneumoniae</i>	25/25	100%	86.7-100	287/288	99.7%	98.1-99.9

<sup>1</sup> Adenovirus was detected in 3/5 FP specimens using an independent molecular method.

<sup>2</sup> The single FP Coronavirus HKU1 specimen was negative when tested using an independent molecular method.

<sup>3</sup> The single FP Coronavirus NL63 specimen was negative when tested using an independent molecular method.

<sup>4</sup> Rhinovirus was detected in 1/2 FN when tested using an independent molecular method. Rhinovirus was detected in 4/10 FP specimens using an independent molecular method.

<sup>5</sup> Influenza H1N1/pdm09 was detected in the single FN specimen.

<sup>6</sup> Parainfluenza Virus 3 was detected in 1/3 FN specimens by an independent molecular method.

<sup>7</sup> Parainfluenza Virus 4 was detected in 1/2 FN specimens by an independent molecular method.

<sup>8</sup> The single FN Respiratory Syncytial Virus was negative for that target by an independent molecular method. The single FP Respiratory Syncytial Virus was negative for that target by an independent molecular method.

<sup>9</sup> *Chlamydomphila pneumoniae* was detected in 4/5 FN specimens by an independent molecular method.

Contrived Specimens

Despite all prospective and retrospective testing efforts, the number of specimens positive for Influenza A H1N1, Parainfluenza 2, Parainfluenza 4, Coronavirus 229E, and *Chlamydomphila pneumoniae* were insufficient to demonstrate system performance. Therefore, contrived specimens were used as surrogate clinical specimens to supplement and test the sensitivity and specificity of the above analytes. Residual negative clinical specimens were spiked with the pathogens at 3x, 5x and 10x LOD levels (50 of each). Contrived samples were provided a unique study identification number and the individual who contrived the samples did not test them; therefore, the status of each contrived specimen was unknown at the time of testing.

**Table 33: QIAstat-Dx Respiratory Panel Performance Compared to FDA-cleared Molecular Comparator – Contrived Specimens**

Analyte	xLOD	TP/(TP+FN)	PPA	95% CI
Influenza A/H1	3	24/24	100%	86.2-100
	5	27/27	100%	87.5-100
	10	24/24	100%	86.2-100
Coronavirus 229E	3	16/16	100%	80.6-100
	5	18/18	100%	82.4-100
	10	16/16	100%	80.6-100
Parainfluenza Virus 2	3	16/16	100%	80.6-100
	5	18/18	100%	82.4-100
	10	16/16	100%	80.6-100
Parainfluenza Virus 4	3	15/16	93.8%	71.7-98.9
	5	18/18	100%	82.4-100
	10	16/16	100%	80.6-100
<i>Chlamydomphila pneumoniae</i>	3	16/16	100%	80.6-100
	5	18/18	100%	82.4-100
	10	16/16	100%	80.6-100

4. Clinical cut-off:

Not applicable.

5. Expected values/Reference range:

In the QIAstat-Dx Respiratory Panel prospective clinical study, a total of 1994 eluted NPS specimens were evaluable by the QIAstat-Dx Respiratory Panel assay. The number and percentage of positive cases per site (Table 34) and per age group (Table 35), as determined by the QIAstat-Dx Respiratory Panel assay, are presented below. One clinical study site consisted of only archived specimens and is not included in the table below.

**Table 34: QIAstat-Dx Respiratory Panel Expected Values by Site**

Organism	Overall (n=1994)		Site 1 (n=418)		Site 2 (n=426)		Site 3 (n=448)		Site 4 (n=303)		Site 5 (n=399)	
	N	Expected value	N	Expected value	N	Expected value	N	Expected value	N	Expected value	N	Expected value
<i>Adenovirus</i>	102	5.1%	44	10.5%	9	2.1%	12	2.7%	30	9.9%	7	1.8%
<i>Coronavirus 229E</i>	8	0.4%	1	0.2%	0	0.0%	0	0.0%	7	2.3%	0	0.0%
<i>Coronavirus HKU1</i>	58	2.9%	4	1.0%	11	2.6%	14	3.1%	12	4.0%	17	4.3%
<i>Coronavirus NL63</i>	42	2.1%	4	1.0%	1	0.2%	15	3.3%	11	3.6%	11	2.8%
<i>Coronavirus OC43</i>	30	1.5%	0	0.0%	5	1.2%	6	1.3%	12	4.0%	7	1.8%
<i>Human Metapneumovirus</i>	120	6.0%	42	10.0%	24	5.6%	14	3.1%	14	4.6%	26	6.5%
<i>Human Rhinovirus/Enterovirus</i>	304	15.2%	59	14.1%	78	18.3%	39	8.7%	53	17.5%	75	18.8%
<i>Influenza A</i>	251	12.6%	120	28.7%	0	0.0%	58	12.9%	38	12.5%	35	8.8%
<i>Influenza A H1</i>	0	0.0%	0	0.0%	0	0.0%	0	0.0%	0	0.0%	0	0.0%
<i>Influenza A H1N1pdm09</i>	85	4.3%	67	16.0%	0	0.0%	4	0.9%	10	3.3%	4	1.0%
<i>Influenza H3</i>	163	8.2%	52	12.4%	0	0.0%	52	11.6%	28	9.2%	31	7.8%
<i>Influenza B</i>	123	6.2%	58	13.9%	0	0.0%	32	7.1%	7	2.3%	26	6.5%
<i>Parainfluenza virus 1</i>	19	1.0%	2	0.5%	1	0.2%	2	0.4%	4	1.3%	10	2.5%
<i>Parainfluenza virus 2</i>	2	0.1%	2	0.5%	0	0.0%	0	0.0%	0	0.0%	0	0.0%
<i>Parainfluenza virus 3</i>	116	5.8%	23	5.5%	19	4.5%	16	3.6%	23	7.6%	35	8.8%
<i>Parainfluenza virus 4</i>	5	0.3%	1	0.2%	0	0.0%	1	0.2%	0	0.0%	3	0.8%
<i>Respiratory Syncytial Virus</i>	217	10.9%	64	15.3%	40	9.4%	35	7.8%	40	13.2%	38	9.5%
<i>Bordetella Pertussis</i>	9	0.5%	2	0.5%	1	0.2%	0	0.0%	6	2.0%	0	0.0%
<i>Chlamydomphila pneumoniae</i>	6	0.3%	2	0.5%	1	0.2%	1	0.2%	1	0.3%	1	0.3%
<i>Mycoplasma pneumoniae</i>	24	1.2%	19	4.5%	0	0.0%	2	0.4%	1	0.3%	2	0.5%

**Table 35: QIAstat-Dx Respiratory Panel Expected Values by Age Group**

Organism	Overall (n=1994)		<=5 yrs (n=626)		6-21 yrs (n=240)		22-49 yrs (n=330)		>49 yrs (n=798)	
	N	Expected value	N	Expected value	N	Expected value	N	Expected value	N	Expected value
<i>Adenovirus</i>	102	5.1%	78	12.4%	7	2.9%	11	3.3%	6	0.8%
<i>Coronavirus 229E</i>	8	0.4%	4	0.6%	4	1.7%	0	0.0%	0	0.0%
<i>Coronavirus HKU1</i>	58	2.9%	29	4.6%	5	2.1%	8	2.4%	16	2.0%
<i>Coronavirus NL63</i>	42	2.1%	25	4.0%	3	1.3%	5	1.5%	9	1.1%
<i>Coronavirus OC43</i>	30	1.5%	20	3.2%	2	0.8%	4	1.2%	4	0.5%
<i>Human Metapneumovirus</i>	120	6.0%	46	7.3%	3	1.3%	17	5.2%	54	6.8%
<i>Human Rhinovirus/Enterovirus</i>	304	15.2%	186	29.7%	35	14.6%	22	6.7%	61	7.6%
<i>Influenza A</i>	251	12.6%	47	7.5%	36	15.1%	64	19.4%	104	13.0%



Organism	Overall (n=1994)		<=5 yrs (n=626)		6-21 yrs (n=240)		22-49 yrs (n=330)		>49 yrs (n=798)	
	N	Expected value	N	Expected value	N	Expected value	N	Expected value	N	Expected value
Influenza A H1	0	0.0%	0	0.0%	0	0.0%	0	0.0%	0	0.0%
Influenza A H1N1 pdm09	85	4.3%	20	3.2%	6	2.5%	30	9.1%	29	3.6%
Influenza H3	163	8.2%	25	4.0%	30	12.6%	35	10.6%	73	9.1%
Influenza B	123	6.2%	11	1.8%	22	9.2%	27	8.2%	63	7.9%
Parainfluenza virus 1	19	1.0%	11	1.8%	0	0.0%	4	1.2%	4	0.5%
Parainfluenza virus 2	2	0.1%	1	0.2%	0	0.0%	0	0.0%	1	0.1%
Parainfluenza virus 3	116	5.8%	70	11.2%	4	1.7%	6	1.8%	36	4.5%
Parainfluenza virus 4	5	0.3%	4	0.6%	0	0.0%	0	0.0%	1	0.1%
Respiratory Syncytial Virus	217	10.9%	135	21.5%	11	4.6%	17	5.2%	54	6.8%
<i>Bordetella Pertussis</i>	9	0.5%	5	0.8%	2	0.8%	0	0.0%	2	0.3%
<i>Chlamydomphila pneumoniae</i>	6	0.3%	1	0.2%	3	1.3%	2	0.6%	0	0.0%
<i>Mycoplasma pneumoniae</i>	24	1.2%	4	0.6%	6	2.5%	11	3.3%	3	0.4%

**N. Instrument Name:**

QIAstat-Dx Analyzer 1.0

**O. System Descriptions:**

1. Modes of Operation:

The QIAstat-Dx Analyzer 1.0 is a portable bench-top unit designed to perform in laboratory and point-of care environments. QIAstat-Dx Respiratory Panel cartridges are keyed to match the corresponding holder on the instrument and in only one direction. Once the cartridge has been inserted into the instrument, the test starts automatically and runs for approximately 74 minutes. When the test is finished, the cartridge is removed by the user and discarded. The QIAstat- Dx Analyzer 1.0 automatically interprets test results and displays a summary on the analyzer display screen. The results can be printed using a connected printer. The Analyzer is equipped with six optical measurement channels employing six fluorescence modules, FAM, NED, ROX, VIC, Cy5, and Cy5.5. Each optical module consists of a LED excitation source and photodiode receiver configured to the fluorescent channel. Fluorescence signals are processed by software running on an on-board computer.

Does the applicant’s device contain the ability to transmit data to a computer, webserver, or mobile device?

Yes  or No

Does the applicant's device transmit data to a computer, webserver, or mobile device using wireless transmission?

Yes \_\_\_\_\_ or No  X

2. Software:

FDA has reviewed applicant's Hazard Analysis and software development processes for this line of product types:

Yes  X  or No \_\_\_\_\_

3. Specimen Identification:

Specimen ID is entered using on screen keyboard or barcode scanner.

4. Specimen Sampling and Handling:

The specimens are manually inserted in the Sample Receiver in the instrument.

5. Calibration:

The QIAstat-Dx Analyzer 1.0 is provided factory calibrated and does not require user calibration. The Analyzer 1.0 includes self-check controls to verify the performance of all sensors and actuators and will alert the user in case of failure. The company will maintain calibration of the instrument through calibration checks performed during QIAGEN technical service visits

6. Quality Control:

Quality control is addressed for each specific FDA-cleared assay to be run on the instrument (separately cleared).

**P. Other Supportive Instrument Performance Characteristics Data Not Covered In The "Performance Characteristics" Section above:**

Not applicable.

**Q. Proposed Labeling:**

The labeling supports the finding of substantial equivalence for this device.

**R. Conclusion:**

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.