

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION  
DECISION SUMMARY  
ASSAY AND INSTRUMENT**

**I Background Information:**

**A 510(k) Number**

K213804

**B Applicant**

Roche Molecular Systems, Inc.

**C Proprietary and Established Names**

cobas SARS-CoV-2 Qualitative for use on the cobas 6800/8800 Systems

**D Regulatory Information**

Product Code(s)	Classification	Regulation Section	Panel
QQX	Class II	21 CFR 866.3981 - Device to detect and identify nucleic acid targets in respiratory specimens from microbial agents that cause the SARS-CoV-2 respiratory infection and other microbial agents when in a multi-target test	MI - Microbiology

**II Submission/Device Overview:**

**A Purpose for Submission:**

New Device

**B Measurand:**

SARS-CoV-2 RNA

**C Type of Test:**

Nucleic acid amplification test with paired reporter and quencher fluorescence labeled probes (TaqMan Technology) using fluorescence resonance energy transfer (FRET)

### **III Intended Use/Indications for Use:**

#### **A Intended Use(s):**

See Indications for Use below.

#### **B Indication(s) for Use:**

cobas SARS-CoV-2 Qualitative for use on the cobas 6800/8800 Systems is a real-time RT-PCR test intended for the qualitative detection of nucleic acids from SARS-CoV-2 in nasal and nasopharyngeal specimens collected from symptomatic individuals suspected of COVID-19 by their healthcare provider.

Results are for the detection of SARS-CoV-2 RNA. Positive results are indicative of the presence of SARS-CoV-2 RNA; clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status. Positive results do not rule out bacterial infection or co-infection with other pathogens.

Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for patient management decisions. Results are meant to be used in conjunction with clinical observations, patient history, recent exposures and epidemiological information, and laboratory data, in accordance with the guidelines provided by the relevant public health authorities. cobas SARS-CoV-2 is intended for use by qualified clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and on the use of the cobas 6800/8800 Systems.

#### **C Special Conditions for Use Statement(s):**

Rx - For Prescription Use Only

#### **D Special Instrument Requirements:**


This test is to be used with the Roche cobas 6800/8800 system only.

### **IV Device/System Characteristics:**


#### **A Device Description:**

cobas SARS-CoV-2 Qualitative is based on fully automated sample preparation (nucleic acid extraction and purification) followed by PCR amplification and detection. The cobas 6800/8800 Systems consists of two separate instruments, the cobas 6800 System and the cobas 8800 System, which both are fully integrated and allow users to perform multiple PCR-based *in vitro* nucleic acid amplification tests. The cobas 6800/8800 Systems consist of the sample supply module, the transfer module, the processing module, and the analytic module. Automated data management is performed by the cobas 6800/8800 Systems software(s), which assigns test results for all tests. Results can be reviewed directly on the system screen and printed as a report.

### **cobas 6600 instrument**



### **cobas 8800 instrument**



The main system functionality is provided by two Software components: the cobas 6800/8800 System Software and Assay Specific Analysis Packages (ASAPs). The cobas 6800/8800 System Software provides basic functionality such as a Graphical User Interface (GUI), instrument management, database functionality, report engines, and LIS interfaces. ASAPs are built using a common Software framework and provides the assay test run conditions (sample preparation and PCR parameters), result analysis functionality (result calculation and algorithms), and result report formatting.

Nucleic acid from patient samples and added internal control RNA (RNA IC) molecules are simultaneously extracted. Nucleic acid is released by addition of proteinase and lysis reagent to the sample. The released nucleic acid binds to the silica surface of the added magnetic glass particles. Unbound substances and impurities, such as denatured protein, cellular debris, and potential PCR inhibitors, are removed with subsequent wash steps and purified nucleic acid is eluted from the magnetic glass particles with elution buffer at elevated temperature. External controls (positive and negative) are processed in the same way.

Selective amplification of target nucleic acid from the sample is achieved using target specific forward and reverse primers for ORF1 a/b non-structural region that is unique to SARS- CoV-2. Additionally, a conserved region in the structural protein envelope E-gene were chosen for pan-Sarbecovirus detection. The pan-Sarbecovirus detection sets will also detect SARS-CoV- 2 virus. Selective amplification of RNA IC is achieved using non-competitive sequence specific forward and reverse primers which have no homology with the coronavirus genome. A thermostable DNA polymerase enzyme is used for amplification.

## B Principle of Operation:

The cobas SARS-CoV-2 Qualitative master mix contains detection probes which are specific for the coronavirus type SARS-CoV-2, members of the Sarbecovirus subgenus, and the RNA IC nucleic acid. The coronavirus and RNA IC detection probes are each labeled with unique fluorescent dyes that act as a reporter. Each probe also has a second dye which acts as a quencher. When not bound to the target sequence, the fluorescent signals of the intact probes are suppressed by the quencher dye. During the PCR amplification step, hybridization of the probes to the specific single-stranded DNA template results in cleavage of the probe by the 5' to 3' exonuclease activity of the DNA polymerase resulting in separation of the reporter and quencher dyes and the generation of a fluorescent signal. With each PCR cycle, increasing amounts of cleaved probes are generated and the cumulative signal of the reporter dye increases concomitantly. Each reporter dye is measured at defined wavelengths, which enables simultaneous detection and discrimination of the amplified coronavirus target and the RNA IC. The master mix includes deoxy uridine triphosphate (dUTP), instead of deoxythymidine triphosphate (dTTP), which is incorporated into the newly synthesized DNA (amplicon). Any contaminating amplicons from previous PCR runs are destroyed by the AmpErase enzyme [uracil-N-glycosylase], which is included in the PCR mix, when heated in the first thermal cycling step. However, newly formed amplicons are not destroyed since the AmpErase enzyme is inactivated once exposed to temperatures above 55°C.

Reagents loaded onto the cobas 6800/8800 Systems are stored at appropriate temperatures and their expiration is monitored by the system. The cobas 6800/8800 Systems allow reagents to be used only if all the conditions are met. The system automatically prevents use of expired reagents via encoded radio frequency identification (RFID) tags that are placed on the reagent cassettes, RMC cassettes and reagent bottles. RFID Labels are used to apply information on test specific and generic reagent cassettes, control tube containers and generic reagent bottles utilized on the cobas 6800/8800 Systems.

### Components of the cobas SARS-CoV-2 Qualitative Kit

Components Name	Formulation
<b>Proteinase Solution (PASE)</b>	Tris buffer, < 0.05% EDTA, calcium chloride, calcium acetate, 8% proteinase, glycerol EUH210: Safety data sheet available on request. EUH208: Contains Subtilisin. May produce an allergic reaction.
<b>RNA Internal Control (RNA IC)</b>	Tris buffer, < 0.05% EDTA, < 0.001% non-Sarbecovirus related armored RNA construct containing primer and probe specific primer sequence regions (non-infectious RNA in MS2 bacteriophage), < 0.1% sodium azide

<b>Elution Buffer (EB)</b>	Tris buffer, 0.2% methyl-4 hydroxybenzoate
<b>Master Mix Reagent 1 (MMX-R1)</b>	Manganese acetate, potassium hydroxide, <0.1% sodium azide
<b>SARS-CoV-2 Master Mix Reagent 2 (SARS-CoV-2 MMX-R2)</b>	Tricine buffer, potassium acetate, < 18% dimethyl sulfoxide, glycerol, < 0.1% Tween 20, EDTA, < 0.12% dATP, dCTP, dGTP, dUTPs, < 0.01% upstream and downstream SARS-CoV-2 and Sarbecovirus primers, < 0.01% Internal Control forward and reverse primers, < 0.01% fluorescent-labeled oligonucleotide probes specific for SARS-CoV-2, Sarbecovirus, and the RNA Internal Control, < 0.01% oligonucleotide aptamer, < 0.1% Z05D DNA polymerase, < 0.10% AmpErase (uracil-N-glycosylase) enzyme (microbial), < 0.1% sodium azide

### Components of the cobas SARS-CoV-2 Qualitative Control Kit

Components Name	Formulation
<b>SARS-CoV-2 Positive Control (SARS-CoV-2 (+) C)</b>	Tris buffer, < 0.05% Sodium azide, < 0.005% EDTA, < 0.003% Poly rA, < 0.01% Non-infectious plasmid DNA (microbial) containing SARS-CoV-2 sequence, <0.01% Non-infectious plasmid DNA (microbial) containing pan-Sarbecovirus sequence

### Components of the cobas Buffer Negative Control Kit

Components Name	Formulation
<b>cobas Buffer Negative Control (BUF (-) C)</b>	Tris buffer, < 0.1% sodium azide, EDTA, < 0.002% Poly rA RNA (synthetic)

Automated data management is performed by the Assay Specific Analysis Package software which assigns test results for all tests. Results can be review directly on the system screen, exported, or printed as a report.

The results interpretation for detecting SARS-CoV-2 is summarized below.

### cobas SARS-CoV-2 Qualitative Results Interpretation

Target 1 (SCoV2)	Target 2 (PanSarb)	Interpretation
<b>Positive</b>	<b>PanSarb Positive</b>	All Target Results were valid. Result for SARS-CoV-2 RNA is Detected.
<b>Positive</b>	<b>Negative</b>	All Target Results were valid. <b>Result for SARS-CoV-2 RNA is Detected.</b> A positive Target 1 result and a negative Target 2 result is suggestive of 1) a sample at concentrations near or below the limit of detection of the test, 2) a mutation in the Target 2, target region, or 3) other factors.

<b>Negative</b>	<b>Positive</b>	All Target Results were valid. <b>Result for SARS-CoV-2 RNA is Presumptive Positive.</b> A negative Target 1 result and a positive Target 2 result is suggestive of 1) a sample at concentrations near or below the limit of detection of the test, 2) a mutation in the Target 1 target region in the oligo binding sites, or 3) infection with some other Sarbecovirus (e.g., SARS-CoV or some other Sarbecovirus previously unknown to infect humans), or 4) other factors. For samples with a Presumptive Positive result, additional confirmatory testing may be conducted, if it is necessary to differentiate between SARS-CoV-2 and SARS-CoV-1 or other Sarbecovirus currently unknown to infect humans, for epidemiological purposes or clinical management.
<b>Negative</b>	<b>Negative</b>	All Target Results were valid. <b>Result for SARS-CoV-2 RNA is Not Detected.</b>
<b>Positive</b>	<b>Invalid</b>	Not all Target Results were valid. <b>Result for SARS-CoV-2 RNA is Detected.</b>
<b>Invalid</b>	<b>Positive</b>	Not all Target Results were valid. <b>Result for SARS-CoV-2 RNA is Presumptive Positive.</b> For samples with a Presumptive Positive result, additional confirmatory testing may be conducted, if it is necessary to differentiate between SARS-CoV-2 and SARS-CoV-1 or other Sarbecovirus currently unknown to infect humans, for epidemiological purposes or clinical management.
<b>Negative</b>	<b>Invalid</b>	Not all Target Results were valid. <b>Sample should be retested.</b> If the result is still invalid, a new specimen should be obtained.
<b>Invalid</b>	<b>Negative</b>	Not all Target Results were valid. <b>Sample should be retested.</b> If the result is still invalid, a new specimen should be obtained.
<b>Invalid</b>	<b>Invalid</b>	All Target Results were invalid. <b>Sample should be retested.</b> If the result is still invalid, a new specimen should be obtained.

**C Instrument Description Information:**

1. Instrument Name: cobas 6800/8800 Systems and software version 1.4.7; cobas SCoV2-QL ASAP Software version 12.1.0
2. Specimen Identification: Specimen identification can be configured in an automated fashion or entered manually. All corresponding user interactions are recorded in the audit trail on the Instrument Gateway
3. Specimen Sampling and Handling: Nasopharyngeal swab (NPS) or Anterior Nasal swab (ANS) specimens collected in transport media.
4. Calibration: Not Applicable

5. Quality Control: The RNA Internal Control (RNA-IC), used to monitor the entire sample preparation and PCR amplification process, is introduced into each specimen during sample processing. In addition, the test utilizes external controls- cobas SARS-CoV-2 Qualitative Control Kit (low titer positive control) and cobas Buffer Negative Control Kit (negative control), provided separately.

**V Substantial Equivalence Information:**

**A Predicate Device Name(s):**

BioFire COVID-19 Test 2

**B Predicate 510(k) Number(s):**

K211079

**C Comparison with Predicate(s):**

<b>Device &amp; Predicate Device(s):</b>	<u>K213804</u>	<u>K211079</u>
<b>Device Trade Name</b>	cobas SARS-CoV-2 Qualitative	BioFire COVID-19 Test 2
<b>Regulation Number</b>	21 CFR 866.3981	same
<b>Regulation Name</b>	Device to detect and identify nucleic acid targets in respiratory specimens from microbial agents that cause the SARS-CoV-2 respiratory infection and other microbial agents when in a multi-target test	same
<b>Product Code</b>	QQX	same
<b>Intended Use</b>	<p>cobas SARS-CoV-2 Qualitative for use on the cobas 6800/8800 Systems is a real-time RT-PCR test intended for the qualitative detection of nucleic acids from SARS-CoV-2 in nasal and nasopharyngeal specimens collected from symptomatic individuals suspected of COVID-19 by their healthcare provider.</p> <p>Results are for the detection of SARS-CoV-2 RNA. Positive results are indicative of the presence of SARS-CoV-2 RNA; clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status. Positive results do not rule out bacterial infection or co-infection with other pathogens.</p> <p>Negative results do not preclude SARS-CoV-2 infection and should</p>	<p>The BioFire COVID-19 Test 2 is a qualitative nested multiplexed RT-PCR in vitro diagnostic test intended for use with the BioFire FilmArray 2.0 and BioFire FilmArray Torch Systems. The BioFire COVID-19 Test 2 detects nucleic acids from severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) in nasopharyngeal swabs (NPS) from symptomatic individuals suspected of COVID-19 by their healthcare provider.</p> <p>Results are for the identification of SARS-CoV-2 RNA. The SARS-CoV-2 RNA is generally detectable in NPS specimens during the acute phase of infection. Positive results are indicative of the presence of SARS-CoV-2</p>



	not be used as the sole basis for patient management decisions. Results are meant to be used in conjunction with clinical observations, patient history, recent exposures and epidemiological information, and laboratory data, in accordance with the guidelines provided by the relevant public health authorities. cobas SARS-CoV-2 is intended for use by qualified clinical laboratory personnel specifically instructed and trained in the techniques of real-time PCR and on the use of the cobas 6800/8800 Systems.	RNA; clinical correlation with patient history and other diagnostic information is necessary to determine patient infection status. Positive results do not rule out co-infection with other pathogens. Results are meant to be used in conjunction with other clinical, epidemiologic, and laboratory data, in accordance with the guidelines provided by the relevant public health authorities. The BioFire COVID-19 Test 2 is intended for use by trained medical and laboratory professionals in a laboratory setting or under the supervision of a trained laboratory professional.
<b>Conditions for use</b>	For prescription use	same
<b>Sample Types</b>	Nasopharyngeal swab specimen Nasal swab specimen	Nasopharyngeal swab specimen
<b>Analyte Targets</b>	SARS-CoV-2	SARS-CoV-2
<b>Sample Preparation Procedure</b>	Automated by cobas 6800/8800 Systems	Automated by BioFire FilmArray 2.0 or BioFire FilmArray Torch systems
<b>Amplification Technology</b>	Real-time PCR	Nested multiplex RT-PCR
<b>Analyte</b>	RNA	RNA
<b>Detection Chemistry</b>	Paired reporter and quencher fluorescence labeled probes (TaqMan Technology) using fluorescence resonance energy transfer (FRET)	Two Step Nested multiplex PCR: Reverse transcription, followed by a multiplexed first stage PCR reaction (PCR1). Multiple simultaneous second-stage PCR reactions (PCR2) to amplify sequences within the PCR1 products using fluorescence double stranded binding dye. Endpoint melting curve data to detect target-specific amplicons
<b>Controls used</b>	1. The RNA Internal Control (RNA-IC) 2. External Positive and negative controls	Two process controls: 1. RNA Process Control (IC) 2. PCR2 Control (A positive result indicates that PCR2 was successful)
<b>Result Analysis</b>	Based on PCR cycle threshold analysis	Endpoint melting curve data to detect target-specific amplicons
<b>Test Interpretation</b>	Automated test interpretation and report generation. User cannot access raw data.	same



## VI Standards/Guidance Documents Referenced:

Class II Special Controls as per 21 CFR 866.3981.

## VII Performance Characteristics:

### A Analytical Performance:

#### 1. Precision/Reproducibility:

##### a) Within-laboratory Precision:

Within-laboratory precision was examined using a panel of heat-inactivated SARS-CoV-2 (USA-WA1/2020) cultures diluted in simulated clinical matrix stabilized in UTM. Sources of variability were examined with a panel consisting of three concentration levels (0.3x, 1x and 3x LoD), using three lots of cobas SARS-CoV-2 Qualitative reagents and three instruments over a course of 15 testing days (2 runs/day x 3 instruments x 5 days/instrument) for a total of 30 runs. A total of 90 replicates per concentration level were tested over the course of the study. A negative panel was tested along with the positive panel members. The tables below summarize the results from the inter-lot precision study, showing the mean Ct values and calculated SD and % CV for each panel member.

#### Qualitative Summary of Within Laboratory Precision

Target	Panel Member	Level (x LoD)	Total Results (N)	Positive Results	Positivity %	Two-sided 95% CI Lower Bound	Two-sided 95% CI Upper Bound
Target 1 (SARS-CoV-2)	High negative	~0.3x	90	9	10%	5%	18%
	Low positive	~1.0x	90	82	91%	83%	96%
	Moderate positive	~3.0x	90	90	100%	96%	100%
Target 2 (pan-Sarbecovirus)	High negative	~0.3x	90	31	34%	25%	45%
	Low positive	~1.0x	90	84	93%	86%	97%
	Moderate positive	~3.0x	90	90	100%	96%	100%
N/A	Negative	Blank	90	0	0%	0%	4%

#### Summary of the Precision Studies based on Ct values

Target	Level (x LoD)	Hit rate	Mean Ct	Within Run		Run-to-Run		Day-to-Day		Instrument -to- Instrument		Lot-to-Lot		Total	
				SD	SD	SD	CV%	CV%	CV%	SD	CV%	SD	CV%	SD	CV%
	~0.3x	10.0%	32.51	0.5	0.0	0.0	0.0	0.0	1.4	0.0	0.0	0.0	0.0	0.5	1.4
	~1.0x	91.1%	32.1	0.6	0.0	0.1	0.3	0.0	1.8	0.0	0.0	0.2	0.6	0.6	1.9

Target	Level (x LoD)	Hit rate	Mean Ct	Within Run		Run-to-Run		Day-to-Day		Instrument -to- Instrument		Lot-to-Lot		Total	
				SD	SD	SD	CV%	CV%	CV%	SD	CV%	SD	CV%	SD	CV%
Target 1 (SARS-CoV-2)	~3.0x	100.0%	31.18	0.3	0.0	0.0	0.0	0.0	0.9	0.0	0.0	0.2	0.7	0.4	1.1
Target 2 (pan-Sarbecovirus)	~0.3x	34.4%	35.36	0.5	0.1	0.3	0.8	0.2	1.5	0.0	0.0	0.5	1.3	0.8	2.2
	~1.0x	93.3%	34.21	0.7	0.0	0.2	0.6	0.0	2	0.0	0.0	0.1	0.3	0.7	2.2
	~3.0x	100.0%	32.9	0.4	0.0	0.0	0.0	0.0	1.1	0.0	0.1	0.0	0.0	0.4	1.1

b) Reproducibility:

The reproducibility of cobas SARS-CoV-2 Qualitative was evaluated at three testing sites, using three reagent lots, with a four-member panel of positive and negative samples resulting in a total number 216 tests per concentration (not including controls). The positive panel members contained the SARS-CoV-2 target in negative simulated clinical matrix stabilized in UTM at three different concentrations (0.3x, 1x, 3x LoD). Each site tested two reagent lots for six days. Two runs were performed each day and three replicates of each panel member were performed for each run. An overall SARS-CoV-2 positive result was determined by a positive detection in either or both SARS-CoV-2 or/and pan-Sarbecovirus channels. The study results are summarized in Table below.

**Summary of Reproducibility Variance Components based on Ct values**

Viral Target	Panel Member Concentration	n <sup>a</sup> /N	Percent Agreement* (%) <sup>b</sup>	Mean Ct	Between Site		Between Lot		Between Day		Between Runs		Within Run		Total SD	Total CV(%)
					SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)	SD	CV(%)		
Negative	0	214/216 <sup>c</sup>	99.1	nc	nc	nc	nc	nc	nc	nc	nc	nc	nc	nc	nc	nc
SARS-CoV-2	~0.3xLoD	45/216	20.8	33.6	0.00	0.0	0.00	0.0	0.11	0.3	0.00	0.0	0.35	1.1	0.37	1.1
SARS-CoV-2	~1xLoD	196/216	90.7	33.2	0.00	0.0	0.09	0.3	0.00	0.0	0.17	0.5	0.37	1.1	0.42	1.3
SARS-CoV-2	~3xLoD	216/216	100.0	32.2	0.05	0.2	0.02	0.1	0.00	0.0	0.03	0.1	0.24	0.8	0.25	0.8
pan-Sarbecovirus	~0.3xLoD	158/216	73.1	36.5	0.18	0.5	0.00	0.0	0.00	0.0	0.00	0.0	0.71	2.0	0.74	2.0
pan-Sarbecovirus	~1xLoD	214/216	99.1	35.4	0.00	0.0	0.00	0.0	0.00	0.0	0.00	0.0	0.67	1.9	0.67	1.9
pan-Sarbecovirus	~3xLoD	216/216	100.0	34.1	0.11	0.3	0.05	0.2	0.00	0.0	0.00	0.0	0.32	0.9	0.34	1.0

Ct = cycle threshold, LoD = limit of detection, SD = standard deviation, CV(%) = percent coefficient of variation, SARS CoV 2 = severe acute respiratory syndrome coronavirus 2, nc = not calculable

Note: SARS-CoV-2 is a dual target assay. Inactivated viral culture material was diluted to ~0.3/1/3x LoD based on the target 2 (SARS-CoV-2) LoD.

<sup>a</sup> n is the number of positive tests which contribute Ct values to the analysis. N is the total number of valid tests for the panel member.

<sup>b</sup> Percent agreement with expected results.

<sup>c</sup> Two negative panel members were tested positive (one for each target). Sequencing showed that one of these samples was positive and the other was negative. The Ct values and the curve analysis of these samples suggest a low level of contamination during specimen handling.

2. Linearity:

Not applicable; this is a qualitative assay.

3. Analytical Specificity/Interference:

a) Cross-Reactivity:

The potential for cross-reactivity (analytical specificity) for cobas SARS-CoV-2 Qualitative assay was evaluated in a study testing samples containing organisms that are commonly found in the respiratory tract and might cause similar clinical symptoms as SARS-CoV-2. A panel of 47 viruses, bacteria, and fungi were included in the study, with each organism tested in three replicates. The organisms were spiked at concentrations of  $1 \times 10^5$  units/mL for viruses and  $1 \times 10^6$  units/mL for bacteria and fungi. Testing was performed with each potential interfering organism in the absence and presence of SARS-CoV-2 target (spiked at  $\sim 3 \times \text{LoD}$ ). No cross reactivity or microbial interference was observed with the organisms at the concentrations tested, as shown below. Additional *in-silico* analysis confirmed no homology was found between any primers or probes for *Bordetella parapertussis*, *Aspergillus sp*, and *Fusobacterium necrophorum*, however no wet testing with these organisms was performed.

**Summary of Analytical Cross-Reactivity for cobas SARS-CoV-2 Qualitative**

Microorganism	Concentration	Positivity Rate	Negativity rate
Human coronavirus 229E	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Human coronavirus OC43	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Human coronavirus HKU1	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Human coronavirus NL63	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
MERS coronavirus	1.0E+05 genomic equivalent/mL	100%	100%
SARS coronavirus	1.0E+05 PFU/mL	100%	100%*
Adenovirus B (Type 34)	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Bocavirus	1.0E+05 cp/mL	100%	100%
Cytomegalovirus	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Epstein Barr virus	1.0E+05 cp/mL	100%	100%
Human Metapneumovirus (hMPV)	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Measles virus	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Mumps virus	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Parainfluenza virus Type 1	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Parainfluenza virus Type 2	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Parainfluenza virus Type 3	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Parainfluenza virus Type 4	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Influenza A (H1N1)	1.0E+05 TCID <sub>50</sub> /mL	100%	100%

Microorganism	Concentration	Positivity Rate	Negativity rate
Influenza A virus (H1N1-2009)	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Influenza A virus (H1N3)	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Influenza A virus (H3N2)	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Influenza B	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Enterovirus E (Type 1)	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Parechovirus	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
Respiratory syncytial virus	1.0E+05 PFU/mL	100%	100%
Rhinovirus	1.0E+05 TCID <sub>50</sub> /mL	100%	100%
<i>Candida albicans</i>	1.0E+06 CFU/mL	100%	100%
<i>Chlamydia pneumoniae</i>	1.0E+06 TCID <sub>50</sub> /mL	100%	100%
<i>Corynebacterium diphtheriae</i>	1.0E+06 CFU/mL	100%	100%
<i>Escherichia coli</i>	1.0E+06 CFU/mL	100%	100%
<i>Haemophilus influenzae</i>	1.0E+06 CFU/mL	100%	100%
<i>Lactobacillus gasseri</i>	1.0E+06 CFU/mL	100%	100%
<i>Legionella pneumophila</i>	1.0E+06 CFU/mL	100%	100%
<i>Legionella jordanis</i> (non-pneumophila)	1.0E+06 CFU/mL	100%	100%
<i>Moraxella catarrhalis</i>	1.0E+06 CFU/mL	100%	100%
<i>Mycobacterium tuberculosis</i>	1.0E+06 cells/mL	100%	100%
<i>Neisseria elongata</i>	1.0E+06 CFU/mL	100%	100%
<i>Neisseria meningitidis</i>	1.0E+06 CFU/mL	100%	100%
<i>Pseudomonas aeruginosa</i>	1.0E+06 CFU/mL	100%	100%
<i>Pneumocystis jirovecii</i>	1:20 of Patient Sample	100%	100%
<i>Staphylococcus aureus</i>	1.0E+06 CFU/mL	100%	100%
<i>Staphylococcus epidermidis</i>	1.0E+06 CFU/mL	100%	100%
<i>Streptococcus pneumoniae</i>	1.0E+06 CFU/mL	100%	100%
<i>Streptococcus pyogenes</i>	1.0E+06 CFU/mL	100%	100%
<i>Streptococcus salivarius</i>	1.0E+06 CFU/mL	100%	100%
<i>Bordetella pertussis</i>	1.0E+06 CFU/mL	100%	100%
<i>Mycoplasma pneumoniae</i>	1.0E+06 CFU/mL	100%	100%

\* Testing of SARS-CoV-1 generated an expected Pan-Sarbecovirus positive result.

b) Interfering substances study:

The potential for interference for the cobas SARS-CoV-2 Qualitative assay was evaluated with substances that may be present in respiratory specimens. A total of 20 potentially interfering substances were tested at or above clinically relevant levels in negative simulated clinical matrix stabilized in UTM in absence and presence of SARS-CoV-2 target (spiked at

~3x LoD). Each sample was tested in 10 replicates. The following table provides the substances and concentrations tested. The FluMist nasal vaccine was not tested.

### Summary of Interfering Substances Study for Cobas SARS-CoV-2 Qualitative

Name	Active Ingredients	Test concentration	Total N	Positivity %	
				Target 1 (SARS-CoV-2)	Target 2 (Pan-Sarbeco)
Control	No Substance	N/A	10	100	100
Human Genomic DNA	-	20 ng/µl	10	100	100
Mucus	-	One sputum swab/mL	10	100	100
Human PBMCs	-	1.0E+03 cells/µL	10	100	100
Human Whole Blood	-	1%, 2% and 5% (v/v)	10	100	100
Afrin Nasal Spray	Oxymetazoline	0.011 mg/mL	10	100	100
Budesonide Nasal spray	Budesonide (glucocorticoid)	0.039 mg/mL	10	100	100
Flonase Nasal Spray	Fluticasone Propionate	0.167 mg/mL	10	100	100
Zicam Nasal Spray	Galphimia glauca, Luffa operculata, Sabadilla	0.023 mg/mL	10	100	100
Saline Nasal Spray with Preservatives	0.65% NaCl, phenylcarbino, Benzalkonium chloride	1.0% (v/v)	10	100	100
Cepacol (Sore Throat Lozenges)	Benzocaine and menthol	5.0 mg/mL	10	100	100
Chloraseptic	Phenol	0.47 mg/mL	10	100	100
Liposomal NUMB520 Spray	Lidocaine and Phenylephrine	2.68 mg/mL	10	100	100
Mupirocin ointment UPS	Mupirocin	0.20 mg/mL	10	100	100
Relenza (Inhalation powder)	Zanamivir	0.0015 mg/mL	10	100	100
Antiviral drug – Tamiflu	Oseltamivir	0.0073 mg/mL	10	100	100
DMSO	N/A	N/A	10	100	100
Tobramycin ophthalmic Solution	Tobramycin	0.018 mg/mL	10	100	100
Vaseline	Petroleum Jelly	1% (w/v)	10	100	100
Snuff Tobacco	Nicotine	1% (w/v)	10	100	100
Analgesic ointment (Vicks@VapoRub)	Camphor-synthetic eucalyptus oil and menthol ointment	1% (w/v)	10	100	100

\* FluMist was not evaluated to assess potential interference.

#### 4. Assay Reportable Range:

Not applicable; this is a qualitative assay.

#### 5. Traceability, Stability, Expected Values (Controls, Calibrators, or Methods):

##### a) Quality Controls

The RNA Internal Control (RNA-IC), used to monitor the entire sample preparation and PCR amplification process, is introduced into each specimen during sample processing. In

addition, the test utilizes external controls- cobas SARS-CoV-2 Qualitative Control Kit (low titer positive control) and cobas Buffer Negative Control Kit (negative control), provided separately. Three lots of Positive Control were included in the clinical evaluation of the cobas SARS-CoV-2 Qualitative conducted at three testing sites. The summary of the Ct values for valid positive and negative control results obtained for each target, combined for sites and lots, is presented below.

**Summary of the Ct values for Controls used in Clinical Studies**

Target Analyte	Total N	Mean Target Ct	Ct SD	Ct %CV
SARS-CoV-2	131	33.66	0.21	0.66
pan-Sarbeco	131	35.56	0.34	1.00
IC (positive samples)	131	33.23	0.22	0.66
IC (negative samples)	131	33.1	0.24	0.7

b) Clinical Specimen stability:

The stability of SARS-CoV-2 viral targets during handling, transportation, and storage, was evaluated using NS and NPS specimens collected in different collection media. Five pools of NS specimens (collected in cobas PCR media, 0.9% saline and copan UTM) and five pools of NPS specimens (collected in UTM) were spiked with the inactivated quantified SARS-CoV-2 virus at 3xLoD. Three replicates per specimen type/storage condition were used. The study data supported storage for specimens collected in UTM and 0.9% physiological saline at 2-25°C for up to 48 hours followed by up to 3 days at 2-8°C and at ≤ -70°C for up to 30 days. Specimens collected in cobas PCR Media can be stored at 2-25°C for up to 24 hours followed by up to 3 days at 2-8°C and at ≤ -70°C for up to 30 days.

6. Detection Limit:

The Limit of Detection (LoD) for cobas SARS-CoV-2 Qualitative was determined using an inactivated quantified SARS-CoV-2 virus (WHO International Standard for SARS-CoV-2, NIBSC code: 20/146). LoD is defined as the lowest concentration of SARS-CoV-2 RNA that can be detected at a rate of at least 95%. A total of 5 concentration levels (500, 250, 125, 62.5 and 31.25 IU/ml) were prepared by diluting the SARS-CoV-2 target in negative simulated clinical matrix stabilized in UTM. Three independent dilution series with three lots of reagents were tested with a total of 24 replicates per concentration.

The concentration level with observed hit rates greater than or equal to 95% was determined to be the LoD for each of the two targets (SARS-CoV-2 and pan-Sarbecovirus) as described below.

### Summary of LoD for SARS-CoV-2 (Target-1)

Viral Strain	Kit lot	Hit rate ≥ 95% [IU/mL] (no. of replicates detected)	Mean Ct at ≥ 95% Hit rate
WHO International Standard for SARS-CoV-2	Lot 1	250 (24/24)	33.2
	Lot 2	125 (23/24)	34.1
	Lot 3	250 (23/24)	33.2

The LoD was confirmed at 250 IU/mL for SARS-CoV-2 (Target 1). For all three reagent lots, at least 23/24 replicates detected the target at 250 IU/ml.

### Summary of LoD for pan-Sarbecovirus (Target-2)

Viral Strain	Kit lot	Hit rate ≥ 95% [IU/mL] (no. of replicates detected)	Mean Ct at ≥ 95% Hit rate
WHO International Standard for SARS-CoV-2	Lot 1	125 (24/24)	35.2
	Lot 2	125 (24/24)	36.0
	Lot 3	125 (23/24)	34.8

The LoD was confirmed at 125 IU/mL for pan-Sarbecovirus (Target 2). For all three reagent lots, at least 23/24 replicates detected the target at 125 IU/ml.

#### 7. Matrix equivalency study:

Because some analytical studies were conducted using simulated matrix, equivalency between simulated and real clinical matrix was evaluated. SARS-CoV-2 target was diluted to make the test panel of four concentration levels (0x, 0.3x, 1x and 3x LoD) in each of three matrices (simulated clinical matrix, nasopharyngeal swab specimen and nasal swab specimen). 25 replicates were tested per each combination.

#### Summary of matrix equivalency study for cobas SARS-CoV-2 Qualitative

Sample Type Matrix (Stabilized in UTM)	Panel Member	Level (x LoD)	Conc. IU/mL	Total Results (N)	Positive Results	Positivity %
Nasal Swab	Above LoD	~3x	750	25	25	100%
	At LoD	~1.0x	250	25	25	100%
	Below LoD	~0.3x	83	25	25	100%
	Negative	0x	0	25	0	0%
	Above LoD	~3x	750	25	25	100%
	At LoD	~1.0x	250	25	25	100%



Nasopharyngeal Swab	Below LoD	~0.3x	83	25	25	100%
	Negative	0x	0	25	0	0%
Simulated Clinical Matrix	Above LoD	~3x	750	25	25	100%
	At LoD	~1.0x	250	25	25	100%
	Below LoD	~0.3x	83	25	24	96%
	Negative	0x	0	25	0	0%

The results showed the performance of the assay with simulated clinical matrix, nasopharyngeal swab, and nasal swab matrices is equivalent.

#### 8. Collection media equivalency (UTM, cobas PCR Media and 0.9% physiological saline)

Equivalency between copan UTM, cobas PCR Media (CPM) and 0.9% physiological saline was evaluated during the development of the subject IVD device. Paired individual negative nasal swab samples were collected and tested before and after spiking heat inactivated SARS-CoV-2 virus (WHO International Standard). At least 20 low positive (1.5-2x LoD) and 10 moderate positive (4x LoD) samples was prepared for each collection media. All negative and positive samples generated expected results for SARS-CoV-2 in the three collection media tested. The mean Ct differences between the three collection media was less than 2.

#### 9. Analytical Reactivity (Inclusivity):

The inclusivity of cobas SARS-CoV-2 Qualitative assay for the detection of SARS-CoV-2 was confirmed by testing nine SARS-CoV-2 strains, including six variant strains. Qualified stocks of different SARS-CoV-2 were used to prepare concentration levels at and around the LoD. Each level was tested with four replicates. The lowest concentration at which all four replicates tested positive are reported in below.

#### **Analytical Reactivity with SARS-CoV-2 Strains**

<b>SARS-CoV-2 Strain</b>	<b>Target Material</b>	<b>Concentration tested with 100% Positivity</b>
Hong Kong/VM20001061/2020	Inactivated Virus	1.06E+02 cp/mL
Italy-INMI1	Inactivated Virus	1.00E+02 cp/mL
USA-WA1/2020	Inactivated Virus	5.03E+01 cp/mL
UK (B.1.1.7)	Inactivated Virus	2.4E+01 cp/mL
Japan / Brazil (P.1)	Extracted RNA	1.9E+02 cp/mL
South Africa (B.1.351)	Inactivated Virus	2.4E+01 cp/mL
US NY (B.1.526)	Inactivated Virus	1.9E+02 cp/mL
India (B.1.617.1)	Inactivated Virus	2.5E+02 cp/mL
India (B.1.617.2)	Inactivated Virus	7.0E+01 cp/mL

*In silico* analysis of additional SARS-CoV-2 sequences indicates that >99.9% of known sequences for SARS-CoV-2 have no changes in primer/probe binding sites at both target regions simultaneously. All sequences are predicted to be detected by at least one of the two target regions.

***In silico* Analysis of cobas SARS-CoV-2 Test Oligo Design as of October 15, 2022**

Assay	Orf1ab (Target 1)				E-gene (Target 2)				Orf1ab & E-gene (Dual Targets)			
	NCBI		GISAID		NCBI		GISAID		NCBI		GISAID	
<b>total</b>	6113594	100%	12695489	100%	6113594	100%	12695489	100%	6113594	100%	12695489	100%
<b>With mismatch</b>	251910	4.12%	464658	3.66%	32124	0.53%	83852	0.66%	1125	0.02%	2971	0.02%
<b>dCp&gt;5 or Tm&lt;65</b>	211	0.00%	493	0.00%	34	0.00%	1025	0.01%	0	0.00%	3	0.00%

***In Silico* Analysis of SARS-CoV-2 Variants of Concern (VOC) as of October 10, 2022**

WHO label	Status	Pango Lineage	Nextstrain Clade	cobas SARS CoV-2 for use on the cobas 6800/8800 Systems
Alpha	de-escalated	B.1.1.7	20I/501Y.V1	detected
Beta	de-escalated	B.1.351	20H/501Y.V2	detected
Gamma	de-escalated	P.1	20J/501Y.V3	detected
Delta	Active VOC	B.1.617.2	21A,I,J	detected
Epsilon	de-escalated	B.1.427/B.1.429	20C/S:452R	detected
Zeta	de-escalated	P.2	20J	detected
Eta	de-escalated	B.1.525	20A/S:484K	detected
Theta	de-escalated	P.3	21E	detected
Iota	de-escalated	B.1.526	20C/S:484K	detected
Kappa	de-escalated	B.1.617.1	20A/S:154K	detected
Lambda	de-escalated	C.37	20D	detected
Mu	de-escalated	B.1.621	21H	detected
Omicron (BA.1 - BA.5, BQ.1, BQ.1.1, BF.7)	Active VOC	B.1.1.529	21K,L,M 22A,B,C,D,E	detected

10. Assay Cut-Off:

The assay cut-off for this device was established based on variety of parameter thresholds taken together by the algorithm to make the qualitative decision whether a curve is to be considered positive or negative. These parameters (e.g., Ct, max. RFI, max. slope, sigmoidality of curve shape) are initially set by the algorithm. As the product development progresses, the data is annotated, and the algorithm is adjusted (“trained”) accordingly. This iterative process in the end leads to a final set of parameter cut-offs. Selected data from testing of clinical specimens or cultures spiked into clinical background were used to confirm adequate separation between the latest Ct value observed in positive specimens in each target detection channel and the assay cut-off.

## 11. Carry-Over:

Carry-over/cross contamination rate for cobas SARS-CoV-2 Qualitative assay was evaluated by testing SARS-CoV-2 high-positive samples and negative samples in a checkerboard configuration. Five runs, with each run containing 48 replicates of SARS-CoV-2 negative samples and 46 replicates of high positive samples were performed. A total of 230 SARS-CoV-2 high positive samples and 239 negative samples were tested. No cross-contamination was observed, as described in the table below.

### Cross Contamination Run Results

Run Number	Total Positive Results from Negative Samples	Total Samples			Cross Contamination Rate
		Positive	Negative	Controls	
Run 1	0	46	47	2	0.0 %
Run 2	0	46	48	2	0.0 %
Run 3	0	46	48	2	0.0 %
Run 4	0	46	48	2	0.0 %
Run 5	0	46	48	2	0.0 %
<b>Total</b>	<b>0</b>	<b>230</b>	<b>239</b>	<b>10</b>	<b>0.0 %</b> (upper one-sided 95% CI of 1.25%)

## B Comparison Studies:

### 1. Method Comparison with Predicate Device:

Not applicable. Refer to the Clinical Studies Section of this document.

## C Clinical Studies:

### Clinical Study Design

The performance of cobas SARS-CoV-2 Qualitative was evaluated in a multi-center study with three external testing sites evaluating prospectively collected clinical specimens in UTM-RT or UVT from individuals with signs and symptoms of respiratory infection. Participants from 12 geographically distributed enrollment centers each provided nasopharyngeal swab (NPS) and nasal swab (NS, anterior nares) specimens as part of a dual collection where (a) the collection order (first specimen collected) was alternated between the NPS and NS specimen, and (b) the collection method for NS specimens was also alternated with 50% of the NS specimens were self-collected on-site with healthcare provider (HCP) instructions while the other 50% were collected by the healthcare provider. The study used a composite comparator method wherein laboratory sites used up to three highly sensitive EUA SARS-CoV-2 molecular assays, testing NPS specimen from each subject. The composite comparator result was defined as the concordant results from two comparator assays (test A and test B). In case of discordance between the initial two comparator assays, the sample was tested by a third assay (test C) and the result of the third test determined the composite comparator result. The composite comparator

result was indetermined when valid results could not be obtained from two assays (i.e., insufficient volume for repeat testing of invalid/failed results).

From March to June 2021, a total of 1,154 participants were enrolled, of which samples from 968 participants were included in the evaluation. Samples from 186 participants were not included: 184 specimens were excluded due to issues associated with specimen shipments and/or being unable to complete testing within the times identified by manufacturer’s instructions, and two subjects were excluded for being previously enrolled in the study (exclusion criteria). When self-reporting COVID-19 vaccination status, 207 (21.4%) of the 968 participants were fully vaccinated.

Of the 968 participants, 961 contributed a NPS specimen which resulted in 942 participants with a confirmed infected status. For NPS, 4 specimens had failed/invalid cobas SARS-CoV-2 Qualitative results, resulting in 938 evaluable NPS results. For NS, 8 specimens were invalid/missing cobas SARS-CoV-2 Qualitative NS results, resulting in 934 evaluable results.

When compared with the NPS composite comparator result, cobas SARS-CoV-2 Qualitative yielded a positive percent agreement (PPA) of 98.7% for NPS and 97.4% for NS specimens. The negative percent agreement (NPA) was 99.7% and 99.9% for NPS and NS specimens, respectively.

**Summary of clinical performance of cobas SARS-CoV-2 Qualitative for nasopharyngeal (NPS) and nasal swabs (NS) versus the NPS composite comparator**

Specimen Type	Total (N)	PPA	PPA 2-sided 95% Score CI	NPA	NPA 2-sided 95% Score CI
Nasopharyngeal (NPS)	938	98.7% (77/78)	(93.1 %, 99.8 %)	99.7% (857/860)	(99.0 %, 99.9 %)
Nasal Swab (NS)*	934	97.4% (76/78)	(91.1 %, 99.3 %)	99.9% (855/856)	(99.3 %, 100 %)

\*Healthcare provider-collected nasal swab specimens and nasal swab specimens self-collected on-site with healthcare provider instructions

**D Clinical Cut-Off:**

Not applicable

**E Expected Values/Reference Range:**

The SARS-CoV-2 positivity determined by the cobas SARS-CoV-2 Qualitative is presented by enrollment site in table below for both NPS and NS specimens. Across the 12 enrollment sites from 8 states, the SARS-CoV-2 positivity ranged from 0% to 16.2% for NPS and from 0% to 15.8% for NS.

**SARS-CoV-2 Prevalence by Collection Site based on cobas SARS-CoV-2 Qualitative**

Site location	NPS results	NPS Prevalence	NS results	NS Prevalence
Oakland, CA	114	0.0% (0/114)	117	0.0% (0/117)

Site location	NPS results	NPS Prevalence	NS results	NS Prevalence
St.Louis, MO	117	11.1% (13/117)	117	10.3% (12/117)
Syracuse, NY	11	9.1% (1/11)	11	9.1% (1/11)
Bronx, NY	64	7.8% (5/64)	64	9.4% (6/64)
Nashville, TN	92	12.0% (11/92)	93	11.8% (11/93)
Easley, SC	127	13.4% (17/127)	128	12.5% (16/128)
Powdersville, SC	117	16.2% (19/117)	114	15.8% (18/114)
Dallas, TX	147	6.8% (10/147)	147	6.8% (10/147)
Dallas, TX	31	0.0% (0/31)	31	0.0% (0/31)
New Orleans, LA	46	4.3% (2/46)	46	2.2% (1/46)
Statesville, NC	33	6.1% (2/33)	34	5.9% (2/34)
Moorseville, NC	42	0.0% (0/42)	42	0.0% (0/42)
ALL	941	8.5% (80/941)	944	8.2% (77/944)

Note: The table included all specimens with valid SARS-CoV-2 results on cobas SARS-CoV-2. A positive SARS-CoV-2 is defined based on the cobas SARS-CoV-2. NPS = nasopharyngeal swab specimen; NS = nasal swab specimen.

#### **F Other Supportive Instrument Performance Characteristics Data:**

Not applicable

#### **VIII Proposed Labeling:**

The labeling supports the finding of substantial equivalence for this device.

#### **IX Conclusion:**

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.